

www.weldoninstruments.com

PROVEEDOR EXCLUSIVO:



Cannulas & I/A Instruments

Eye Speculas


Forceps & Scissors

Needle Holders & Retractors

661 - 752

Ophthalmology

Ophthalmology
Oftalmologia
Augenheilkunde



WS-661-001
Fixier Measuring Instrument
4.5 cm / 1 3/4"



Castroviejo
WS-661-002
Marking Caliper
15 mm



Castroviejo
WS-661-003
Marking Caliper
20 mm



Mendez
WS-661-004
11.5 cm / 4 1/2"
Measuring Ring



Maloney
WS-661-005
Keratometer



Modell Giessen
WS-661-006
Scale 1 - 30 mm



Scale
WS-661-007
1-20 cm

Cannulas | Cánulas | Kanülen

Kratz
WS-662-001

Capsule Polishing Cannula


Kratz
WS-662-002

 Irrigating and Polishing Cannula
 20 Gauge / 0.90 x 34 mm

Teping
WS-662-003

 Aspirating Polishing Cannula
 23 Gauge / 0.64 mm

WS-662-004

 Capsule Polishing Cannula
 With Bulbous Tip
 18 Gauge / 1.20 mm

Kyu
WS-662-005

 Capsule Polishing Cannula
 22 Gauge / 0.70 mm

Kyu
WS-662-006

 Capsule Polishing Cannula
 23 Gauge / 0.64 mm

Greite
WS-662-007

 Capsule Polishing Cannula With Spoon
 22 Gauge / 0.70 mm

Brauweiler
WS-662-008

 Irrigating Cannula, Angled
 With Bulbous Tip, Rough
 21 Gauge / 0.80 mm

Busse
WS-662-009

 Cyclodialysis Spatula Cannula
 Length of Spatula 12 mm
 26 Gauge / 0.45 mm

WS-662-010

 Cyclodialysis Spatula Cannula
 Extremely Delicate
 Length of Spatula 10 mm
 26 Gauge / 0.45 mm

Mainz Model
WS-662-011

 Cyclodialysis Spatula Cannula
 Length of Spatula 12 mm
 2 Lateral Irrigating Openings
 20 Gauge / 0.90 mm

WS-662-012


 Irrigating Cystotome Cannula
 26 Gauge / 0.45 mm




Enzyme Cannula
WS-663-013
24 Gauge / 0.60 mm



Sautter
WS-663-014
Hydrodissection Cannula
27 Gauge / 0.40 mm




Sautter
WS-663-015
Hydrodissection Cannula
26 Gauge / 0.45 mm



Tuebingen Model
WS-663-016
Hydrodissection Cannula
30 Gauge / 0.30 mm



Frankfurt Model
WS-663-017
Hydrodissection Cannula
27 Gauge / 0.40 mm



Müller
WS-663-018
Irrigating Cannula
With Probe Tip
26 Gauge / 0.45 x 19 mm
Right



Müller
WS-663-019
Irrigating Cannula
With Prob Tip
26 Gauge / 0.45 x 19 mm
Left



Rentsch
WS-663-020
Irrigating Cannula
26 Gauge / 0.45 x 24 mm



Binkhorst
WS-663-021
Irrigating Cannula
22 Gauge / 0.70 x 36 mm
Right



Binkhorst
WS-663-022
Irrigating Cannula
22 Gauge / 0.7 x 36 mm
Left



Irrigating Cannula
WS-663-023
With Bulbous Tip
Extremely Delicate




Simcoe
WS-663-024
Irrigating/ Aspirating Cannula
Irrigating And Aspirating Tube
23 Gauge / 0.64 mm Each



Irrigating/ Aspirating Cannula
WS-663-025
Round Curve
Irrigating And Aspirating Tube
23 Gauge / 0.64 mm Each




Irrigating/ Aspirating Cannula
WS-663-026 22 Gauge / 0.70 mm Each
WS-663-027 26 Gauge / 0.45 mm Each
Irrigating And Aspirating Tube




Engels
WS-664-028
 Iris Retractor
 Irrigating Cannula Tube
 21 Gauge / 0.8 mm




Brauweiler
WS-664-029
 Aspirating Cannula
 Angled, Bevelled, Rough
 21 Gauge / 0.80 mm




H.-R. Koch
WS-664-030
 Aspirating Cannula
 Angled, Rough
 21 Gauge / 0.80 mm




Brauweiler
WS-664-031
 Irrigating Cannula
 Angled
 21 Gauge / 0.80 mm



Brauweiler
WS-664-032
 Irrigating Cannula
 Angled
 With Bulbous Tip, Rough
 21 Gauge / 0.80 mm



Agarwal
WS-664-033
 Irrigating Nucleus Chopper
 Tube 20 Gauge / 0.90 mm



Bonn Model
WS-664-034
 Lacrimal Cannula
 With Probe Tip



Bangerter
WS-664-035 21 Gauge / 0.80 x 60 mm
WS-664-036 23 Gauge / 0.64 x 60 mm
WS-664-037 21 Gauge / 0.80 x 40 mm
WS-664-038 23 Gauge / 0.64 x 40 mm
 Lacrimal Probe Cannula



Jünemann
 Lacrimal Probe Cannula
WS-664-039
 To Insert Suture Material
 21 Gauge / 0.80 x 35 mm



Jünemann
 Lacrimal Probe Cannula
WS-664-040
 To Insert Suture Material
 21 Gauge / 0.80 x 35 mm



Lacrimal Cannula
 WS-665-041 Silver
 WS-665-042 Stainless Steel
 Conical, Straight



Lacrimal Cannula
 WS-665-043 Silver
 WS-665-044 Stainless Steel
 Conical, Round Curve



Lacrimal Cannula
 WS-665-045 Silver
 WS-665-046 Stainless Steel
 Conical, Short Angle



Lacrimal Cannula
 WS-665-047 Silver
 WS-665-048 Stainless Steel
 Conical, Wide Angle



Lacrimal Cannula
 WS-665-049
 Cylindrical, Straight
 20 Gauge / 0.90 x 35 mm



Anel
 Lacrimal Cannula
 WS-665-050
 Straight
 23 Gauge / 0.64 mm



Anel
 Lacrimal Cannula
 WS-665-051
 Round Curve
 23 Gauge / 0.64 mm



Zirm
 WS-665-052
 Lasik Dissector Cannula
 25 Gauge / 0.5 mm



Zirm
 WS-665-053
 Lasik Irrigating Cannula
 20 Gauge / 0.9 mm



Nagata
 WS-665-054
 Sub-Tenon Injection Cannula
 24 Gauge / 0.60 mm



Dardenne
 WS-665-055
 Positioning Needle Cannula
 30 Gauge / 0.30 mm



Waldemar Bauer
 WS-665-056
 Lasik Cannula
 Double-Armed
 6 Irrigating Ports
 25 Gauge / 0.5 mm



Retrobulbar
 WS-665-057 24 Gauge / 0.55 x 60 mm
 WS-665-058 25 Gauge / 0.50 x 35 mm
 WS-665-059 24 Gauge / 0.55 x 35 mm
 WS-665-060 25 Gauge / 0.50 x 50 mm
 WS-665-061 25 Gauge / 0.50 x 40 mm
 Injection Cannula



WS-666-062
Irrigating Vectis
Cannula



WS-666-063
Irrigating Vectis
Cannula
Small Model



Klöckner
WS-666-064
Expressor Cannula



PH.Jacobi
WS-666-065
Trabecular Aspirating Cannula



Burk
WS-666-066
Irrigating
Trabeculotomy Cannula
Right



Burk
WS-666-067
Irrigating
Trabeculotomy Cannula
Left



Irrigating Cannula
WS-666-068
For Visco Canalostomy
25 Gauge / 0.49 mm
7.0 mm Angled
Tapered Tube 5.0 mm Curved
0.19 mm



Aspirating Cannula
WS-666-069
Delicate
19 Gauge / 1.10 mm



Aspirating Cannula
WS-666-070
17 Gauge / 1.50 mm



Healon Cannula
WS-666-071
21 Gauge / 0.80 x 30 mm



Aspirating Cannula
WS-666-072
Inside: 24 Gauge / 0.55 mm
Outside: 20 Gauge / 0.9 mm



Aspirating Cannula
WS-666-073
Tapered Tube / 3.5 mm
25 Gauge / 0.5 mm



Aspirating Cannula
WS-666-074
Lateral Aspirating Port / 0.5 mm
20 Gauge / 0.9 mm



Laqua
WS-666-075
Aspirating Cannula
Inside: 21 Gauge / 0.8 mm
Outside: 19 Gauge / 1.1 mm



Kampik
WS-667-076
Aspirating Cannula
Inside: 22 Gauge / 0.7 mm
Outside: 20 Gauge / 0.9 mm



Aspirating Cannula
WS-667-077
With Silicone Brush
20 Gauge / 0.9 mm




Aspirating Cannula
WS-667-078
20 Gauge / 0.9 mm
Tapered Tube 7.0 mm
3.5 mm Angled
24 Gauge / 0.60 mm



Aspirating Tip
WS-667-079
20 Gauge / 0.9 mm



Aspirating Tip
WS-667-080
3.0 mm Angled
20 Gauge / 0.9 mm



Aspirating Tip
WS-667-081
With Silicone Tube at The Tip
20 Gauge / 0.9 mm



Aspirating Tip
WS-667-082
3.0 mm Angled
With Silicone Tube at The Tip
20 Gauge / 0.9 mm



Aspirating Tip
WS-667-083
Curved
Lateral Aspirating Port 0.3 mm
24 Gauge / 0.60 mm



Aspirating Tip
WS-667-084
Tapered Tube
5.0 mm Angled
30 Gauge / 0.3 mm
Aspirating Port Oval



Aspirating Tip
WS-667-085
20 Gauge / 0.9 mm
Tapered Tube 2.0 mm
24 Gauge / 0.6 mm



Weber
WS-667-086 5.0 mm
WS-667-087 10.0 mm
Aspirating Cannula
For Viscous Fluid
Inside: 22 Gauge / 0.7 mm
Outside: 20 Gauge / 0.9 mm



Roider
WS-667-088
Aspirating Cannula
For Viscous Fluid
Lateral Aspiration Port 0.7 mm
19 Gauge / 1.0 mm



Pfcl-Cannula
WS-667-089
With Two Tubes
24 Gauge / 0.55 mm
Lateral Tube
25 Gauge / 0.5 mm



Bartz-Schmidt
WS-667-090
20 Gauge / 0.9 mm
2 Lateral Ports 0.5 mm
90° Transposed
Tapered Tube 19 mm
25 Gauge / 0.5 mm
With Probe Tip 0.9 mm
Aspirating Port 0.15 mm



Heidelberg Model
WS-667-091
Injection Cannula
For Viscous Fluid,
1.20 mm, 30°



WS-668-092
Irrigating Cannula Handpiece



Brauweiler WS-668-093
Irrigating Cannula Handpiece




Brauweiler WS-668-094
Aspirating Cannula Handpiece



D.D Koch-T. Kohnen WS-668-095
9 cm / 3 1/2"
Aspirating Handpiece
Aspirating Port 0.30 mm
Rough Tip
Tube 22 Gauge / 0.7 mm



D.D Koch-T. Kohnen WS-668-096
9 cm / 3 1/2"
Irrigating Handpiece
Tube 22 Gauge / 0.7 mm



Aspirating Handpiece WS-668-097
9 cm / 3 1/2"
Aspirating Port 0.25 mm
Conical And Rough Tip
Oval Tube 0.95 x 0.55 mm



Irrigating Handpiece WS-668-098
9 cm / 3 1/2"
Irrigating Port 0.50 mm
Conical Tip
Oval Tube 0.95 x 0.55 mm



WS-668-099
9 cm / 3 1/2"
Aspirating Handpiece
Aspirating Port 0.30 mm
Tube 21 Gauge / 0.80 mm
Conical And Rough Tip



WS-668-100
9 cm / 3 1/2"
Irrigating Handpiece
Irrigating Ports 2 x 0.50 mm
Tube 21 Gauge / 0.80 mm
Conical And Rough Tip



Engels WS-668-101
9 cm / 3 1/2"
Iris Retractor-Irrigating Handpiece
Tube 21 Gauge / 0.8 mm



R. Koch
WS-669-102
 9 cm / 3½"
 Aspirating Port 0.3 mm
 Tube 19 Gauge / 1.1 mm
 Rough Tip



Brauweiler
WS-669-103
 9 cm / 3½"
 Angled
 Declined And Rough
 21 Gauge / 0.80 mm



H.-R. Koch
WS-669-104
 9 cm / 3½"
 Aspirating Handpiece
 Angled, Rough
 21 Gauge / 0.80 mm



Brauweiler
WS-669-105
 9 cm / 3½"
 Irrigating Handpiece
 Angled
 21 Gauge / 0.80 mm



Brauweiler
WS-669-106
 9 cm / 3½"
 Irrigating Handpiece, Angled
 With Bulbous Tip, Rough
 21 Gauge / 0.80 mm



R. Koch
WS-669-107
 9 cm / 3½"
 Irrigating Handpiece With Chopper
 Angled, Irrigating Ports 2 x 0.75 mm
 Tube 19 Gauge / 1.1 mm



Hoffmann
WS-669-108
 9 cm / 3½"
 Irrigating Handpiece With Chopper
 Angled, Irrigating Ports 2 x 0.5 mm
 Oval Tube 1.1 mm x 0.9 mm



Wilhelm
WS-669-109
 9 cm / 3½"
 Irrigating Handpiece With Chopper
 Angled, Irrigating Ports 2 x 0.75 mm
 Tube 19 Gauge / 1.1 mm



Neuhann/Sauder
WS-669-110
 9 cm / 3½"
 Irrigating Handpiece With Chopper
 Angled, Irrigating Ports 2 x 0.75 mm
 Tube 19 Gauge / 1.1 mm



HöH
WS-669-111
 9 cm / 3½"
 Chopper
 18 Gauge / 1.2 mm



Anis
WS-670-112
11.5 cm / 4½"
Ball Capsule Polisher



Rentsch
WS-670-113
11.5 cm / 4½"
Capsule Polishing Curette
Angled, 135°



Rentsch
WS-670-114
11.5 cm / 4½"
Capsule Polishing Curette
Angled, 135°



Meyhöfer
WS-671-001
11.5 cm / 4½"
Chalazion Curette
1 mm



Meyhöfer
WS-671-002
11.5 cm / 4½"
Chalazion Curette
2 mm



Meyhöfer
WS-671-003
11.5 cm / 4½"
Chalazion Curette
2.5 mm



Meyhöfer
WS-671-004
11.5 cm / 4½"
Chalazion Curette
3 mm



Meyhöfer
WS-671-005
11.5 cm / 4½"
Chalazion Curette
4 mm



Meyhöfer
WS-671-006
11.5 cm / 4½"
Chalazion Curette
4.5 mm

Curettes, Spoons | Curetas, Cucharas | Küretten, Löffel



Hebra
WS-672-007
Chalazion Curette
Small



Hebra
WS-672-008
Chalazion Curette
Medium



Hebra
WS-672-009
Chalazion Curette
Large



Axenfeld
WS-672-010
Chalazion Curette
Double Ended



Bunge
WS-672-011
Exenteration Curette
Medium



Bunge
WS-672-012
Exenteration Curette
Large



Rycroft A.C. Cannula
WS-673-013
 20 Gauge / 0.90 x 22 mm
WS-673-014
 25 Gauge / 0.50 x 25 mm
 5 Per Box, Sterile




Kratz Capsule Polisher
WS-673-015
 20 Gauge / 0.60 mm
WS-673-016
 21 Gauge / 0.80 mm
 5 Per Box, Sterile



Cystotome Cannula
WS-673-017
 Side Cutting
 25 Gauge / 0.50 mm
WS-673-018
 25 Gauge / 0.50 mm
 5 Per Box, Sterile



Sautter
WS-673-019
 Hydro Dissection Cannula
 27 Gauge / 0.40 x 22 mm
 5 Per Box, Sterile



A.C. Maintainer
WS-673-020
 20 Gauge / 0.90 mm
 5 Per Box, Sterile



Hydro Dissection Cannula
WS-673-021
 With Flat Tip
 25 Gauge / 0.50 x 20 mm
 5 Per Box, Sterile




Lacrimal Cannula
WS-673-022
 Curved
 26 Gauge / 0.45 x 30 mm
 5 Per Box, Sterile



Irrigating vectis
Cannula
WS-673-023
 Double Hole
 25 Gauge / 0.50 Mm
 5 Per Box, Sterile



Irrigating / Aspirating Cannula
WS-673-024
 21 Gauge / 0.80 X 35 mm
 5 Per Box, Sterile



Simcoe
Irrigating / Aspirating
Cannula
WS-673-025
 Aspiration Port, 0.30 mm
 23/23 Gauge / 0.64 / 0.64 mm
 5 Per Box, Sterile



Retrolbulbar Cannula
WS-673-026
 25 Gauge / 0.50 x 40 mm
 5 Per Box, Sterile



Retrolbulbar cannula
WS-673-027
 Extra Long
 25 Gauge / 0.50 x 50 mm
 5 Per Box, Sterile



Atkinson
WS-673-028
 Retrolbulbar Cannula
 25 Gauge / 0.50 x 40 mm
 5 Per Box, Sterile



Atkinson
WS-673-029
 Retrolbulbar Cannula
 23 Gauge / 0.64 x 40 mm
 5 Per Box, Sterile



Atkinson
WS-673-030
 Peribulbar Cannula
 25 Gauge / 0.50 x 30 mm
 5 Per Box, Sterile



Barraquer
WS-674-001
Eye Speculum
Fenestrated



Barraquer
WS-674-002
Eye Speculum
Solid Blades



Barraquer
WS-674-003
Eye Speculum
Fenestrated
For Premature Infants



Barraquer
WS-674-004
Eye Speculum
Fenestrated, For Infants



Barraquer
WS-674-005
4 cm / 1½"
Eye Speculum
Fenestrated, Round Blades



Barraquer
WS-674-006
Eye Speculum
Fenestrated, Round Blades



Barraquer
WS-674-007
Eye Speculum
Fenestrated, Round Blades
For Infants



Barraquer
WS-674-008
Eye Speculum
Fenestrated, Open Blades



Barraquer
WS-674-009
Frankfurt Model
Eye Speculum
Fenestrated, Large



Kratz-Barraquer
WS-675-010
Eye Speculum
Fenestrated, Round Blades




OHNO
WS-675-011 Right
WS-675-012 Left
Eye Speculum
Fenestrated




H.R.KOCH
WS-675-013
Eye Speculum
For Nasal Placement
Fenestrated
Open Blades




Eye Speculum
WS-675-014
For Nasal Placement
Fenestrated
Round Blades



Kershner
WS-675-015
Eye Speculum
Open Blades
For Nasal And Temporal Placement



Kershner
WS-675-016
Eye Speculum
Solid Blades
For Nasal And Temporal Placement



Liebermann-Tennant
WS-675-017
Eye Speculum
Fenestrated, Open Blades
For Temporal Placement

Eye Speculas, Adjustable | Especulos Oculares, Adjustable | Augenspekula, Justierbar


Barraquer-Oosterhuis
WS-676-018
 Eye Speculum
 Adjustable, Fenestrated



Barraquer-Oosterhuis
WS-676-019
 Eye Speculum
 Adjustable, Fenestrated
 For Infants



Barraquer-Oosterhuis
WS-676-020
 Eye Speculum
 Adjustable, Solid Blades



Barraquer-Oosterhuis
WS-676-021
 Eye Speculum
 Adjustable, Solid Blades
 For Infants



Barraquer-Oosterhuis
WS-676-022
 Eye Speculum
 Adjustable, Fenestrated
 Round Blades



Barraquer-Oosterhuis
WS-676-023
 Eye Speculum
 Adjustable, Fenestrated
 Round Blades, For Infants



Neuhann
WS-676-024 Solid Blades
WS-676-025 Fenestrated
 Eye Speculum
 Adjustable (Not Illustrated)



Barraquer-Oosterhuis-Kratz
WS-676-026
 Eye Speculum
 Adjustable, Fenestrated
 Open Blades




Kratz-Barraquer
WS-676-027
 Eye Speculum
 For Adults



Tazawa
WS-676-028 Right, Adjustable, Fenestrated
WS-676-029 Left, Adjustable, Fenestrated
 Eye Speculum
 Round Blades




Weiss
WS-677-030
Eye Speculum
Fenestrated



Weiss
WS-677-031
Eye Speculum
Fenestrated, Open Blades



Weiss
WS-677-032
Eye Speculum
Fenestrated, For Infants



Weiss
WS-677-033
Eye Speculum
Solid Blades



Weiss
WS-677-034
Eye Speculum
Solid Blades, For Infants

Eye Speculas, Adjustable | Especulos Oculares, Adjustable | Augenspekula, Justierbar



Cook
WS-678-035
Eye Speculum
Large



Cook
WS-678-036
Eye Speculum
Medium



Cook
WS-678-037
Eye Speculum
Small



Murdock
WS-678-038
Eye Speculum
Right



Murdock
WS-678-039
Eye Speculum
Left




Murdock
WS-678-040
Eye Speculum
Small, Right



Murdock
WS-678-041
Eye Speculum
Small, Left



Cook
WS-679-042
Eye Speculum
For Infants



Bangerter
WS-679-043
Eye Speculum



Pretori
WS-679-044 Left
WS-679-045 Right
Eye Specula



Sauer
WS-679-046
Eye Speculum




Takeda
WS-679-047
Eye Speculum
Fenestrated, For Infants




Castroviejo
WS-679-048 Small (14 mm Blades)
WS-679-049 Medium (16 mm Blades)
WS-679-050 Large (18 mm Blades)
Lidsperre

Eye Speculas,Adjustable | Especulos Oculares,Adjustable | Augenspekula,Justierbar




Mellinger
WS-680-051
Eye Speculum




Mellinger
WS-680-052
Eye Speculum
For Infants




Mellinger-Beard
WS-680-053
Eye Speculum




Mellinger-Beard
WS-680-054
Eye Speculum
For Infants



Schott
WS-681-055
Eye Speculum
Fenestrated
Weight 17gm




Schott
WS-681-056
Eye Speculum
Fenestrated, For Infants
Weight 17gm



Schott
WS-681-057
Eye Speculum
Fenestrated, For Infants
Weight 13 gm



Schott
WS-681-058
Eye Specula
Solid Blades
Weight 17gm



Schott
WS-681-059
Eye Speculum
Fenestrated With Aspiration
For Infants
Weight 17gm

Fixation Rings | Anillos de Fijación | Fixierungs-Ringe


Thornton
WS-682-001
10 cm / 4"
Fixation Ring
Inside Diameter 16 mm
12 Teeth



Thornton
WS-682-002
10 cm / 4"
Fixation Ring
Inside Diameter 16 mm
Swivelable 12 Teeth



Clear-Cornea-Incision
WS-682-003
10 cm / 4"
Fixation Ring
Inside Diameter 16 mm
Swivelable 12 Teeth



Clear-Cornea-Incision
WS-682-004
10 cm / 4"
Fixation Ring
Inside Diameter 16 mm
12 Teeth



Clear-Cornea-Incision
WS-682-005
10 cm / 4"
Fixation Ring
Inside Diameter 13 mm
9 Teeth



Fixation Ring
WS-682-006
10 cm / 4"
With Concentric Grooves
In The Interior Part Of The Ring
Inside Diameter 12 mm



Flieringa Fixation Ring
WS-682-007 15 mm
WS-682-008 16 mm
WS-682-009 17 mm
WS-682-010 18 mm
WS-682-011 19 mm
WS-682-012 20 mm
10 cm / 4"



Eye Dressing Forceps
WS-683-001AA
 Delicate, 0.7 mm Wide At Tip
WS-683-001A
 0.8 mm Wide At Tip
 10.0 cm / 4"
 Serrated, Straight



Eye Dressing Forceps
WS-683-001BB
 Delicate, 0.7 mm Wide At Tip
WS-683-001B
 0.8 mm Wide At Tip
 9.9 cm / 3³/₁₆"
 Serrated, Half Curved



Eye Dressing Forceps
WS-683-001CC
 Delicate, 0.7 mm Wide At Tip
WS-683-001C
 0.8 mm Wide At Tip
 9.7 cm / 3³/₁₆"
 Serrated, Full Curved




Iris Forceps
WS-683-001DD
 Delicate, 0.7 mm Wide At Tip
WS-683-001D
 0.8 mm Wide At Tip
 10 cm / 4"
 Straight, 1 x 2 Teeth




Iris Forceps
WS-683-001EE
 Delicate, 0.7 mm Wide At Tip
WS-683-001E
 0.8 mm Wide At Tip
 9.9 cm / 3³/₁₆"
 Half Curved, 1 x 2 Teeth



Iris Forceps
WS-683-001FF
 Delicate, 0.7 mm Wide At Tip
WS-683-001F
 0.8 mm Wide At Tip
 9.7 cm / 3³/₁₆"
 Full Curved, 1 x 2 Teeth



Arruga
WS-683-001
 10.5 cm / 4¹/₂"
 Capsule Forceps
 Calibrated Delicate




WS-683-002
 10.5 cm / 4¹/₂"
 Capsulorhexis Forceps
 Without Ridge
 Extremely Delicate




H.R. Koch
WS-683-003
 10.5 cm / 4¹/₂"
 Capsulorhexis Forceps
 With 2 - 13 mm Gauge



WS-683-004
 10.5 cm / 4¹/₂"
 Capsulorhexis Forceps
 Extremely Delicate




Masket
WS-683-005
 10.5 cm / 4¹/₂"
 Capsulorhexis Forceps



WS-683-006
 10.5 cm / 4¹/₂"
 Capsulorhexis Forceps
 Extremely Delicate



WS-683-007
 10.5 cm / 4¹/₂"
 Capsulorhexis Forceps
 Strong Model



WS-683-008
 10.5 cm / 4¹/₂"
 Capsulorhexis Forceps
 Extremely Delicate
 Round Handle

Forceps | Pinzas | Pinzetten


Lambert
WS-684-009
 10 cm / 4"
 Chalazion Forceps



Ayer
WS-684-010
 10 cm / 4"
 Chalazion Forceps




Desmarres
WS-684-011 20 mm
WS-684-012 25 mm
WS-684-013 31 mm
 9 cm / 3½"
 Chalazion Forceps



Cauer
WS-684-014
 Chalazion Forceps




Tobias Neuhann
WS-684-015
 12.5 cm / 5"
 Multifunction Micro Forceps
 23 Gauge / 0.65 mm




N. Inamura
WS-684-016
 11.5 cm / 4½"
 Capsulorhexis Forceps



Berke
WS-684-017
 9.5 cm / 3¾"
 Ptosis Forceps



Henry
WS-685-018
10 cm / 4"
Cilia Forceps



Gradle
WS-685-019
10 cm / 4"
Cilia Forceps



Barraquer
WS-685-020
10 cm / 4"
Cilia Forceps



Saupe
WS-685-021
10 cm / 4"
Cilia Forceps



Beer
WS-685-022
10 cm / 4"
Cilia Forceps



Bonaccolto
WS-685-023
11 cm / 4 1/8"



Barraquer
WS-685-024
10 cm / 4"
Cilia Forceps




Hamburg Model
WS-685-025
10 cm / 4"
Cilia And Foreign Body Forceps

Forceps | Pinzas | Pinzetten



Von Barsewisch
WS-686-026
8 cm / 3"
Anterior Chamber Forceps



Colibri Forceps
WS-686-027
8 cm / 3"
With Tying Platform




Micro Colibri Forceps
WS-686-028
8 cm / 3"
With Tying Platform



Micro Colibri Forceps
WS-686-029
8 cm / 3"
With Tying Platform-Extremely Delicate



Zuerich Model
WS-686-030
8 cm / 3"
Suturing Forceps
Straight Delicate



Zuerich Model
WS-686-031
8 cm / 3"
Suturing Forceps
Angled Delicate



Von Graefe
WS-687-032A
11.4 cm / 4½",
Plain, 4.5 mm Jaws



Von Graefe
WS-687-032B
11.4 cm / 4½",
With Catch, 4.5 mm Jaws
Fixation Forceps



Von Graefe
WS-687-032CC Serrated
WS-687-032C 1 x 2 Teeth
7.2 cm / 2¾"
Curved, 0.7 mm Tip
Iris Forceps



Von Graefe
WS-687-032DD Serrated
WS-687-032D 1 x 2 Teeth
7.2 cm / 2¾"
Straight, 0.7 mm Tip
Iris Forceps



Von Graefe
WS-687-032EE Serrated
WS-687-032E 1 x 2 Teeth
7.2 cm / 2¾"
Curved, 0.7 mm Tip
Iris Forceps



Krefeld Model
WS-687-032
10 cm / 4"
Conjunctiva Forceps



Freiburg Model
WS-687-033
10 cm / 4"
Conjunctiva Forceps



Kaufmann
WS-687-034
9 cm / 3½"
Conjunctiva Forceps
With Tying Platform



Draeger
WS-687-035
10 cm / 4"
Conjunctiva Forceps
Non-Traumatic



Förster
WS-687-036
10 cm / 4"
Conjunctiva Forceps

Forceps | Pinzas | Pinzetten


Remy
WS-688-037
 8 cm / 3"
 Keratoplasty Forceps



Neuhann
WS-688-038
 8 cm / 3"
 Keratoplasty Forceps



Fixation Forceps
WS-688-039
 10 cm / 4"
 Straight, Width Of Prong 3 mm
 1 x 2 Teeth Each Side, 0.12 mm



Fixation Forceps
WS-688-040
 10 cm / 4"
 Curved, Width Of Prong 3 mm
 1 x 2 Teeth Each Side, 0.12 mm



Kremer
WS-688-041
 11.5 cm / 4½"
 Fixation Forceps
 Curved With Lock, Width of Prong 13 mm
 1 x 2 Teeth Each Side 0.12 mm



Kremer
WS-688-042
 11.5 cm / 4½"
 Fixation Forceps
 Straight With Lock, Width of Prong 13 mm
 1 x 2 Teeth Each Side 0.12 mm



Dressing Forceps-Ophthalmical

- WS-689-043 0.5 mm
- WS-689-044 0.7 mm
- WS-689-045 1 mm
- WS-689-046 1.2 mm
- WS-689-047 1.4 mm
- 10 cm / 4"
- Serrated



Dressing Forceps-Ophthalmical

- WS-689-048 0.5 mm
- WS-689-049 0.7 mm
- WS-689-050 1 mm
- WS-689-051 1.2 mm
- WS-689-052 1.4 mm
- 10 cm / 4"
- 1 x 2 Teeth



Graefe

- WS-689-053 Without Lock
- WS-689-054 With Lock
- 10 cm / 4"
- Fixation Forceps



Heidelberg Model

- WS-689-055
- 10 cm / 4"
- Fixation Forceps



Jess

- WS-689-056
- 10 cm / 4"
- Fixation Forceps



Barraquer

- WS-689-057
- 10 cm / 4"
- Fixation Forceps



Twist Fixation Hook

- WS-689-058 Right
- WS-689-059 Left

Forceps | Pinzas | Pinzetten



Iris Forceps
WS-690-060
9 cm / 3½"
Curved, Serrated



Iris Forceps
WS-690-061
9 cm / 3½"
Straight, Serrated



Bonn Model Iris Forceps
WS-690-062
9 cm / 3½"
Straight, Serrated



Bonn Model Iris Forceps
WS-690-063
9 cm / 3½"
Angled 6 mm, Serrated



Bonn Model Iris Forceps
WS-690-064
9 cm / 3½"
Angled 8 mm, Serrated



Bonn Model Iris Forceps
WS-690-065
9 cm / 3½"
Angled 6 mm, 1 x 2 Teeth



Bonn Model Iris Forceps
WS-690-066
9 cm / 3½"
Angled 8 mm, 1 x 2 Teeth



Bonn Model
WS-691-067
10 cm / 4"
Iris Forceps
Straight, 1 x 2 Teeth




Bonn Model
WS-691-068
10 cm / 4"
Iris Forceps
Long, Straight, 1 x 2 Teeth



Bonn Model
WS-691-069
10 cm / 4"
Iris Forceps
Long, Straight, 1 x 2 Teeth
With Tying Platform



Swiss Model
WS-691-070
10 cm / 4"
Iris Forceps
Curved, 1 x 2 Teeth



Heidelberg Model
WS-691-071
9 cm / 3 1/2"
Iris Forceps
Light Curve, 1 x 2 Teeth



Heidelberg Model
WS-691-072
9 cm / 3 1/2"
Iris Forceps
Strong Curve, 1 x 2 Teeth

Forceps | Pinzas | Pinzetten



Iris Forceps
WS-692-073
8 cm / 3"
Straight, 1 x 2 Teeth




Iris Forceps
WS-692-074
7.5 cm / 3"
Straight, 1 x 2 Teeth, Strong



Hamburg Model Iris Forceps
WS-692-075
9 cm / 3½"
Straight, 1 x 2 Teeth



Swiss Model Iris Forceps
WS-692-076
8 cm / 3"
Extremely Delicate, 1 x 2 Teeth



Elze
WS-692-077
7.5 cm / 3"
Iris Forceps



Hess-Engelking Iris Forceps
WS-692-078
1 x 2 Teeth



Neuhann
WS-693-079
8.5 cm / 3 1/4"
Implantation Forceps



Fechner
WS-693-080
8.5 cm / 3 1/4"
Implantation Forceps



Shepard
WS-693-081 Without Stop
WS-693-082 With Stop
10 cm / 4"
Implantation Forceps



Tenner
WS-693-083
10 cm / 4"
Implantation Forceps



Doden
WS-693-084
10 cm / 4"
Implantation Forceps

Forceps | Pinzas | Pinzetten




Elze
WS-694-085
10 cm / 4"
Implantation Forceps
For Retropupillary Lenses Without Holes



Holding Forceps
WS-694-086
10 cm / 4"
For Foldable IOLs



T.Kohnen
WS-694-087
10 cm / 4"
Implantation Forceps
Silicone Tube Coated




Folder
WS-694-088
10 cm / 4"
For Foldable IOLs




Spaleck
WS-694-089
10 cm / 4"
Folder For Silicone And
Acrylic IOLs Without Lock



Thomas Neuhann
WS-695-090
10 cm / 4"
Folder For Foldable IOLs



WS-695-091
10 cm / 4"
Folder
With Notches For Acrylic IOLs




Nucleus Splitter Forceps
WS-695-092
10 cm / 4"



H.-R.Koch
WS-695-093
10 cm / 4"
For Foldable 3-Piece IOLs



H.-R.Koch
WS-695-094
10 cm / 4"
Plate Haptic Iol Folder
Folding Technique For Plate Haptic



Akahoshi
WS-695-095
10 cm / 4"
Nucleus-Spreading Forceps
Pre-Chopper

Forceps | Pinzas | Pinzetten



Folder
WS-696-096
10 cm / 4"
With Lock
For Pharmacia CEEoN Edge



Explantation Forceps
WS-696-097
10 cm / 4"
For Soft IOLs




H.R.Koch
WS-696-098
12.5 cm / 5"
Universal Prefolder
For Foldable Silicone, Hydrogel Or Acrylic IOLs




H.R.Koch
WS-696-099
12.5 cm / 5"
Prefolder
For Pharmacia CEEoN Edge



Holding Forceps
WS-697-100
10 cm / 4"
For Foldable IOLs



Strobel
WS-697-101
10 cm / 4"
Holding Forceps
For Pharmacia CEEon Edge



Lasik Corneal Flap
WS-697-102
10 cm / 4"
Forceps



Burrato
WS-697-103
10 cm / 4"
Lasik Corneal Flap
Forceps




Sattler
WS-697-104
10 cm / 4"
Muscle Forceps
Right




Sattler
WS-697-105
10 cm / 4"
Muscle Forceps
Left

Forceps | Pinzas | Pinzetten



Remky
WS-698-106
10 cm / 4"
Superior Rectus Forceps



Elschnig
WS-698-107
10 cm / 4"
Superior Rectus Forceps



Elschnig
WS-698-108
10 cm / 4"
Superior Rectus
Forceps, Delicate



Leyhecker
WS-698-109
10 cm / 4"
Superior Rectus Forceps




Castrovieio
WS-698-110
10 cm / 4"
Superior Rectus Forceps
With Tying Platform
Angled, 0.9 mm



Dardenne
WS-698-111
10 cm / 4"
Nucleus Forceps



Walser
 WS-699-112 0.5 mm
 WS-699-113 0.8 mm
 10 cm / 4"
 Suturing Forceps, Straight
 With Tying Platform



Castroviejo (New Model)
 WS-699-114 0.12 mm
 WS-699-115 0.3 mm
 WS-699-116 0.5 mm
 WS-699-117 0.9 mm
 WS-699-118 1.2 mm
 10 cm / 4"
 Suturing Forceps, Straight
 With Tying Platform



Munich Model
 WS-699-119
 10 cm / 4"
 Suturing Forceps, Straight, Conical
 With Tying Platform
 0.35 mm



Castroviejo
 WS-699-120 0.3 mm
 WS-699-121 0.6 mm
 WS-699-122 0.9 mm
 WS-699-123 1.2 mm
 10 cm / 4"
 Suturing Forceps, Straight
 With Tying Platform



Tuebingen Model
 WS-699-124
 10 cm / 4"
 Tying Forceps, Straight



Tuebingen Model
 WS-699-125
 10 cm / 4"
 Tying Forceps, Curved

Forceps | Pinzas | Pinzetten


Tuebingen Model
WS-700-126
 10 cm / 4"
 Tying Forceps
 Delicate



Bonn Model
WS-700-127
 10 cm / 4"
 Tying Forceps
 Extremely Delicate



Tying Forceps
WS-700-128
 10 cm / 4"
 Angled, 7.5 mm



Dardenne
WS-700-129
 10 cm / 4"
 Tying Forceps
 Angled, 6 mm




Kelman-Mcpherson
WS-700-130
 10 cm / 4"
 Fragment Forceps
 Angled, 7.5 mm




Gaskin
WS-700-131
 10 cm / 4"
 Fragment Forceps
 Angled, 9 mm



Bishop Harmon
 8.9 cm / 3 1/2",
WS-700-131A 0.6mm, 1 x 2 Teeth
WS-700-131B 0.4mm, 1 x 2 Teeth, Extra Delicate
WS-700-131C 0.3mm, 1 x 2 Teeth, Micro Surgery
WS-700-131D 0.6mm, Serrated
 Iris Forceps



Gaskin
WS-701-132
10 cm / 4"
Fragment Forceps
Angled, 11 mm



Barraquer
WS-701-133
6 cm / 2½"
Suture Removing Forceps



Paufigue
WS-701-134
8 cm / 3"
Tying Forceps



Paufigue
WS-701-135
8 cm / 3"
Tying Forceps
Delicate



Bangarter
WS-701-136
10 cm / 4"
Tying Forceps



WS-701-137
11.5 cm / 4½"
Tying Forceps

Forceps Vitreoretinal | Pinzas Vitreoretinal | Pinzetten Vitreoretinal


Essen Model
WS-702-001
 10 cm / 4"
 Cerclage Forceps




Holz
WS-702-002
 Double Muscle Hook
 For Encircling Bands



Holz
WS-702-003
 Cerclage Forceps



Watzke
WS-702-004
 11.5 cm / 4½"
 Spreading Forceps



Hoerauf
WS-702-005
 9 cm / 3½"
 Spreading Forceps
 For Silicon Oil Removal



Pfennigwerth
WS-702-006
 10 cm / 4"
 Holding Forceps
 For Vitrectomy Lenses



Holding Forceps
WS-702-007
 For Sclera Pins



Insertion Forceps
WS-702-008
 9 cm / 3½"
 For Sclera Pins



WS-703-009
 Vitrectomy Forceps
 Straight, Serrated Tips
 20 Gauge / 0.9 mm
 WS-703-010
 Titanium
 12.5 cm / 5"



WS-703-011
 Vitrectomy Forceps
 Straight, Rough Tips
 20 Gauge / 0.9 mm
 WS-703-012
 Titanium
 12.5 cm / 5"



WS-703-013
 Vitrectomy Forceps
 Straight, Silicone Coated Tips
 20 Gauge / 0.9 mm
 WS-703-014
 Titanium
 12.5 cm / 5"



WS-703-015
 Vitrectomy Forceps
 Straight, Grip Tips
 20 Gauge / 0.9 mm
 WS-703-016
 Titanium
 12.5 cm / 5"



WS-703-017
 Vitrectomy Forceps
 45° Angled
 20 Gauge / 0.9 mm
 WS-703-018
 Titanium
 12.5 cm / 5"

Forceps Vitreoretinal | Pinzas Vitreoretinal | Pinzetten Vitreoretinal


WS-704-019
 Vitrectomy Membrane Peeler Forceps
 45° Angled
 20 Gauge / 0.9 mm
WS-704-020
 Titanium
 12.5 cm / 5"



WS-704-021
 Vitrectomy Forceps
 Horizontal, Subretinal
 45° Angled, 20 Gauge / 0.9 mm
WS-704-022
 Titanium
 12.5 cm / 5"



WS-704-023
 Vitrectomy Membrane Peeler Forceps
 90° Angled
 20 Gauge / 0.9 mm
WS-704-024
 Titanium
 12.5 cm / 5"



WS-704-025
 Vitrectomy Forceps
 Straight, Annular Tips
 20 Gauge / 0.9 mm
WS-704-026
 Titanium
 12.5 cm / 5"



WS-704-027
 Vitrectomy Forceps
 45° Angled, Vertical Tips, Subretinal
 20 Gauge / 0.9 mm
WS-704-028
 Titanium
 12.5 cm / 5"



WS-705-029
 Vitrectomy Forceps
 With Membrane Peeler Tips
 20 Gauge / 0.9 mm
WS-705-030
 Titanium
 12.5 cm / 5"



WS-705-031
 Vitrectomy Membrane Dissector
 Straight, Serrated Tips
 20 Gauge / 0.9 mm
WS-705-032
 Titanium
 12.5 cm / 5"



WS-705-033
 Vitrectomy Foreign Body Forceps
 20 Gauge / 0.9 mm
WS-705-034
 Titanium
 12.5 cm / 5"



WS-705-035
 Vitrectomy Membrane Dissector
 Curved
 20 Gauge / 0.9 mm
 12.5 cm / 5"



WS-705-036
 Vitrectomy Membrane Dissector
 45° Angled
 20 Gauge / 0.9 mm
WS-705-037
 Titanium
 12.5 cm / 5"

Forceps Vitreoretinal | Pinzas Vitreoretinal | Pinzetten Vitreoretinal


WS-706-038 0.2 mm x 0.05 mm
WS-706-039 0.3 mm x 0.05 mm
WS-706-040 0.2 mm x 0.1 mm
WS-706-041 0.3 mm x 0.1 mm
 Vitrectomy Forceps, Straight
 End Gripping, Extremely Delicate
 20 Gauge / 0.9 mm
 12.5 cm / 5"



WS-706-042
 Vitrectomy Forceps
 90° Angled, Pointed Tips
 20 Gauge / 0.9 mm
WS-706-043
 12.5 cm / 5"



WS-706-044
 12.5 cm / 5"
 Vitrectomy Forceps, Straight
 End Gripping, Extremely Delicate
 0.2 mm x 0.05 mm
 24 Gauge / 0.6 mm



WS-706-045
 12.5 cm / 5"
 Vitrectomy Forceps, Straight
 Grip Tips
 24 Gauge / 0.6 mm



Shaaf
WS-707-001
9 cm / 3½"
Cilia & Foreign Body Forceps



Hamburg Model
WS-707-002
9 cm / 3½"
Cilia & Foreign Body Forceps



WS-707-003
11.5 cm / 4½"
Foreign Body Needle
Straight



WS-707-004
11.5 cm / 4½"
Foreign Body Needle
Curved




Plange
WS-707-005
11.5 cm / 4½"
Foreign Body Spud




WS-707-006
11.5 cm / 4½"
Foreign Body Gouge



Hamburg Model
WS-708-001
 11.5 cm / 4½"
 Keratome
 Angled



Bonn Model
WS-708-002
 11.5 cm / 4½"
 Scleral Knife
 7 mm Blade



Bonn Model
WS-708-003
 11.5 cm / 4½"
 Scleral Knife
 4 mm Blade



Phaco Slit Knife
 Angled

WS-708-004	1.50 mm	WS-708-008	2.50 mm
WS-708-005	1.70 mm	WS-708-009	2.85 mm
WS-708-006	1.80 mm	WS-708-010	3 mm
WS-708-007	2 mm	WS-708-011	3.20 mm
		WS-708-012	3.50 mm

11.5 cm / 4½"



Neuhann
WS-708-013
 11.5 cm / 4½"
 lol- Knife
 Angled, 3.5 mm



Brauweiler
WS-708-014 5 mm, For 5 mm Lens Ø
WS-708-015 5.5 mm, For 6 mm Lens Ø
WS-708-016 6.0 mm, For 6.5 mm Lens Ø
WS-708-017 6.5 mm, For 7 mm Lens Ø
 11.5 cm / 4½"
 lol- Knife, Angled



Chalazion Knife
WS-709-018
11.5 cm / 4½"
Delicate



Sato Knife
WS-709-019
11.5 cm / 4½"
Long Blade



Sato Knife
WS-709-020
11.5 cm / 4½"



Dardenne
WS-709-021
10 cm / 4"
Sato Knife



Foerster
WS-709-022
10 cm / 4"
Sato Knife



Micro Round Knife
WS-709-023
10 cm / 4"
For Visco Canalostomy
1 mm


Stiletto

For Instruments With

WS-710-024	24 Gauge / 0.6 mm
WS-710-025	22 Gauge / 0.7 mm
WS-710-026	20 Gauge / 0.9 mm
WS-710-027	19 Gauge / 1.1 mm
WS-710-028	16 Gauge / 1.6 mm



Frankfurt Model
WS-710-029
11.5 cm / 4½"
Crescent Knife




Paufrage
WS-710-030
11.5 cm / 4½"
Keratoplasty Knife



Kuhnt
WS-710-031
11.5 cm / 4½"
Corneal Scarifier
Fig. 1



Kuhnt
WS-710-032
11.5 cm / 4½"
Corneal Scarifier
Fig. 2



Kuhnt
WS-710-033
11.5 cm / 4½"
Corneal Scarifier
Fig. 3



Lohmann
WS-711-001
11.5 cm / 4½"
Lasek- Epi- Peeler



Lasek-Trephine
WS-711-002 8 mm
WS-711-003 9 mm
11.5 cm / 4½"
With Crosshairs



Lasek-Funnel
WS-711-004 8.5 mm
WS-711-005 9.5 mm
11.5 cm / 4½"
With Crosshairs



Lohmann
WS-711-006
11.5 cm / 4½"
A-Lasek Trephine
8 x 10 mm
With Crosshairs



Lohmann
WS-711-007
11.5 cm / 4½"
A-Lasek Funnel
8.9 x 10.7 mm
With Crosshairs



Weber
WS-712-008
12.5 cm / 5"
Lens Loop
Small



Weber
WS-712-009
12.5 cm / 5"
Lens Loop
Large



Snellen
WS-712-010
12.5 cm / 5"
Lens Loop




Art
WS-712-011
12.5 cm / 5"
Lens Loop
Serrated



Lucky
WS-713-012
12.5 cm / 5"
Localizer



Laqua
WS-713-013
12.5 cm / 5"
Localizer



Laqua
WS-713-014
12.5 cm / 5"
Localizer
Double Ended



Meyer / Schwickerath
WS-713-015
12.5 cm / 5"
Localizer



Kreissig
WS-713-016
12.5 cm / 5"
Localizer
Double Ended



Sclera Depressor
WS-713-017
12.5 cm / 5"
Double Ended



Sclera Depressor
WS-713-018
12.5 cm / 5"
Double Ended
Roughened Ends



Foerster
WS-713-019
12.5 cm / 5"
Inside Diameter 6.5 mm
Fixation Ring With Localizer



Thomas Neuhann
WS-714-001
11.5 cm / 4½"
Nucleus Dividing Hook



Zirm
WS-714-002
11.5 cm / 4½"
Nucleus Dividing Hook



Tobias Neuhann
WS-714-003
11.5 cm / 4½"
Chopper



Agarwal
WS-714-004
11.5 cm / 4½"
Phaco-Karate-Chopper



Haefliger
WS-714-005 For Right Hand
WS-714-006 For Left Hand
11.5 cm / 4½"
Nucleus Dividing Hook



Nagahara
WS-714-007
11.5 cm / 4½"
Nucleus Dividing Hook



Höh
WS-715-008
Chopper
With Irrigation
18 Gauge / 1.2 mm



Agarwal
WS-715-009
Irrigation Nucleus Chopper
20 Gauge / 0.9 mm



R.Koch
WS-715-010
10 cm / 4"
Irrigation Handpiece With Chopper
Angled
Irrigation Ports 2 x 0.75 mm
Tube 19 Gauge / 1.1 mm



Hoffmann
WS-715-011
10 cm / 4"
Irrigation Handpiece With Chopper
Angled
Irrigation Ports 2 x 0.5 mm
Oval Tube 1.1 x 0.9 mm



Thomas Neuhann/Sauder
WS-715-012
10 cm / 4"
Irrigation Handpiece With Chopper
Irrigation Ports 2 x 0.75 mm
Tube 19 Gauge / 1.1 mm



Tobias Neuhann
WS-716-013
11.5 cm / 4½"
Nucleus Rotator And Chopper



Thomas Neuhann / Nagahra
WS-716-014
11.5 cm / 4½"
Spatula and
Nucleus Dividing Hook



Thomas Neuhann
WS-716-015
11.5 cm / 4½"
Nucleus Rotator and
Nucleus Dividing Hook



Thomas Neuhann
WS-716-016
11.5 cm / 4½"
Iris Hook And
Nucleus Dividing Hook



Thomas Neuhann
WS-716-017
11.5 cm / 4½"
Spatula And
Nucleus Dividing Hook



Mainz Model
WS-717-018
11.5 cm / 4½"
Manipulator



Mainz Model
WS-717-019
11.5 cm / 4½"
Manipulator, New Model



Bremen Model
WS-717-020
11.5 cm / 4½"
Lens Loop Manipulator



Jacobi Manipulator
WS-717-021
WS-717-022
11.5 cm / 4½"
For 0.30 mm Lol Hole
For 0.20 mm Lol Hole



Sinsky
WS-717-023
11.5 cm / 4½"
Manipulator
For 0.20 mm Lol Hole



Lester
WS-717-024
11.5 cm / 4½"
Lens Loop Manipulator



Tandogan
WS-717-025
11.5 cm / 4½"
Explantation Hook For Lols



Frankfurt ,Model
WS-718-026
16.5 cm / 6½"
Iris Hook



Frankfurt Model
WS-718-027
16.5 cm / 6½"
Manipulator



Jacobi
WS-718-028
11.5 cm / 4½"
Iris Hook



Bonn Model
WS-718-029
11.5 cm / 4½"
Iris Hook
Extremely Delicate



Lund
WS-718-030
11.5 cm / 4½"
Iris Hook



Dardenne
WS-718-031
11.5 cm / 4½"
(Push-Pull)
Iris Hook



Neuhann
WS-719-032
10 cm / 4"
Nucleus Rotator



Agarwal
WS-719-033
10 cm / 4"
Globe Stabilizer And
Nucleus Manipulator



Liebermann
WS-719-034 Right Hand
WS-719-035 Left Hand
10 cm / 4"
Micro Finger
Nucleus Manipulator



Drysdale
WS-719-036
10 cm / 4"
Nucleus Manipulator
And Anterior Capsule Polisher



Henning-Friedrich
WS-719-037
10 cm / 4"
Lens Extraction Hook



Muenster Model
WS-719-038
15.5 cm / 6"
Suture Hook



Graefe
WS-720-039
12.5 cm / 5"
Strabismus Hook
Short


Hamburg Model
WS-720-043
12.5 cm / 5"
Strabismus Hook



Malbran
WS-720-044
12.5 cm / 5"
Strabismus Hook

Bonn Model
WS-720-045
12.5 cm / 5"
Strabismus Hook

Jameson
WS-720-048
12.5 cm / 5"
Strabismus Hook



Wecker
WS-721-001
12.5 cm / 5"
Iris Spatula
Wide



Wecker
WS-721-002
12.5 cm / 5"
Iris Spatula
Slim



Wecker
WS-721-003
12.5 cm / 5"
Iris Spatula
Conical




Frankfurt Model
WS-721-004
12.5 cm / 5"
Iris Spatula



Frankfurt Model
WS-721-005
12.5 cm / 5"
Iris Spatula
Extremely Delicate



Iris Spatula
WS-721-006
12.5 cm / 5"
Angled, Conical
Extremely Delicate



Bangerter
WS-721-007 1.00 x 0.25 mm
WS-721-008 1.50 x 0.25 mm
WS-721-009 2.00 x 0.25 mm
12.5 cm / 5"
Iris Spatula, Angled



Bangerter
WS-721-010 1.00 x 0.25 mm
WS-721-011 1.00 x 0.25 mm
WS-721-012 2.00 x 0.25 mm
12.5 cm / 5"
Iris Spatula, Round Handle, Angled



Elschnig
WS-722-013
 12.5 cm / 5"
 Spatula
 10 mm Long

Castroviejo
WS-722-014
 12.5 cm / 5"
 Spatula
 10 mm Long

Castroviejo
WS-722-015
 12.5 cm / 5"
 Spatula
 15 mm Long

Castroviejo
WS-722-016
 12.5 cm / 5"
 Spatula
 Double Ended
 10 / 15 mm Long

Bangerter
WS-722-017 6 mm
WS-722-018 10 mm
WS-722-019 12 mm
 12.5 cm / 5"
 Spatula



Castroviejo
WS-723-020
10 cm / 4"
Spatula
7 mm long



Spatula
WS-723-021
10 cm / 4"
Angled, 0.5 x 0.25 mm



Spatula
WS-723-022
10 cm / 4"
Angled, 1.0 x 0.25 mm



Hirschmann
WS-723-023
10 cm / 4"
Spatula



Steinmtz
WS-723-024
10 cm / 4"
Spatula
For Manual Lens Fragmentation
For Left-Handed Persons



Binkhorst
WS-723-025
12.5 cm / 5"
10 mm Long
Extremely Delicate



Rentsch
WS-723-026
12.5 cm / 5"
Spatula

Manipulators, Spatulas | Manipulantes, Espátulas | Manipulators, Spachtel



Essen Model
WS-724-027
11.5 cm / 4½"
Spatula



Wells
WS-724-028
11.5 cm / 4½"
Enucleation Spoon



Engels
WS-724-029
Phaco Spoon



D.D Koch
WS-724-030
9 cm / 3½"
Phaco Spatula



Amon
WS-724-031
9 cm / 3½"
Phaco Spatula



Zirm
WS-724-032
10 cm / 4"
Lasik Cornea Protector



Zirm
WS-724-033
10 cm / 4"
Lasik Cornea Surface
Cleaner



Foerster
WS-725-034
 12.5 cm / 5"
 Sato Knife
 Conical Shaft
 Front 30 Gauge / 0.3 mm
 Rear 20 Gauge / 0.9 mm



Membrane Spatula
WS-725-035
 12.5 cm / 5"
 1 mm Angled
 Long Rakes
 20 Gauge / 0.9 mm



Membrane Spatula
WS-725-036
 12.5 cm / 5"
 1 mm Angled
 Short Rakes
 20 Gauge / 0.9 mm



Membrane Spatula Knife
WS-725-037
 12.5 cm / 5"
 5 mm Angled
 All Edges Sharp
 20 Gauge / 0.9 mm



Membrane Spatula
WS-725-038
 12.5 cm / 5"
 0.5 mm Curved
 Short Rakes, Round Edge
 Rear 20 Gauge / 0.9 mm



Membrane Spatula
WS-725-039
 12.5 cm / 5"
 1.5 mm Angled
 20 Gauge / 0.9 mm



Membrane Spatula
WS-725-040
 12.5 cm / 5"
 2 and 4 mm Angled
 20 Gauge / 0.9 mm

Manipulators, Spatulas Vitreoretinal | Manipulantes, Espátulas Vitreoretinal
 Manipulators, Spachtel Vitreoretinal



Membrane Spatula
 WS-726-041
 12.5 cm / 5"
 With Pointed Tip
 5 mm Angled
 20 Gauge / 0.9 mm



Membrane Spatula
 WS-726-042
 12.5 cm / 5"
 Subretinal Pic, Flat
 4 mm Angled
 0.3 x 0.2 mm
 20 Gauge / 0.9 mm



Membrane Spatula
 WS-726-043
 12.5 cm / 5"
 Subretinal Pic, Flat
 4 mm Angled
 0.3 x 0.15 mm
 20 Gauge / 0.9 mm



Membrane Spatula
 WS-726-044
 12.5 cm / 5"
 2 mm Angled
 20 Gauge / 0.9 mm



Lucke
 WS-726-045
 12.5 cm / 5"
 Membrane Spatula
 Round Curved
 20 Gauge / 0.9 mm



Membrane Spatula
 WS-726-046
 12.5 cm / 5"
 5 mm Angled
 20 Gauge / 0.9 mm



Membrane Spatula
 WS-726-047
 12.5 cm / 5"
 2.5 mm Angled
 20 Gauge / 0.9 mm



Membrane Spatula
 WS-727-048
 12.5 cm / 5"
 Rough, Angled
 20 Gauge / 0.9 mm For MIS



Membrane Peeler
 WS-727-049
 12.5 cm / 5"
 1 mm Angled
 20 Gauge / 0.9 mm



Morris-Kroll
 WS-727-050
 12.5 cm / 5"
 Membrane Peeler
 1.5 mm Angled
 20 Gauge / 0.9 mm



Membrane Peeler
 WS-727-051
 12.5 cm / 5"
 Blunt
 With Peeling Edges
 20 Gauge / 0.9 mm



Membrane Peeler
 WS-727-052
 12.5 cm / 5"
 Rough, Angled
 20 Gauge / 0.9 mm For MIS



Membrane Peeler
 WS-727-053
 12.5 cm / 5"
 1 mm Angled
 24 Gauge / 0.6 mm



Neuhann
WS-728-001
 10 cm / 4"
 Suture Marker
 For Keratoplasty To Mark The Sutures,



Grene
WS-728-002
 10 cm / 4"
 Axis Marker
 Diameter 0.5 mm
 Concave Tip



Sinsky
WS-728-003
 10 cm / 4"
 Marking Hook
 Diameter 0.2 mm
 Round Tip



Hoffer
 Inside Diameter

WS-728-004	3 mm	WS-728-010	5 mm
WS-728-005	3.25 mm	WS-728-011	5.50 mm
WS-728-006	3.50 mm	WS-728-012	6 mm
WS-728-007	3.75 mm	WS-728-013	6.50 mm
WS-728-008	4 mm	WS-728-014	7 mm
WS-728-009	4.50 mm	WS-728-015	7.50 mm


10 cm / 4"
 Corneal Marker For Optic Zone



Thomas Neuhann
WS-728-016
 10 cm / 4"
 Corneal Marker
 4 Blades
 Inside Diameter 3 mm
 Outside Diameter 12 mm



Thomas Neuhann
WS-728-017
 10 cm / 4"
 Corneal Marker
 6 Blades
 Inside Diameter 3 mm
 Outside Diameter 12 mm



Thomas Neuhann
WS-728-018
 10 cm / 4"
 Corneal Marker
 8 Blades
 Inside Diameter 3 mm
 Outside Diameter 12 mm



Barraquer
WS-729-001
Needle Holder
Straight, Without Lock
Handle 6 mm
Total Length 10.5 cm



Barraquer
WS-729-002
Needle Holder
Straight, With Lock
Handle 6 mm
Total Length 10.5 cm



Barraquer
WS-729-003
Needle Holder
Curved, Without Lock
Handle 6 mm
Total Length 10.5 cm



Barraquer
WS-729-004
Needle Holder
Curved, With Lock
Handle 6 mm
Total Length 10.5 cm



Barraquer-Trounman
WS-729-005
Needle Holder
Straight, Without Lock
Handle 6 mm
Total Length 10 cm



Barraquer-Trounman
WS-729-006
Needle Holder
Straight, With Lock
Handle 6 mm
Total Length 10 cm

Needle Holders | Porta Agujas | Nadelhalter


Barraquer-Troutman
WS-730-007
 Needle Holder
 Curved, Without Lock
 Handle 6 mm
 Total Length 10 cm



Barraquer-Troutman
WS-730-008
 Needle Holder
 Curved, With Lock
 Handle 6 mm
 Total Length 10 cm



Castroviejo
WS-730-009
 Needle Holder
 Straight, Without Lock
 Total Length 9 cm



Castroviejo
WS-730-010
 Needle Holder
 Straight, With Lock
 Total Length 9 cm



Barraquer
WS-730-011
 Needle Holder
 Curved, Without Lock
 Jaw 0.55 mm
 Handle 8 mm
 Total Length 12 cm



Barraquer
WS-730-012
 Needle Holder
 Curved, With Lock
 Jaw 0.55 mm
 Handle 8 mm
 Total Length 12 cm



Barraquer
WS-731-013
 Needle Holder
 Curved, Without Lock
 Jaw 0.45 mm
 Handle 8 mm
 Total Length 12 cm



Barraquer
WS-731-014
 Needle Holder
 Curved, With Lock
 Jaw 0.45 mm
 Handle 8 mm
 Total Length 12 cm



Jacobi
WS-731-015
 Needle Holder
 Straight, Without Lock
 Handle 7 mm
 Total Length 11.5 cm



Jacobi
WS-731-016
 Needle Holder
 Straight, With Lock
 Handle 7 mm
 Total Length 11.5 cm



Jacobi
WS-731-017
 Needle Holder
 Curved, Without Lock
 Handle 7 mm
 Total Length 11.5 cm



Jacobi
WS-731-018
 Needle Holder
 Curved, With Lock
 Handle 7 mm
 Total Length 11.5 cm

Needle Holders | Porta Agujas | Nadelhalter


Barraquer
WS-732-019
 Needle Holder
 Curved, Without Lock
 Jaw 0.8 mm, Handle 8 mm
 Total Length 14 cm



Barraquer
WS-732-020
 Needle Holder
 Curved, With Lock
 Jaw 0.8 mm, Handle 8 mm
 Total Length 14 cm



Barraquer
WS-732-021
 Needle Holder
 Curved, Without Lock
 Jaw 0.55 mm, Handle 8 mm
 Total Length 14 cm



Barraquer
WS-732-022
 Needle Holder
 Curved, With Lock
 Jaw 0.55 mm, Handle 8 mm
 Total Length 14 cm



Barraquer
WS-732-023
 Needle Holder
 Curved, Without Lock
 Handle 10 mm
 Total Length 14 cm



Barraquer
WS-732-024
 Needle Holder
 Curved, With Lock
 Handle 10 mm
 Total Length 14 cm



Barraquer
WS-733-025
Needle Holder
Original Model
Curved, Without Lock
Handle 10 mm
Total Length 13 cm



Barraquer
WS-733-026
Needle Holder
With Strong Jaw
Curved, Without Lock
Handle 10 mm
Total Length 14 cm



Barraquer
WS-733-027
Needle Holder
With Strong Jaw
Curved, With Lock
Handle 10 mm
Total Length 14 cm



Castroviejo
WS-733-028
Needle Holder
Straight, Without Lock
Handle 10 mm
Total Length 13.5 mm



Castroviejo
WS-733-029
Needle Holder
Straight, With Lock
Handle 10 mm
Total Length 13.5 cm

Needle Holders | Porta Agujas | Nadelhalter



Castroviejo
WS-734-030
Needle Holder
Straight, Without Lock
Total Length 14 cm



Castroviejo
WS-734-031
Needle Holder
Straight, With Lock
Total Length 14 cm



Castroviejo
WS-734-032
Needle Holder
Curved, Without Lock
Total Length 14 cm



Castroviejo
WS-734-033
Needle Holder
Curved, With Lock
Total Length 14 cm



Castroviejo
WS-734-034
Needle Holder
Strong Curved, Without Lock
Total Length 14 cm



Castroviejo
WS-734-035
Needle Holder
Strong Curved, With Lock
Total Length 14 cm



Kalt
WS-735-036
13.5 cm / 5¼"
Needle Holder
Straight




Kalt
WS-735-037
13.5 cm / 5¼"
Needle Holder
Curved



Arruga
WS-735-038
13.5 cm / 5¼"
Needle Holder
Straight



Arruga
WS-735-039
13.5 cm / 5¼"
Needle Holder
Curved

Needle Holders | Porta Agujas | Nadelhalter


Boynton
WS-736-040
 Needle Holder



Combined Needle Holder / Scissors
WS-736-041 Without Lock
WS-736-042 With Lock
 12.5 cm / 5"



Axenfeld
WS-736-043
 Needle Holder



Rohrschneider
WS-736-044
 9 cm / 3½"
 Needle Holder
 Total Length 100 mm



Holding Forceps
WS-736-045
 10 cm / 4"
 For Sclera Pins




Kellnar
 WS-737-001 Right
 WS-737-002 Left
 Retrograde Probe



Hamburg Model
 WS-737-003 Right
 WS-737-004 Left
 Retrograde Probe



**Mackensen
 Trabecular Probe**
 WS-737-005 Right
 WS-737-006 Left
 Radius 7 mm



Neuhann
 WS-737-007 0.30 mm
 WS-737-008 0.25 mm
 U-Shaped Probe
 With Holding Plate
 Radius 7 mm



**Harms
 Trabecular Probe**
 WS-737-009 Right
 WS-737-010 Left
 Radius 7 mm



**Jünemann
 Lacrimal Probe**
 WS-737-011 0.60 mm, With Handle
 WS-737-012 0.80 mm, With Handle
 WS-737-013 1 mm, With Handle



**Neuhann
 WS-737-014**
 11.5 cm / 4½"
 Holder For U-Shaped Probe
 With Lock



Pigtail Probe
 WS-737-015 With Suture Holes
 WS-737-016 With Hook



Wilder
WS-738-017
10 cm / 4"
Lacrimal Probe
Delicate



Wilder
WS-738-018
10 cm / 4"
Lacrimal Probe
Medium



Wilder
WS-738-019
10 cm / 4"
Lacrimal Probe
Strong



Wilder
WS-738-020
10 cm / 4"
Lacrimal Probe
Slim, Delicate



Wilder
WS-738-021
10 cm / 4"
Lacrimal Probe
Slim, Medium



Wilder
WS-738-022
10 cm / 4"
Lacrimal Probe
Slim, Strong



Giessen Model
WS-738-023
8 cm / 3"
Lacrimal Probe



Bowman
 WS-739-024 Size 0000/000
 WS-739-025 Size 00/0
 WS-739-026 Size 1 / 2
 WS-739-027 Size 3 / 4
 WS-739-028 Size 5 / 6
 12.5 cm / 5"
 Lacrimal Probe
 Cylindrical, Silver



Bowman
 WS-739-029 Size 0000/000
 WS-739-030 Size 00/0
 WS-739-031 Size 1 / 2
 WS-739-032 Size 3 / 4
 WS-739-033 Size 5 / 6
 12.5 cm / 5"
 Lacrimal Probe
 Cylindrical, Stainless Steel



Bowman
 WS-739-034 Size 0000/000
 WS-739-035 Size 00/0
 WS-739-036 Size 1 / 2
 WS-739-037 Size 3 / 4
 WS-739-038 Size 5 / 6
 12.5 cm / 5"
 Lacrimal Probe
 Bulbous, Silver



Bowman
 WS-739-039 Size 0000/000
 WS-739-040 Size 00/0
 WS-739-041 Size 1 / 2
 WS-739-042 Size 3 / 4
 WS-739-043 Size 5 / 6
 12.5 cm / 5"
 Lacrimal Probe
 Bulbous, Stainless Steel



Mohr
WS-740-001
9 cm / 3½"
Glaucoma Punch



Jaeger
WS-740-002
10 cm / 4"
Lid Plate
Metal



Giessen Model
WS-740-003
12.5 cm / 5"
Lid Retractor



American Model
WS-740-004 Blunt
WS-740-005 Sharp
10 cm / 4"
Retractor, 4 x 4 Prongs



Müller
WS-740-006
7 cm / 2¼"
Retractor



Goldstein
WS-740-007
Retractor




Axenfeld Retractor
WS-741-008
12.5 cm / 5"
Four Pronged
Blunt




Axenfeld Retractor
WS-741-009 Blunt
WS-741-010 Sharp
12.5 cm / 5"
Three Pronged



Axenfeld Retractor
WS-741-011 Blunt
WS-741-012 Sharp
12.5 cm / 5"
Two Pronged



Rolle Retractor
WS-741-013
12.5 cm / 5"
Four Pronged
Sharp



Guthrie Retractor
WS-741-014
12.5 cm / 5"
Sharp, Large



Guthrie Retractor
WS-741-015
12.5 cm / 5"
Sharp, Small

Scissors | Tijeras | Scheren



Tuebingen Model
WS-742-001
8 cm / 3¼"
Capsulotomy Scissors
Straight




Tuebingen Model
WS-742-002
8 cm / 3¼"
Capsulotomy Scissors
Curved



Vannas
WS-742-003
8 cm / 3¼"
Capsulotomy Scissors
Extremely Delicate
Straight



Vannas
WS-742-004
8 cm / 3¼"
Capsulotomy Scissors
Extremely Delicate
Bent



Vannas
WS-742-005
8 cm / 3¼"
Capsulotomy Scissors
Extremely Delicate
Angled on Flat




Vannas
WS-742-006
8 cm / 3¼"
Capsulotomy Scissors
Extremely Delicate
Angled to Left



Capsulotomy Scissors
WS-743-007 9 mm Blade
WS-743-008 11 mm Blade
8 cm / 3 1/4"
Angled on Flat



Capsulotomy Scissors
WS-743-009
8 cm / 3 1/4"
13 mm Blade
Tip Lightly Curved



Capsulotomy Scissors
WS-743-010
8 cm / 3 1/4"
13 mm Blade
Tip Strongly Curved



Capsulotomy Scissors
WS-743-011
8 cm / 3 1/4"
10 mm Blade
Tip Strongly Curved



Capsulotomy Scissors
WS-743-012
8 cm / 3 1/4"
11 mm Blade
Tip Lightly Curved



Capsulotomy Scissors
WS-743-013
8 cm / 3 1/4"
7 mm Blade
Tip Lightly Curved

Scissors | Tijeras | Scheren



Sutter
WS-744-014
10 cm / 4"
Corneal Scissors
Light Curved
Sharp / Sharp



Castroviejo
WS-744-015
10 cm / 4"
Corneal Scissors
Blunt / Blunt
Bent




Troutman-Castroviejo
WS-744-016 Right
WS-744-017 Left
10.5 cm / 4"
Corneal Scissors, With Stop



Troutman
WS-744-018 Right
WS-744-019 Left
10 cm / 4"
Corneal Scissors
Delicate



Castroviejo
WS-744-020 Right
WS-744-021 Left
10.5 cm / 4"
Corneal Scissors



Scissors
WS-745-022
10 cm / 4"
Strong Curve
Blunt / Blunt



Scissors
WS-745-023
10 cm / 4"
Light Curve
Blunt / Blunt



Scissors
WS-745-024
10 cm / 4"
Light Curve
Sharp / Sharp



Scissors
WS-745-025
10 cm / 4"
Straight
Sharp / Sharp



Scissors
WS-745-026
11 cm / 4¼"
Straight
Sharp / Sharp


Scissors | Tijeras | Scheren



Scissors
WS-746-027
10 cm / 4"
Bent
Blunt / Blunt



Scissors
WS-746-028
10 cm / 4"
Angled On Flat
Blunt / Blunt



Castroviejo
WS-746-029
10 cm / 4"
Corneal Scissors
Strong Curve
Blunt / Blunt



Castroviejo
WS-746-030
10 cm / 4"
Corneal Scissors
With Stop, Strong Curved
Blunt / Blunt



Scissors
WS-746-031
10 cm / 4"
Strong Curve
Blunt / Blunt



Katzin
WS-747-032
10 cm / 4"
Keratoplasty Scissors
Curved Right



Katzin
WS-747-033
10 cm / 4"
Keratoplasty Scissors
Curved Left



Scissors
WS-747-034
11 cm / 4 1/4"
Straight
Sharp / Sharp



Castroviejo
WS-747-035
11.5 cm / 4 1/2"
Anterior Synechia Scissors



Castroviejo
WS-747-036
10.5 cm / 4"
Anterior Synechia Scissors
Delicate



Noyes
WS-748-037
 11 cm / 4¼"
 Iris Scissors, Straight
 Sharp / Sharp



Noyes
WS-748-038
 11 cm / 4¼"
 Iris Scissors, Straight
 Sharp / Blunt



Noyes
WS-748-039
 11 cm / 4¼"
 Iris Scissors, Straight
 Blunt / Blunt



Barraquer
WS-748-040
 6 cm / 2½"
 Iris Scissors
 7 mm Blade
 Blunt / Blunt



Wecker
WS-748-041 Sharp / Sharp
WS-748-042 Sharp / Blunt, Right
WS-748-043 Sharp / Blunt, Left
WS-748-044 Blunt / Blunt
 9 cm / 3½"
 Iris Scissors, 11 mm Blade



Wecker
WS-748-045 Sharp / Sharp
WS-748-046 Sharp / Blunt, Right
WS-748-047 Sharp / Blunt, Left
WS-748-048 Blunt / Blunt
 9 cm / 3½"
 Iris Scissors, 7 mm Blade



Trabeculum Scissors
WS-749-049
8 cm / 3"
Straight



Trabeculum Scissors
WS-749-050
8 cm / 3"
Curved



Combined
WS-749-051 Without Lock
WS-749-052 With Lock
12 cm / 4½"
Needle Holder / Scissors



Westcott Scissors
WS-749-053
Delicate
Blunt / Blunt
10 cm / 3¾"



Westcott
WS-749-054 Blunt / Blunt
WS-749-055 Sharp / Sharp
12 cm / 4½"
Scissors, Delicate



Bonn Model
WS-750-056
9 cm / 3½"
Tenotomy Scissors
Curved



Bonn Model
WS-750-057
9 cm / 3½"
Strabismus Scissors
Straight



Bonn Model
WS-750-058
9 cm / 3½"
Strabismus Scissors
Curved



Bonn Model
WS-750-059
9 cm / 3½"
Eye Scissors
Curved, Extremely Delicate




Bonn Model
WS-750-060
9 cm / 3½"
Eye Scissors
Straight, Extremely Delicate



Eye Scissors
WS-750-061
11.5 cm / 4½"
Bent
Sharp / Blunt



Eye Scissors
WS-750-062
11.5 cm / 4½"
Bent
Blunt / Blunt



Kroll
 WS-751-063
 12.5 cm / 5"
 Vitreous Body Scissors
 With Membrane Peeler
 20 Gauge / 0.9 mm
 WS-751-064
 Titanium, Angled

WS-751-065
 12.5 cm / 5"
 Vitreous Body Scissors
 For MIS
 20 Gauge / 0.9 mm
 WS-751-066
 Titanium
 Angled, Blunt / Blunt

WS-751-067
 12.5 cm / 5"
 Vitreous Body Scissors
 20 Gauge / 0.9 mm
 WS-751-068
 Titanium
 45° Angled

WS-751-069
 12.5 cm / 5"
 Vitreous Body Scissors
 With Membrane Peeler
 20 Gauge / 0.9 mm
 WS-751-070
 Titanium
 Bent

WS-751-071
 12.5 cm / 5"
 Vitreous Body Scissors
 20 Gauge / 0.9 mm
 WS-751-072
 Titanium
 Round, Curve

Scissors | Tijeras | Scheren



Heimann
WS-752-076
12.5 cm / 5"
Vitreous Body Scissors
20 Gauge / 0.9 mm
WS-752-077
Titanium
Angled, Sharp / Sharp



Heimann
WS-752-073
12.5 cm / 5"
Vitreous Body Scissors
20 Gauge / 0.9 mm
WS-752-074
Titanium
Straight, Sharp / Sharp



Heimann
WS-752-078
12.5 cm / 5"
Vitreous Body Scissors
20 Gauge / 0.9 mm
WS-752-079
Titanium
Straight, Blunt / Blunt



Heimann
WS-752-075
12.5 cm / 5"
Vitreous Body Scissors
24 Gauge / 0.6 mm
Angled, Blunt / Blunt