




WELDON
Instruments




Dental Instruments

Catalogue
Edition 2016

Extracting Forceps	1 - 52
Bone Rongeurs	53 - 58
Bone Cutting Forceps	59 - 61
Root Elevators	62 - 76
Root-Tip Picks	77
Syndesmotomes	78 - 79
Instruments Set For Oral Surgery	80
Scissors	81 - 92
Scissors With Tungsten Carbide Inserts	93 - 94
Scalpel Handles & Scalpel Blades	95 - 96
Operating & Abscess Knives	97 - 98
Reverdin Needles & Awls	99
Needle Holders	100 - 106
Needle Holders With Tungsten Carbide Inserts	107 - 113
Haemostatic Forceps	114 - 115
Towel Clamps	116
Dressing Forceps	117
Tongue Depressors	118
Mouth Gags & Tongue Forceps	119 - 121
Retractors	122 - 130
Nasal Specula	131
Laryngoscopes	132 - 133
Tissue & Suture Pliers	134 - 137
Cotton & Dressing Pliers	138 - 141
Mouth Mirrors & Handles For Mouth Mirrors	142
Instruments Set For Implant Dentistry	143
Explorers & Double-Ended Explorers	144 - 146
Periodontal Pocket Probes	147 - 148
Combination Explorer - Probes	149
Heidemann Separating Spatulas - Excavators	150 - 153
Endodontic Instruments	154 - 157
Rubber Dam Instruments & Clamps	158 - 166
Amalgam Carriers	167
Filling Instruments	168 - 184
Operative Instruments	185


Cavity Preparation Instruments	186
Excavators	187 - 190
Periodontia Instruments	191 - 212
Scalers	213 - 214
Bone Files	215 - 216
Bone Curettes	217 - 219
Crown Instruments	220 - 225
Matrix Retainers & Bands	226 - 227
Saliva Ejectors	228 - 229
Instrument Trays - Stainless Steel	230 - 232
Endodontic Boxes	233 - 234
Alcohol Burner	235
Sterilizing Forceps	236
Impression Trays	237 - 262
Articulators	263 - 266
Wax & Modeling Carvers	267 - 272
Porcelain Sculpturing Brushes	273
Ceramic Sculpturing Instruments	274
Waxing Instruments	275
Plaster Knives	276
Spatulas For Plaster & Alginate	277 - 278
Plaster Cutting Pliers & Flexible Spatula	279
Implantology	280 - 283
Mallets	284
Fret - Saw & Blades	285
Measuring Instruments	286 - 289
Crucible Tongs	290
Laboratory Tweezers	291
Pliers For Orthodontics & Prosthetics	292 - 331
Ortho Boxes	332
Orthodontic Instruments	333 - 334
Dental Syringes	335 - 340
Index	341 - 356
Instrument Care	357 - 363
Promotional Instruments	364 - 366



Most of the extracting forceps are available in Titanium & Super light models.



Super light




Titanium



Conventional



WD-001-001
Fig. 1
Upper Laterals & Canines



WD-001-002
Fig. 2
Upper Laterals & Canines




WD-001-004
Fig. 4
Lower Incisors & Canines




WD-001-003
Fig. 7
Upper Premolars



WD-001-005
Fig. 13
Lower Premolars






WD-003-012
Fig. 45
Lower Roots




WD-003-013
Fig. 33
Lower Roots




WD-003-016
Fig. 79
Lower Third Molars



WD-003-014
Fig. 51A
Upper Roots



WD-003-015
Fig. 67A
Upper Third Molars







WD-005-021
Fig. 349
Upper Roots



WD-005-022
Fig. 151 1/2
Upper Roots



WD-005-023
Fig. 346
Lower Roots




WD-005-024
Fig. 145
Lower Roots




WD-005-025
Fig. 359
Lower Roots



WD-006-026
Fig. 3
Upper Incisors & Canines
Crowded Internally, Or Externally.




WD-006-027
Fig. 8
Lower Premolars



WD-006-028
Fig. 18A
Upper Molars




WD-006-029
Fig. 19
Upper Third Molars




WD-006-030
Fig. 23
Lower Molars, Right




WD-006-031
Fig. 24
Lower Molars, Left




WD-007-032
Fig. 20
Lower Third Molars




WD-007-033
Fig. 29
Upper Roots



WD-007-034
Fig. 29N
Upper Roots & Incisors



WD-007-035
Fig. 30
Upper Roots



WD-007-036
Fig. 31
Lower Roots




WD-007-037
Fig. 33A
Lower Roots




WD-008-038
Fig. 33L
Lower Roots




WD-008-039
Fig. 38A
Lower Roots



WD-008-040
Fig. 43A
Upper Roots




WD-008-041
Fig. 41
Upper Roots



WD-008-042
Fig. 44
Upper Roots



WD-008-043
Fig. 48
Upper Roots





WD-010-050
Fig. 52
Upper Roots & Canines




WD-010-051
Fig. 59
Upper Roots



WD-010-052
Fig. 60
Lower Roots



WD-010-053
Fig. 54
For Separating Upper Molars




WD-010-054
Fig. 55
For Separating Upper Molars



WD-010-055
Fig. 56
For Separating Lower Molars



WD-011-056
Fig. 59L
Lower Incisor Roots




WD-011-057
Fig. 65L
Upper Molars Left



WD-011-058
Fig. 65R
Upper Molars Right



WD-011-059
Fig. 66R
Upper Molars Right



WD-011-060
Fig. 66L
Upper Molars Left



WD-012-061
Fig. 67
Upper Third Molars



WD-012-062
Fig. 67 1/2R
Upper Third Molars, Right



WD-012-063
Fig. 67 1/2L
Upper Third Molars, Left



WD-012-064
Fig. 73
Lower Molars



WD-012-065
Fig. 74M
Lower Roots




WD-012-066
Fig. 68
Lower Roots



WD-013-067
Fig. 74N
Lower Roots




WD-013-068
Fig. 75
Lower Canines, Premolars



WD-013-069
Fig. 79A
Lower, Third Molars



WD-013-070
Fig. 79C
Lower, Third Molars



WD-013-071
Fig. 87
Lower Molars With Decayed or
Broken-Down Crowns




WD-013-072
Fig. 99
Broken-Down Lower Molars



WD-014-073
Fig. 107
Upper Canines
Ogden-Felsch



WD-014-074
Fig. 121
Upper Third Molars



WD-014-075
Fig. 166
Trotter
Universal Pattern For Upper Theet



WD-014-076
Fig. 167
Trotter
Universal Pattern For Lower Teeth



WD-014-077
Fig. 168
Trotter
Upper Canines




WD-014-078
Fig. 169
Trotter
Lower Canines & Premolars




WD-015-079
Fig. 1
Universal Tooth Forcep
For Upper Incisors & Canines.



WD-015-080
Fig. 2
Upper Molars




WD-015-081
Fig. 3
Lower Incisors, Premolars & Roots



WD-015-082
Fig. 4
Lower Molars



WD-016-083
Fig. 1
Witzel
Lower Roots
Root Splinter Forceps




WD-016-084
Fig. 2
Witzel
Upper Roots




WD-016-085
Fig. 1
Stieglitz
Root Fragments




WD-016-086
Fig. 2
Stieglitz
Root Fragments




WD-017-087
Fig. 1
Upper Centrals & Canines




WD-017-088
Fig. 2
Upper Laterals & Canines



WD-017-089
Fig. 4
Lower Incisors & Canines



WD-017-090
Fig. 7
Upper Premolars




WD-017-091
Fig. 13
Lower Premolars

Extracting Forceps (Anatomically Shaped Handle)


Fórceps De Extracción Con Mango Adaptado | Zahnzangen Mit Paßformgriff




WD-018-092
Fig. 17
Upper Molars, Right




WD-018-093
Fig. 18
Upper Molars, Left



WD-018-094
Fig. 22
Lower Molars




WD-018-095
Fig. 21
Lower Molars




WD-019-096
Fig. 45
Lower Roots




WD-019-097
Fig. 33
Lower Roots



WD-019-100
Fig. 79
Lower Third Molars



WD-019-098
Fig. 51A
Upper Roots




WD-019-099
Fig. 67A
Upper Third Molars

Extracting Forceps (Anatomically Shaped Handle)

Fórceps De Extracción Con Mango Adaptado | Zahnzangen Mit Paßformgriff



WD-020-101
Fig. 74
Lower Roots




WD-020-103
Fig. 89
Upper Molars, Right




WD-020-104
Fig. 90
Upper Molars, Left



WD-020-102
Fig. 86C
Lower Molars



WD-021-106
Fig. 151 1/2
Upper Roots



WD-021-107
Fig. 359
Lower Roots

WD-021-105
Fig. 349
Upper Roots


WD-021-108
Fig. 145
Lower Roots

Extracting Forceps (Anatomically Shaped Handle)

Fórceps De Extracción Con Mango Adaptado | Zahnzangen Mit Paßformgriff




WD-022-109
Fig. 3
Upper Incisors & Canines



WD-022-110
Fig. 31
Upper Roots



WD-022-111
Fig. 34
Upper Incisors & Canines



WD-022-112
Fig. 36
Lower Roots Incisive Canines



WD-022-113
Fig. 36N
Lower Roots Incisive Canines



WD-023-114
Fig. 46L
Lower Roots



WD-023-115
Fig. 46L
Lower Fragments



WD-023-116
Fig. 49
Upper Roots



WD-023-117
Fig. 51
Upper Roots



WD-023-118
Fig. 51B
Upper Roots

Extracting Forceps (Anatomically Shaped Handle)

Fórceps De Extracción Con Mango Adaptado | Zahnzangen Mit Paßformgriff



WD-024-119
Fig. 51LX
Upper Roots




WD-024-120
Fig. 52




WD-024-121
Fig. 56
Separating Forceps For Molars



WD-024-122
Fig. 59
Lower Roots & Fragments




WD-024-123
Fig. 66L
Upper Molars Left



WD-024-124
Fig. 66R
Upper Molars Right



WD-025-125
 Fig. 67L
 Upper Molars Left




WD-025-126
 Fig. 67R
 Upper Molars Right



WD-025-127
 Fig. 73
 Lower Molars



WD-025-128
 Fig. 74N
 Lower Roots Incisive Canines




WD-025-129
 Fig. 79C
 Lower Third Molars



WD-025-130
 Fig. 80
 Upper Third Molars

Extracting Forceps (Anatomically Shaped Handle)


Fórceps De Extracción Con Mango Adaptado | Zahnzangen Mit Paßformgriff




WD-026-131
Fig. 87
Lower Molars




WD-026-132
Fig. 97
Upper Roots & Fragments



WD-026-133
Fig. 127R
Upper Incisors & Premolars Right



WD-026-134
Fig. 127L
Upper Incisors & Premolars Left



WD-027-135
Fig. 1
Standard
Upper Lateral & Canines

WD-027-136
Fig. 6
Separating Forceps

WD-027-137
Fig. 10S
Upper Molars

WD-027-138
Fig. 10H
Upper Molars

Extracting Forceps (American Pattern) | Fórceps De Extracción (Forma Americana)
Zahnzangen (Amerikanische Form)



WD-028-139
Fig. 16
Lower Molars



WD-028-140
Fig. 16S
Lower Molars (For Children)




WD-028-141
Fig. 15
Harris
Lower Molars



WD-028-142
Fig. 23
Lower Molars




WD-028-143
Fig. 17
Lower Molars




Harris
WD-029-144
Fig. 18R
Upper Molars, Right



Harris
WD-029-145
Fig. 18L
Upper Molars, Left




Physick
WD-029-146
Fig. 5
Lower Third Molars



Parmly
WD-029-147
Fig. 32
Upper Premolars & Molars

Extracting Forceps (American Pattern) | Fórceps De Extracción (Forma Americana)
Zahnzangen (Amerikanische Form)




WD-030-149
Fig. 53L
Upper Molars, Left



WD-030-148
Fig. 53R
Upper Molars, Right



WD-030-150
Fig. 62
Incisors, Premolars, Deciduous Teeth



WD-030-151
Fig. 65
Upper Incisors & Roots



Nevius
WD-031-153
Fig. 88R
Upper Molars, Right



Tomes
WD-031-152
Fig. 69
Upper & Lower Roots



Kells
WD-031-154
Fig. 99C
Upper Incisors, Canines, Premolars



Nevius
WD-031-155
Fig. 88L
Upper Molars, Left

Extracting Forceps (American Pattern) | Fórceps De Extracción (Forma Americana)
Zahnzangen (Amerikanische Form)



Hull
WD-032-156
Fig. 101
Deciduous Teeth



Cryer
WD-032-157
Fig. 150
Upper Incisors, Premolars, Roots



Cryer
WD-032-158
Fig. 150A
Upper Incisors, Premolars, Roots



Cryer
WD-032-159
Fig. 150S
Upper Teeth (For Children)



Cryer
WD-033-160
Fig. 151
Lower Incisors, Premolars, Roots



Cryer
WD-033-161
Fig. 151A
Lower Incisors, Premolars, Roots



Cryer
WD-033-162
Fig. 151S
Lower Teeth (For Children)

Extracting Forceps (American Pattern) | Fórceps De Extracción (Forma Americana)
Zahnzangen (Amerikanische Form)



WD-034-163
Fig. 222
Lower Third Molars



WD-034-166
Fig. 210H
Upper Third Molars



WD-034-164
Fig. 210S
Upper Third Molars



WD-034-165
Fig. 286
Upper Incisors & Roots



WD-035-174
Set of 7, Including
W-035-167
W-035-168
W-035-169
W-035-170
W-035-171
W-035-172
W-035-173



WD-035-167
13.5 cm / 5 1/4"
Fig. 3
Upper Molars



WD-035-168
14 cm / 5 1/2"
Fig. 2
Upper Premolars



WD-035-169
14.2 cm / 5 1/2"
Fig. 4
Upper Roots



WD-035-170
13.5 cm / 5 1/4"
Fig. 5
Lower Incisors



WD-035-171
13.5 cm / 5 1/4"
Fig. 1
Upper Incisors



WD-035-172
13.5 cm / 5 1/4"
Fig. 6
Lower Molars



WD-035-173
13 cm / 5"
Fig. 7
Lower Roots



WD-036-175
Fig. 37
Upper Incisors



WD-036-176
Fig. 39
Upper Molars




WD-036-177
Fig. 39R
Upper Molars, Right



WD-036-178
Fig. 39L
Upper Molars, Left




WD-036-179
Fig. 51S
Upper Roots, Right




WD-037-185
Set of 10, including
WD-036-175
WD-036-176
WD-036-177
WD-036-178
WD-036-179
WD-037-180
WD-037-181
WD-037-182
WD-037-183
WD-037-184




WD-037-180
Fig. 13S
Lower Premolars




WD-037-181
Fig. 22S
Lower Molars



WD-037-182
Fig. 30S
Upper Incisors & Roots



WD-037-183
Fig. 33S
Lower Roots



WD-037-184
Fig. 38
Lower Incisors & Canines

Extracting Forceps For Children (Anatomically Shaped Handle)


Fórceps De Extracción Con Mango Adaptado | Zahnzangen Mit Paßformgriff Für Kinder



WD-038-191
Set of 5, including
WD-038-186
WD-038-187
WD-038-188
WD-038-189
WD-038-190




WD-038-186
Fig. 22S
Lower Molars



WD-038-187
Fig. 39L
Upper Molars, Left




WD-038-188
Fig. 73S
Lower Molars




WD-038-189
Fig. 39R
Upper Molars, Right



WD-038-190
Fig. 39S
Upper Molars, Left




WD-039-192
Fig. 18A
Molars




WD-039-193
Fig. 29
Roots & Incisors




WD-039-194
Fig. 34N
Incisors & Bicuspsids,
Deep-Gripping



WD-039-195
Fig. 35N
Incisors & Bicuspsids,
Deep-Gripping




WD-039-196
Fig. 51C
Roots




WD-039-197
Fig. 67N
Third Molars, Deep-Gripping



Lawrence-Read
WD-040-198
Fig. 76
Roots




Lawrence-Read
WD-040-199
Fig. 76N
Roots




WD-040-200
Fig. 97
Very Fine Roots



WD-040-201
Fig. 151
Very Fine Roots



WD-040-202
Fig. 153
Root Fragments, Deep-Gripping



WD-041-203
Fig. 4
Incisors & Cuspids



WD-041-204
Fig. 13
Bicusps




WD-041-205
Fig. 22
Molars



WD-041-206
Fig. 23
Molar, Right



WD-041-207
Fig. 24
Molar, Left



WD-042-208
Fig. 33
Roots




WD-042-209
Fig. 33A
Roots, Slender Jaw



WD-042-210
Fig. 33C
Roots, With Anti-Slip Gripping Surface



WD-042-211
Fig. 33M
Roots, Slender Jaw



WD-043-212
Fig. 46N
Bicuspid, Deep-Gripping



WD-043-213
Fig. 74
Incisors & Roots



WD-043-214
Fig. 74XN
Incisors & Roots, Slender Jaw




WD-043-215
Fig. 75
Cuspid & Bicuspid




WD-043-216
Fig. 133
Very Fine Roots




WD-043-217
Fig. 146
Root Fragments, Deep-Gripping



WD-044-219
Fig. 79A
Third Molars



WD-044-220
Fig. 79N
Third Molars



WD-044-221
Fig. 86A
Molars With Carious
Or Broken Caps



WD-044-218
Fig. 79
Third Molars




WD-044-222
Fig. 86B
Molars With Carious Or Broken Caps



WD-044-223
Fig. 86C
Molars With Carious Or Broken Caps



WD-045-225
Fig. 233
Roots, With Parallel Beak




WD-045-226
Fig. 51
Roots



WD-045-227
Fig. 74N
Incisors & Roots, With Fine Beak



WD-045-224
Fig. 39
Molars



WD-045-228
Fig. 39A
Molars




WD-045-229
Fig. 39L
Molars, Left




WD-045-230
Fig. 39R
Molars, Right



WD-046-231
Fig. 30
Incisors & Roots




WD-046-232
Fig. 30S
Incisors & Roots




WD-046-233
Fig. 73
Molars



WD-046-234
Fig. 37
Incisors & Cuspids



Klein
WD-047-236
Fig. 137
Incisors & Cuspids



Klein
WD-047-237
Fig. 139
Bicuspid

WD-047-235
Fig. 51S
Roots

Lawrence-Read
WD-047-238
Fig. 76S
Roots



Klein
WD-048-240
Fig. 5
Incisors



Klein
WD-048-241
Fig. 6
Molars



WD-048-239
Fig. 40
Molars



Klein
WD-048-242
Fig. 7
Roots



WD-049-243
Fig. 150Sk
Teeth & Roots, Universal



WD-049-244
Fig. 17Sk
Molars



WD-049-245
Fig. 23Sk
Molars



WD-049-246
Fig. 151Sk
Teeth & Roots, Universal



WD-050-247
Fig. 150AS
Bicuspid, Incisors & Roots



Henahan
WD-050-248
Fig. 1A
Incisors




Parmly
WD-050-249
Fig. 32A
Molars, Bicuspid & Roots




WD-050-250
Fig. 62
Teeth & Roots, Universal, Also For Children



WD-050-251
Fig. 150X
Bicuspid, Incisors & Roots



Martin Ergo
WD-051-252
Fig. 74N Diam
Roots




Witzel
WD-051-253
Fig. 501
Root Fragments, Universal



Cryer
WD-051-254
Fig. 151X
Bicuspid, Incisors & Roots





WD-051-255
Fig. 151AS
Bicuspid, Incisors & Roots








WD-052-256
Fig. 222
Third Molar







Witzel
WD-052-257
Fig. 502
Root Fragments, Universal

















1/1

WD-059-001
14 cm / 5½"



1/1

WD-059-002
17 cm / 6¾"



1/1

WD-059-003
20 cm / 8"



1/1

WD-059-004
22 cm / 8¾"



1/1

WD-059-005
14 cm / 5½"



1/1

WD-059-006
17 cm / 6¾"



1/1

WD-059-007
20 cm / 8"



1/1

WD-059-008
22 cm / 8¾"



Liston




Littauer
WD-059-009
15 cm / 6"



1/1



1/1



Niro
WD-059-010
14 cm / 5½"



1/1




Cleveland
WD-060-011
17 cm / 6¾"




Cleveland
WD-060-012
15 cm / 6"




Cleveland
WD-060-013
12 cm / 4¾"



Cleveland
WD-060-014
16 cm / 6¼"




Cottle-Kazanjian
WD-061-015
19 cm / 7½"



Böhler
15 cm / 6"




Böhler
WD-061-018
15 cm / 6"



WD-061-019 18.5 cm / 7¼"
WD-061-021 23 cm / 9"
Straight
WD-061-020
WD-061-022
Curved

Ruskin-Liston



Root Elevators Are Available In Following Handles



Small Handle



Regular Handle



SB
Blasted Short Handle



SP
Special Pattern



Kopp
WD-062-012
WD-062-012-SB
WD-062-012-SP
Fig. 1

Kopp
WD-062-013
WD-062-013-SB
WD-062-013-SP
Fig. 2



Kopp
WD-062-018
WD-062-018-SB
WD-062-018-SP
Fig. 1

Kopp
WD-062-019
WD-062-019-SB
WD-062-019-SP
Fig. 2



Kopp
WD-062-014
WD-062-014-SB
WD-062-014-SP
Fig. 3

Kopp
WD-062-015
WD-062-015-SB
WD-062-015-SP
Fig. 4



Kopp
WD-062-020
WD-062-020-SB
WD-062-020-SP
Fig. 3

Kopp
WD-062-021
WD-062-021-SB
WD-062-021-SP
Fig. 4



Kopp
WD-062-016
WD-062-016-SB
WD-062-016-SP
Fig. 5

Kopp
WD-062-017
WD-062-017-SB
WD-062-017-SP
Fig. 6



Kopp



Kopp



Bein
WD-063-007
WD-063-007-SB
WD-063-007-SP
Fig. 2B



Bein
WD-063-008
WD-063-008-SB
WD-063-008-SP
Fig. 3B



Bein
WD-063-009
WD-063-009-SB
WD-063-009-SP
Fig. 4B



Bein
WD-063-010
WD-063-010-SB
WD-063-010-SP
Fig. 5B



Bein
WD-063-011
WD-063-011-SB
WD-063-011-SP
Fig. 6B




Bein
WD-063-006
WD-063-006-SB
WD-063-006-SP
Fig. 1B



Bein
WD-063-012
WD-063-012-SB
WD-063-012-SP
Fig. 1



Bein
WD-063-013
WD-063-013-SB
WD-063-013-SP
Fig. 2



Bein
WD-063-014
WD-063-014-SB
WD-063-014-SP
Fig. 3



Bein
WD-063-015
WD-063-015-SB
WD-063-015-SP
Fig. 4



Bein



Flohr
WD-064-016
WD-064-016-SB
WD-064-016-SP
Fig. 1F



Flohr
WD-064-017
WD-064-017-SB
WD-064-017-SP
Fig. 2F



Flohr
WD-064-018
WD-064-018-SB
WD-064-018-SP
Fig. 3F



Schmeckeber
WD-064-019
WD-064-019-SB
WD-064-019-SP
Fig. 4B SA
Concave Point



Schmeckeber
WD-064-020
WD-064-020-SB
WD-064-020-SP
Fig. 5B SA
Concave Point



Apexo
WD-064-021
WD-064-021-SB
WD-064-021-SP
Fig. 4 SA
Flat Point



Apexo
WD-064-022
WD-064-022-SB
WD-064-022-SP
Fig. 5 SA
Flat Point



Apical
WD-064-022A
WD-064-022A-SB
WD-064-022A-SP
Fig. 301



Apical
WD-064-022B
WD-064-022B-SB
WD-064-022B-SP
3 mm, Fig. 302





Apical
WD-064-022C
WD-064-022C-SB
WD-064-022C-SP
3 mm, Fig. 303




Apical
WD-064-022D
WD-064-022D-SB
WD-064-022D-SP
Fig. 304











Cryer
WD-068-056
WD-068-056-SB
WD-068-056-SP
Fig. 44



Cryer
WD-068-057
WD-068-057-SB
WD-068-057-SP
Fig. 45



Cryer
WD-068-058
WD-068-058-SB
WD-068-058-SP
Fig. 25



Cryer
WD-068-059
WD-068-059-SB
WD-068-059-SP
Fig. 26



Cryer
WD-068-060
WD-068-060-SB
WD-068-060-SP
Fig. 31



Cryer
WD-068-061
WD-068-061-SB
WD-068-061-SP
Fig. 32



WD-068-062
WD-068-062-SB
WD-068-062-SP
Fig. 80



WD-068-063
WD-068-063-SB
WD-068-063-SP
Fig. 81



WD-068-064
WD-068-064-SB
WD-068-064-SP
Fig. 81A



WD-068-065
WD-068-065-SB
WD-068-065-SP
Fig. 81E



Coupland
WD-069-066
WD-069-066-SB
WD-069-066-SP
Fig. A1C



Coupland
WD-069-067
WD-069-067-SB
WD-069-067-SP
Fig. A2C



Coupland
WD-069-068
WD-069-068-SB
WD-069-068-SP
Fig. A3C



WD-069-069
WD-069-069-SB
WD-069-069-SP
Fig. 190



WD-069-070
WD-069-070-SB
WD-069-070-SP
Fig. 191



Woodward
WD-069-071
WD-069-071-SB
WD-069-071-SP
Fig. 15



Woodward
WD-069-072
WD-069-072-SB
WD-069-072-SP
Fig. 16



WD-069-073
WD-069-073-SB
WD-069-073-SP
Fig. 23



WD-069-074
WD-069-074-SB
WD-069-074-SP
Fig. 24



WD-069-075
WD-069-075-SB
WD-069-075-SP
Fig. 98



Coupland
WD-070-076
WD-070-076-SB
WD-070-076-SP
Fig. 1C



Coupland
WD-070-077
WD-070-077-SB
WD-070-077-SP
Fig. 2C



Coupland
WD-070-078
WD-070-078-SB
WD-070-078-SP
Fig. 3C



Mac Millan
WD-070-079
WD-070-079-SB
WD-070-079-SP
Fig. 11M



Mac Millan
WD-070-080
WD-070-080-SB
WD-070-080-SP
Fig. 12M



Mac Millan
WD-070-081
WD-070-081-SB
WD-070-081-SP
Fig. 13M



Crane
WD-070-082
WD-070-082-SB
WD-070-082-SP
Fig. 8



University Of Michigan
WD-070-083
WD-070-083-SB
WD-070-083-SP
Fig. 21



Woodward
WD-070-084
WD-070-084-SB
WD-070-084-SP
Fig. 69W



Friedman
WD-070-085
WD-070-085-SB
WD-070-085-SP
Fig. 31F



Friedman
WD-070-086
WD-070-086-SB
WD-070-086-SP
Fig. 32F



WD-071-087
WD-071-087-SB
WD-071-087-SP
Fig. 1



WD-071-088
WD-071-088-SB
WD-071-088-SP
Fig. 2



Flohr
WD-071-089
WD-071-089-SB
WD-071-089-SP
Fig. 2



Flohr
WD-071-090
WD-071-090-SB
WD-071-090-SP
Fig. 3



Cryer
WD-071-091
WD-071-091-SB
WD-071-091-SP
Fig. 7



Cryer
WD-071-092
WD-071-092-SB
WD-071-092-SP
Fig. 8



Bein
WD-071-094
WD-071-094-SB
WD-071-094-SP
Fig. 2



Bein
WD-071-095
WD-071-095-SB
WD-071-095-SP
Fig. 3



Bein
WD-071-096
WD-071-096-SB
WD-071-096-SP
Fig. 4



Bein
WD-071-097
WD-071-097-SB
WD-071-097-SP
Fig. 5



Bein
WD-071-098
WD-071-098-SB
WD-071-098-SP
Fig. 6



Bein
WD-071-093
WD-071-093-SB
WD-071-093-SP
Fig. 1

Root Elevators (Anatomically Shaped Handle)

Elevadores Para Raíces Con Mango Adaptado | Wurzelheber Mit Paßformgriff



Bein
WD-072-099
Fig. 1B



Bein
WD-072-100
Fig. 2B



Bein
WD-072-101
Fig. 3B



Bein
WD-072-102
Fig. 4B



Bein
WD-072-103
Fig. 5B



Bein
WD-072-104
Fig. 6B



Flohr
WD-072-105
Fig. 1F



Flohr
WD-072-106
Fig. 2F



Flohr
WD-072-107
Fig. 3F



Schmeckebier-Apexo
WD-072-108
Fig. 4SA




Schmeckebier-Apexo
WD-072-109
Fig. 5SA



Miller-Apexo
WD-073-110
Fig. 71



Miller-Apexo
WD-073-111
Fig. 72



Miller-Apexo
WD-073-112
Fig. 73



Miller-Apexo
WD-073-113
Fig. 74



WD-073-114
Fig. 19



WD-073-115
Fig. 20



Seldin
WD-073-116
Fig. 1L



Seldin
WD-073-117
Fig. 1R



Seldin
WD-073-118
Fig. 4L



Seldin
WD-073-119
Fig. 4R



Seldin
WD-073-120
Fig. 301W



Seldin
WD-073-121
Fig. 304W

Root Elevators (Anatomically Shaped Handle)

Elevadores Para Raíces Con Mango Adaptado | Wurzelheber Mit Paßformgriff



Apexo
WD-074-123
Fig. 301



Apexo
WD-074-123
Fig. 302



Apexo
WD-074-124
Fig. 303



Heidbrink
WD-074-125
Fig. 1H



Heidbrink
WD-074-126
Fig. 2H



Heidbrink
WD-074-127
Fig. 40H



Heidbrink
WD-074-128
Fig. 41H



Cryer
WD-074-129
Fig. 44



Cryer
WD-074-130
Fig. 45



Cryer
WD-074-131
Fig. 25



Cryer
WD-074-132
Fig. 26



Cryer
WD-074-133
Fig. 31



Cryer
WD-074-134
Fig. 32



Walter F. Barry
WD-075-135
Stainless Steel Handle



Walter F. Barry
WD-075-136 Fig. 320



Walter F. Barry
WD-075-137 Fig. 321



Walter F. Barry
WD-075-138 Fig. 322



Walter F. Barry
WD-075-139 Fig. 323



Walter F. Barry
WD-075-140 Fig. 350



Walter F. Barry
WD-075-141 Fig. 351



Walter F. Barry
WD-075-142 Fig. 370



Walter F. Barry
WD-075-143 Fig. 371



Lecluse
WD-075-144
Stainless Steel Handle



Lecluse
WD-075-145
Fig. 1LC



Winter
WD-076-146
Stainless Steel Handle



Winter
WD-076-147
Fig. 11LW



Winter
WD-076-148
Fig. 11RW



Winter
WD-076-149
Fig. 11LW



Winter
WD-076-150
Fig. 11RW



Winter
WD-076-151
Fig. 12LW




Winter
WD-076-152
Fig. 12RW



Winter
WD-076-153
Fig. 14LW



Winter
WD-076-154
Fig. 14RW





WD-078-002
Fig. 0

WD-078-003
Fig. 1




WD-078-004
Fig. 2

WD-078-005
Fig. 3



WD-078-006
Fig. 4

WD-078-007
Fig. 5



WD-078-001
11 cm / 4 1/4"
Handle Only



WD-078-008
Fig. 6



WD-078-009
Fig. 7



WD-078-010
Fig. 8



WD-078-011
Fig. 9



WD-078-012
Fig. 11


WD-078-013
Complete Set
Consisting Of 1 Handle & 11 Interchangeable Tips
Fig.0 - Fig.11

WD-078-014
Complete Set
Consisting Of 1 Handle & 20 Tips
Fig.0 - Fig.11 & Fig.21 - Fig.19




WD-078-014
Complete Set
Consisting Of 1 Handle & 20 Tips
Fig.0 - Fig.11 & Fig.21 - Fig.19





WD-080-001
Contents


Curette Gracey x 3
 Scissor Goldman-Fox x 2
 Bone Curette Lucas x 2
 Knife Kirkland
 Knife Orban
 Periosteals Molt
 Periosteals Prichard
 Scalpel Handles Straight
 Scalpel Handles Curved
 Chisel Ochsenbein x 2
 Needle Holder - Castroviejo
 Needle Holder - Mathieu
 Tray


WD-080-002
Contents


Handle For Mouth Mirror
 Mouth Mirror
 Explorer
 Cutting And Dressing Pliers - London College
 Tissue Plier - Semken - Taylor
 Knives Goldman - Fox
 Knives Orban
 Periosteals Molt
 Periosteals Prichard
 Retractor Farabeuf
 Retractor Sternberg
 Retractor Senn-Miller
 Needle Holder Crile-wood
 Needle Holder Mathieu
 Scissor Goldman - Fox
 Haemostatic Forceps - Halstead Mosquito, Straight
 Haemostatic Forceps - Halstead Mosquito, Curved
 Tray


WD-080-003
Contents

Cotton And Dressing Plier
 Tissue Plier
 Suture Plier
 Retractor
 Needle Holder Mathieu
 Needle Holder Castroviejo
 Scalpel Handle
 Surgical Scissor - Spencer
 Surgical Scissor Curved - Metzenbaum
 Handle For Mouth Mirror
 Mouth Mirror
 Retractor - Buchwald
 Tray




WD-081-001
10.5 cm / 4 1/4"




WD-081-002
11.5 cm / 4 1/2"




WD-081-003
13 cm / 5"




WD-081-004
14 cm / 5 1/2"



WD-081-005
15.5 cm / 6"




WD-081-006
16.5 cm / 6 1/2"



WD-081-007
17.5 cm / 7"




Operating
Sharp / Blunt
Straight



1/1

WD-082-008
10.5 cm / 4¼"



1/1

WD-082-009
11.5 cm / 4½"




1/1

WD-082-010
13 cm / 5"




Operating
Sharp / Sharp
Straight



1/1

WD-082-011
14 cm / 5½"



1/1

WD-082-012
15.5 cm / 6"




1/1


WD-082-013
17.5 cm / 7"




WD-083-014
11.5 cm / 4½"



WD-083-015
13 cm / 5"




WD-083-016
14 cm / 5½"



WD-083-017
15.5 cm / 6"




WD-083-018
16.5 cm / 6½"



WD-083-019
17.5 cm / 7"




Operating
Sharp / Blunt
Curved




1/1

WD-084-020
11.5 cm / 4½"



1/1

WD-084-021
13 cm / 5"




1/1

WD-084-022
14 cm / 5½"




Operating
Sharp / Sharp
Curved



1/1

WD-084-023
15.5 cm / 6"






1/1

WD-084-024
16.5 cm / 6½"




WD-084-025
17.5 cm / 7"






Straight

Mayo
WD-086-036 14.5 cm / 5¾"
WD-086-037 17 cm / 6¾"




Curved

Mayo
WD-086-038 14.5 cm / 5¾"
WD-086-039 17 cm / 6¾"




WD-086-040
Straight


WD-086-041
Curved



WD-086-042
Blunt, Straight




WD-086-043
Blunt, Curved




WD-086-044
Straight

WD-086-045
Curved




Stevens
11.5 cm / 4½"
Dissecting & Strabismus



Stevens
11.5 cm / 4½"
Dissecting & Strabismus



Metzenbaum
WD-087-046 14.5 cm / 5 3/4"
WD-087-048 16 cm / 6 1/4"
WD-087-050 18 cm / 7 1/8"
Straight



Metzenbaum
WD-087-047 14.5 cm / 5 3/4"
WD-087-049 16 cm / 6 1/4"
WD-087-051 18 cm / 7 1/8"
Curved



Metzenbaum-Delicate
WD-087-052
14.5 cm / 5 3/4"
Straight



Metzenbaum-Delicate
WD-087-053
14.5 cm / 5 3/4"
Curved



Lister
WD-087-054
18 cm / 7 1/8"
Bandage Scissors



Lister
WD-087-055
14.5 cm / 5 3/4"
Bandage Scissors




Universal
WD-087-056 15 cm / 6"
WD-087-057 17 cm / 6 3/4"
Autoclavable Up To 143° C



Universal Scissors
are available in
above colors




Barraquer
WD-088-058
14 cm / 5 1/2"
Surgical Scissors
Straight



Barraquer
WD-088-059
Surgical Scissors
14 cm / 5 1/2"
Curved



Castro
WD-088-060
Surgical Scissors
45°



Fine Operating Scissors
12 cm / 4¾"

WD-089-061
Fig. 1
Straight



WD-089-062
Fig. 2
Curved



Goldman-Fox
WD-089-063
13 cm / 5"
Straight, Saw Edge




Goldman-Fox
WD-089-064
13 cm / 5"
Curved, Saw Edge



Iris
11.5 cm / 4½"

WD-089-065
Straight



WD-089-066
Curved




WD-089-067
Side Curved



Fine Operating Scissors
11.5 cm / 4½"



WD-089-068
Angular



WD-089-069
Angular, Saw Edge



Dean
WD-090-070
17 cm / 6 3/4"




Neumann
WD-090-071
12.5 cm / 5"




La Grange
WD-090-072
Curved, Saw Edge
10.5 cm / 4 1/4"



Goldman-Fox (Special)
WD-090-073
13 cm / 5"
Curved, Saw Edge




Quinby
WD-090-074
12.5 cm / 5"



Quinby
WD-090-075
12.5 cm / 5"



Kelly
WD-091-076
16 cm / 6¼"
Straight



Kelly
WD-091-077
16 cm / 6¼"
Curved



Locklin
WD-091-078
16.5 cm / 6½"
Saw Edge



Locklin
WD-091-079
16.5 cm / 6½"
Saw Edge



Sullivan
WD-091-080
14.5 cm / 5¾"
Straight, Saw Edge



Sullivan
WD-091-081
14.5 cm / 5¾"
Curved, Saw Edge



Noyes - Castroviejo
WD-092-082
12 cm / 4¾"
Very Fine Points



Westcott
WD-092-083
11 cm / 4¼"
Sharp Pointed



Westcott
WD-092-084
11 cm / 4¼"
Blunt, Pointed



WD-092-085
14.5 cm / 5¾"
Very Fine Points



Goldman-Fox
WD-093-086
13 cm / 5"
Straight, One Blade Serrated
TC



Goldman-Fox
WD-093-087
13 cm / 5"
Curved, One Blade Serrated
TC



Iris
WD-093-088
11 cm / 4 1/4"
Straight
TC



Iris
WD-093-089
11 cm / 4 1/4"
Curved
TC

Gum Scissors With Tungsten Carbide Inserts | Tijeras Para Encías Con Carburo De Tungsteno
Zahnfleischscheren Mit Hartmetalleinlagen



Dean
WD-094-090
17 cm / 6 3/4"
TC



Locklin
WD-094-091
16.5 cm / 6 1/2"
TC




Locklin
WD-094-092
16.5 cm / 6 1/2"
TC





Kelly
WD-094-093
16 cm / 6 1/4"
Straight
TC



Kelly
WD-094-094
16 cm / 6 1/4"
Curved
TC








Dieffenbach




Ash
WD-098-033
Fig. 1



Ash
WD-098-034
Fig. 2



Schroeder
WD-098-035
22 cm / 8 3/4"



WD-098-036
13 cm / 5 1/8"



Reverdin



Awl For Lower Jaw



WD-099-001
Fig. 1
14 cm / 5½"




WD-099-002
Fig. 2
14 cm / 5½"




WD-099-003
Fig. 3
14 cm / 5½"




WD-099-004
Fig. 4
14 cm / 5½"




WD-099-005
Fig. 1
19.5 cm / 7¾"



WD-099-006
Fig. 2
19.5 cm / 7¾"



WD-099-007
Fig. 3
19.5 cm / 7¾"



WD-099-008
Fig. 4
19.5 cm / 7¾"



Awl For Upper Jaw



WD-099-010
18 cm / 7"



Kilner
WD-099-011
Skin Hook & Retractors
16 cm / 6¼"

WD-099-009
13 cm / 5"



Derf - Wright
WD-100-001
12 cm / 4¾"



Derf
WD-100-002
12 cm / 4¾"



Collier
WD-100-003
12.5 cm / 5"



Hegar-Baumgartner
WD-100-004
14 cm / 5½"



Halsey-Smooth
WD-100-005
13 cm / 5"



Halsey
WD-100-006
13 cm / 5"



Baumgartner
WD-100-007
14.5 cm / 5¾"



Mini-Ryder
WD-101-008
18 cm / 7 $\frac{1}{8}$ "



Mini-Ryder
WD-101-009
15 cm / 6"



Fine Swedish
WD-101-010
15 cm / 6"



Micro-Vascular
WD-101-011
15 cm / 6"



Crile-Wood
WD-102-012
15 cm / 6"




Adson
WD-102-013
17.5 cm / 7"



Mayo-Hegar
WD-102-014 14 cm / 5½"
WD-102-015 15 cm / 6"
WD-102-016 16 cm / 6¼"



Mayo-Hegar
WD-102-017
18 cm / 7½"




Gillies
WD-103-018 14.5 cm / 5¾"
WD-103-019 16 cm / 6¼"



Olsen-Hegar
WD-103-020 14 cm / 5½"
WD-103-021 17 cm / 6¾"
Needle Holder & Scissors Combined



Toennis
WD-103-022
18 cm / 7⅞"
Light Pattern



Eiselberg-Mathieu
WD-103-023
19.5 cm / 7¾"



Mathieu
WD-104-024
17 cm / 6 3/4"



Mathieu
WD-104-025
17 cm / 6 3/4"



Mathieu
WD-104-026
14 cm / 5 1/2"




Mathieu
WD-104-027
14 cm / 5 1/2"



Boynton
WD-105-028
12 cm / 4 3/4"



Lichtenberg
WD-105-029
17 cm / 6 3/4"



Arruga
WD-105-030
17 cm / 6 3/4"



Castroviejo
WD-106-031
14 cm / 5½"
Serrated, Straight Jaws



Castroviejo
WD-106-032
14 cm / 5½"
Serrated, Curved Jaws



Castroviejo
WD-106-033
13 cm / 5"
Smooth Jaws



Castroviejo
WD-106-034
13 cm / 5"
Serrated Jaws



Castroviejo-Barraquer
WD-106-035
12 cm / 4¾"
Smooth Jaws



Derf
WD-107-001
12 cm / 4¾"
TC



Collier
WD-107-002
12.5 cm / 5"
TC



Hegar-Baumgartner
WD-107-003
14 cm / 5½"
TC



Halsey-Smooth
WD-107-004
13 cm / 5"
TC



Halsey
WD-107-005
13 cm / 5"
TC

Needle Holders With Tungsten Carbide Inserts | Porta Agujas Con Carburo De Tungsteno
Nadelhalter Mit Hartmetalleinlagen





Needle Holders With Tungsten Carbide Inserts | Porta Agujas Con Carburo De Tungsteno
Nadelhalter Mit Hartmetalleinlagen



Gillies
WD-110-016 14.5 cm / 5¾"
WD-110-017 16 cm / 6¼"
TC



Baumgartner
WD-110-018
14.5 cm / 5¾"
TC



Toennis
WD-110-019
18 cm / 7½"
Light Pattern
TC



Mathieu
WD-111-020
17 cm / 6 3/4"
TC



Mathieu
WD-111-021
17 cm / 6 3/4"
TC



Mathieu
WD-111-022
14 cm / 5 1/2"
TC



Mathieu
WD-111-023
14 cm / 5 1/2"
TC


Needle Holders With Tungsten Carbide Inserts | Porta Agujas Con Carburo De Tungsteno
Nadelhalter Mit Hartmetalleinlagen



Boynton
WD-112-024
12 cm / 4 3/4"
TC



Lichtenberg
WD-112-025
17 cm / 6 3/4"
TC



Arruga
WD-112-026
16 cm / 6 1/4"
TC




Castroviejo
WD-113-027
14 cm / 5½"
TC




Castroviejo
WD-113-028
14 cm / 5½"
TC




Castroviejo
WD-113-029
13 cm / 5"
TC




Hartmann
 WD-114-001 10 cm / 4" WD-114-003
 WD-114-002 12 cm / 4¾" WD-114-004
 Straight Curved




Halsted-Mosquito
 WD-114-005 10 cm / 4" WD-114-007
 WD-114-006 12 cm / 4¾" WD-114-008
 Straight Curved




Halsted-Mosquito
 WD-114-009 10 cm / 4" WD-114-011
 WD-114-010 12 cm / 4¾" WD-114-012
 Straight Curved
 1:2 Teeth




Kelly
 WD-114-013 12.5 cm / 5" WD-114-015
 WD-114-014 14 cm / 5½" WD-114-016
 Straight Curved



Crile
WD-115-017 14 cm / 5½" WD-115-019
WD-115-018 15 cm / 6" WD-115-020
Straight Curved




Pean
WD-115-021 14 cm / 5½" WD-115-023
WD-115-022 16 cm / 6¼" WD-115-024
Straight Curved







Kocher
WD-115-025 14 cm / 5½" WD-115-027
WD-115-026 16 cm / 6¼" WD-115-028
Straight Curved
1:2 Teeth



Spencer-Wells
WD-115-029 13 cm / 5" WD-115-032
WD-115-030 14 cm / 5½" WD-115-033
WD-115-031 16 cm / 6¼" WD-115-034
Straight Curved





Lip Retractors


Columbia
WD-118-001
14 cm / 5½"



Sternberg
WD-118-002
14 cm / 5½"



Sternberg
WD-118-003
16 cm / 6¼"




Sternberg
WD-118-004
14 cm / 5½"

Cheek Retractors


Cawood Minnesota
WD-118-005
15 cm / 6"
20 mm



WD-118-006
14.5 cm / 5¾"



WD-118-007
14.5 cm / 5¾"



Buchwald
WD-118-008
17 cm / 6¾"



Brünings
WD-118-009
19 cm / 7½"



Collin
WD-119-010
16 cm / 6¼"



Young
WD-119-011
17 cm / 6¾"




Heister
WD-119-012
12 cm / 4¾"



Roser-Koenig
WD-119-013
19 cm / 7½"



Roser-Koenig
WD-119-014
16 cm / 6¼"



Denhart
WD-119-015
12 cm / 4¾"



Ferguson-Ackland
WD-120-016
14.5 cm / 5¾"



Molt
WD-120-017
14 cm / 5½"
For Adults




Molt
WD-120-018
11 cm / 4¼"
For Children




Doyen-Collin
WD-120-019
12 cm / 4¾"




Doyen-Jansen
WD-120-020
14 cm / 5½"




Whitehead
WD-121-021
13 cm / 5"



Whitehead
WD-121-022 11 cm / 4¼"
WD-121-023 13 cm / 5"
WD-121-024 15 cm / 6"



Jennings
WD-121-025 10.5 cm / 4"
WD-121-026 12.5 cm / 5"
WD-121-027 14.5 cm / 5¾"
WD-121-028 16.5 cm / 6½"




Senn-Miller
WD-122-001
16 cm / 6¼"



Senn-Miller
WD-122-002
16 cm / 6¼"



Ragnell
WD-122-003
16 cm / 6¼"




Straight
WD-122-004
13 cm / 5"




Barsky
WD-122-005
15 cm / 6"
Rigid



Barsky
WD-122-006
15 cm / 6"
Flexible




Gillies
19 cm / 7½"



WD-122-007
Fig. 1

WD-122-008
Fig. 2



Guthrie
WD-123-009
13 cm / 5"



Guthrie
WD-123-010
13 cm / 5"



Knapp
WD-123-011
13 cm / 5"



Knapp
WD-123-012
16 cm / 6 1/4"



Knapp
WD-123-013
16 cm / 6 1/4"
Flexible



Knapp
WD-123-014
16 cm / 6 1/4"
Flexible




Machatt
WD-123-015
16 cm / 6 1/4"




Machatt
WD-123-016
16 cm / 6 1/4"









Williger
13 cm / 5"
Straight, Sharp




Williger
13 cm / 5"
Left, Sharp




Williger
13 cm / 5"
Right, Sharp




Williger
13 cm / 5"
Straight, Blunt



Williger
13 cm / 5"
Left, Blunt



Williger
13 cm / 5"
Right, Blunt



WD-127-048
40 x 11 mm

WD-127-049
55 x 11 mm



WD-127-050
40 x 11 mm


WD-127-051
55 x 11 mm



Ramus
WD-127-052
72 x 11 mm
21 cm / 8 1/4"







Ramus
WD-127-053
35 x 11 mm
21 cm / 8 1/4"


Middeldorpf

WD-128-054
21.5 cm / 8½"
14 x 17 mm




WD-128-055
22 cm / 8½"
20 x 22 mm

WD-128-056
23.5 cm / 9¼"
26 x 30 mm

WD-128-057
20 cm / 8"

Langenbeck
21 cm / 8¼"

WD-128-058
30 x 11 mm

WD-128-059
30 x 14 mm

WD-128-060
30 x 16 mm

WD-128-061
40 x 11 mm

WD-128-062
50 x 11 mm


Langenbeck
21 cm / 8¼"

WD-128-063
30 x 11 mm

WD-128-064
30 x 14 mm

WD-128-065
30 x 16 mm

WD-128-066
40 x 11 mm

WD-128-067
50 x 11 mm




WD-129-068
Sharp




WD-129-069
Blunt



Weitlaner
10.5 cm / 4 1/4"
2 : 3 Teeth



WD-129-070
Sharp



WD-129-071
Blunt




Weitlaner-Wullstein
13 cm / 5"
3 : 4 Teeth




WD-129-072
13 cm / 5"
Sharp



WD-129-073
16.5 cm / 6 1/2"
Sharp



WD-129-074
13 cm / 5"
Blunt



WD-129-075
16.5 cm / 6 1/2"
Blunt



Weitlaner
3 : 4 Teeth




Weitlaner
WD-130-076
13 cm / 5"
Sharp



Weitlaner
WD-130-077
13 cm / 5"
Blunt



Alm
WD-130-078
7 cm / 2 3/4"
Sharp



Alm
WD-130-079
7 cm / 2 3/4"
Blunt



Alm
WD-130-080
10 cm / 4"
Sharp



Alm
WD-130-081
10 cm / 4"
Blunt



Langenbeck
WD-131-001
22 mm Ø
Finger Protectors



Langenbeck
WD-131-002
24.5 mm Ø
Finger Protectors



Langenbeck
WD-131-003
26.5 mm Ø
Finger Protectors



WD-131-004



WD-131-005



WD-131-006



Hartmann-Halle
15 cm / 6"



Antrum Probes

Bowman
WD-131-007
Fig. 1/2, 0.9/1.1 mm Ø
WD-131-008
Fig. 7/8, 1.8/1.9 mm Ø
13 cm / 5"



WD-131-009
Fig. 1, 1 mm Ø
WD-131-010
Fig. 2, 2 mm Ø
14 cm / 5 1/2"


Cotton Applicators




Farrell
(Lathbury-Uebe)
WD-131-011
Fig. 1, 0.9 mm Ø
12 cm / 4 3/4"
WD-131-012
Fig. 2, 0.9 mm Ø
14 cm / 5 1/2"
WD-131-013
Fig. 3, 1.1 mm Ø
16.5 cm / 6 1/2"




WD-132-004
110 mm
Fig. 4



WD-132-003
90 mm
Fig. 3



WD-132-002
75 mm
Fig. 2



WD-132-001
55 mm
Fig. 1



WD-132-005
Conventional Handle



WD-132-006
Spare Bulbs



Mc Intosh
WD-132-007
Laryngoscopes In Case
Consists of
WD-132-001 x 1
WD-132-003 x 1
WD-132-004 x 1
WD-132-005 x 1
WD-132-006 x 3



McIntosh Mega Blade Fiber Optic		
	Fig.	Working Length (mm)
WD-133-008	0	55
WD-133-009	1	70
WD-133-010	2	78
WD-133-011	3	90
WD-133-012	4	108
WD-133-013	5	133
WD-133-014	6	150



McIntosh Fiber Optic - (Cold Light Guide Fixed) Re-Useable		
	Fig.	Working Length (mm)
WD-133-015	0	55
WD-133-016	1	70
WD-133-017	2	78
WD-133-018	3	90
WD-133-019	4	108
WD-133-020	5	133
WD-133-021	6	150



WD-133-122
Fiber Optic Handle

Tissue Pliers - Suture Pliers | Pinzas De Disección - Pinzas Para Suture
Chirurgische Pinzetten - Nahtpinzette



Suture Pliers
WD-134-001
15 cm / 6"



WD-134-002
Very Delicate
10 cm / 4"
1:2 Teeth



Semken-Taylor
WD-134-003
12.5 cm / 5"



Semken-Taylor
WD-134-004
12.5 cm / 5"



Semken-Taylor
WD-134-005
12.5 cm / 5"
1:2 Teeth



Semken-Taylor
WD-134-006
12.5 cm / 5"
1:2 Teeth



Adson
WD-135-007 12 cm / 4 3/4"
WD-135-008 15 cm / 6"



Adson
WD-135-009 12 cm / 4 3/4"
WD-135-010 15 cm / 6"
1:2 Teeth




Micro-Adson
WD-135-011
12 cm / 4 3/4"




Micro-Adson
WD-135-012
15 cm / 6"



Micro-Adson
WD-135-013
12 cm / 4 3/4"
1:2 Teeth




Micro-Adson
WD-135-014
15 cm / 6"
1:2 Teeth



Adson-Brown
Non-Traumatic
WD-135-015 12 cm / 4 3/4"
WD-135-016 15 cm / 6"



WD-136-017
11.5 cm / 4½"
Fig. 4



WD-136-018
11.5 cm / 4½"
Fig. 5




WD-136-019
12.5 cm / 5"
Fig. 1




WD-136-020
12.5 cm / 5"
Fig. 2



WD-136-021
11.5 cm / 4½"
Fig. 3



Dressing Forceps
Narrow



Tissue Forceps
1:2 Teeth



Potts-Smith



Goldman-Fox
(Crane-Kaplan)
WD-138-004
15 cm / 6"
Fig. 1



Goldman-Fox
(Crane-Kaplan)
WD-138-005
15 cm / 6"
Fig. 2



Locking
WD-138-006
16 cm / 6 1/4"
Fig. 18L
Self-Locking



WD-138-007
15 cm / 6"
Fig. 1
Self-Locking



WD-138-008
15 cm / 6"
Fig. 2
Self-Locking



London-College
15 cm / 6"



WD-139-009
Serrated



WD-139-010
Serrated



WD-139-011
Serrated



WD-139-012
Smooth



WD-139-013
Smooth



WD-139-014
Smooth



London-College
15 cm / 6"
With lock



WD-139-015
NO. 5



WD-139-016
NO. 8
For Silver Points



Perry
WD-139-017
13 cm / 5"
Serrated



Perry
WD-139-018
13 cm / 5"
Smooth



Meriam
16 cm / 6 1/4"



WD-140-019
NO. 1



WD-140-020
NO. 3



WD-140-021
NO. 5



Flagg
16 cm / 6 1/4"



WD-140-022
NO. 1



WD-140-023
NO. 2




WD-140-024
NO. 3



WD-140-025
NO. 4



College
WD-141-026
16 cm / 6¼"
Fig. 2




College
WD-141-027
16 cm / 6¼"
Fig. 3



College
WD-141-028
16 cm / 6¼"
Fig. 17



College
WD-141-029
16 cm / 6¼"
Fig. 18



College
WD-141-030
16.5 cm / 6½"
Fig. 19



Magnifying



WD-142-001
Plane

WD-142-002
Magnifying



WD-142-003
Plane

WD-142-004
Magnifying

WD-142-005
Rodium




WD-142-006
Plane

WD-142-007
Magnifying

WD-142-008
Rodium




Plane



Cone Socket
WD-142-009
Stainless Steel



Cone Socket
WD-142-010
Stainless Steel



Cone Socket
WD-142-011
Stainless Steel




WD-142-012
Patient Bib Chains (Plastic)




WD-142-013
Parallelometer Mouth Mirror With Handle

WD-142-014
Parallelometer Mouth Mirror




Conservative Treatment Instruments Set
WD-143-001

- Consists of :
- Periosteals Molt
 - Periosteals Prichard
 - Periosteals Goldman-Fox
 - Knife Kirkland
 - Mallets-Mead
 - Needle Holder Castroviejo
 - Needle Holder Mathieu
 - Scalpel Handle
 - Instrument Tray



Conservative Treatment Instruments Set
WD-143-002

- Consists of :
- Scalpel Handle
 - Scalpel Blades
 - Retractor Sternberg
 - Scissor Goldman-Fox
 - Periosteals Molt
 - Periosteals Prichard
 - Needle Holder Crile Wood
 - Tissue Plier Semken Taylor
 - Explorer
 - Handle For Mouth Mirror
 - Mouth Mirror
 - Instrument Tray



Conservative Treatment Instruments Set
WD-143-003

- Consists of :
- Instrument Tray
 - Filling Instrument
 - Placement Instruments
 - Amalgam Carver
 - Excavator
 - Mouth Mirror
 - Handle For Mouth Mirror
 - Cotton & Dressing Plier
 - Cavity Preparation Instruments
 - Explorer-Probes
 - Cement Spatula
 - Heidemann Spatula
 - Gum Scissor



HH
Hollow Handle



HH
Hollow Handle



HH
Hollow Handle



SD
Solid Handle



SD
Solid Handle



SD
Solid Handle



WD-145-001-SD
WD-145-001-HH
Fig. 16



WD-145-002-SD
WD-145-002-HH
Fig. 9



WD-145-003-SD
WD-145-003-HH
Fig. 10



WD-145-004-SD
WD-145-004-HH
Fig. 12



WD-145-005-SD
WD-145-005-HH
Fig. 2



WD-145-006-SD
WD-145-006-HH
Fig. 11



WD-145-007-SD
WD-145-007-HH
Fig. 6



WD-145-008-SD
WD-145-008-HH
Fig. 8



WD-145-009-SD
WD-145-009-HH
Fig. 5



WD-145-014-SD
WD-145-014-HH
Fig. 23



WD-145-010-SD
WD-145-010-HH
Fig. 17



WD-145-011-SD
WD-145-011-HH
Fig. 17A



WD-145-012-SD
WD-145-012-HH
Fig. 3



WD-145-013-SD
WD-145-013-HH
Fig. 3A



WD-146-016-SD
WD-146-016-HH
Fig. 23/17


WD-146-017-SD
WD-146-017-HH
Fig. 23/17A

WD-146-018-SD
WD-146-018-HH
Fig. 2

WD-146-019-SD
WD-146-019-HH
Fig. 2A

WD-146-020-SD
WD-146-020-HH
Fig. 3

Root Canal Explorers



Cattoni
WD-146-015-SD
WD-146-015-HH
Fig. UT.1/2


WD-146-021-SD
WD-146-021-HH
Fig. 16/23

WD-146-022-SD
WD-146-022-HH
Fig. 16/3

WD-146-023-SD
WD-146-023-HH
Fig. 6/1

WD-146-024-SD
WD-146-024-HH
Fig. 6/17

WD-146-025-SD
WD-146-025-HH
Fig. 6/23



Cattoni
WD-147-026-SD
WD-147-026-HH
Fig. UT.1/2


Goldman-Fox / Williams
WD-147-027-SD
WD-147-027-HH

Nabers
WD-147-028-SD
WD-147-028-HH
Fig. 1N

Nabers
WD-147-029-SD
WD-147-029-HH
Fig. 2N

Color-Coded







WD-150-001-SD WD-150-002-SD WD-150-003-SD WD-150-004-SD WD-150-005-SD WD-150-006-SD WD-150-007-SD WD-150-008-SD
 WD-150-001-HH WD-150-002-HH WD-150-003-HH WD-150-004-HH WD-150-005-HH WD-150-006-HH WD-150-007-HH WD-150-008-HH
 Fig. 1 Fig. 2 Fig. 00 Fig. 0 Fig. 1 Fig. 2 Fig. 3 Fig. 4



WD-151-017-SD
 WD-151-017-HH
 Set Of 5 Pieces

WD-151-009-SD
 WD-151-009-HH
 Fig. 1

WD-151-010-SD
 WD-151-010-HH
 Fig. 2

WD-151-011-SD
 WD-151-011-HH
 Fig. 3

WD-151-012-SD
 WD-151-012-HH
 Fig. 0

WD-151-013-SD
 WD-151-013-HH
 Fig. 1

WD-151-014-SD
 WD-151-014-HH
 Fig. 2

WD-151-015-SD
 WD-151-015-HH
 Fig. 3

WD-151-016-SD
 WD-151-016-HH
 Fig. 4



WD-152-018-SD
WD-152-018-HH
Fig. 1

WD-152-019-SD
WD-152-019-HH
Fig. 2

WD-152-020-SD
WD-152-020-HH
Fig. 3




WD-152-021-SD
WD-152-021-HH
Fig. 1

WD-152-022-SD
WD-152-022-HH
Fig. 2

WD-152-023-SD
WD-152-023-HH
Fig. 3





Spreader
WD-154-001-SD
WD-154-001-HH
Fig. 1

Spreader
WD-154-002-SD
WD-154-002-HH
Fig. 2


Spreader
WD-154-003-SD
WD-154-003-HH
Fig. 3

Spreader
WD-154-004-SD
WD-154-004-HH
Fig. 4

Spreader
WD-154-005-SD
WD-154-005-HH
Fig. 5

Spreader
WD-154-006-SD
WD-154-006-HH
Fig. 6

Spreader
WD-154-007-SD
WD-154-007-HH
Fig. 7



WD-155-016-SD
WD-155-016-HH
Endodontic Set
With Measuring Plate



Luks
WD-155-008-SD WD-155-009-SD WD-155-010-SD WD-155-011-SD
WD-155-008-HH WD-155-009-HH WD-155-010-HH WD-155-011-HH
Fig. 1 Fig. 2 Fig. 3 Fig. 4
Root-Canal Pluggers



Vertical Condenser WD-155-012-SD WD-155-013-SD WD-155-014-SD WD-155-015-SD
WD-155-012-HH WD-155-013-HH WD-155-014-HH WD-155-015-HH
Fig. 1 Fig. 2 Fig. 3 Fig. 4



WD-156-017-SD
WD-156-017-HH
Fig. 6/1

WD-156-018-SD
WD-156-018-HH
Fig. 6/17
Root-Canal Explorers

WD-156-019-SD
WD-156-019-HH
Fig. 6/23




WD-156-020-SD
WD-156-020-HH
Fig. 31L

WD-156-021-SD
WD-156-021-HH
Fig. 32L
Root-Canal Excavators

WD-156-022-SD
WD-156-022-HH
Fig. 33L




WD-157-023
15 cm / 6"
Fig. 2
Self-Locking



Stieglitz
WD-157-024
14 cm / 5½"
Fig. 1
Silver Point Forceps
Straight Tip




Stieglitz
WD-157-025
14 cm / 5½"
Fig. 2
Angled Tip



WD-157-026
13 cm / 5"
Fig. 1

Pliers For Removing Broken Broaches




WD-157-027
13 cm / 5"
Fig. 2



Ainsworth
WD-158-001
Punch Forceps



Ainsworth
WD-158-002
Punch Forceps




Ivory
WD-158-003
With 1 Row Of Holes
Punch Forceps




Ivory
WD-158-004
With 2 Row Of Holes
Punch Forceps



Ivory
WD-158-005
Light Weight




Brewer
WD-158-006



Stocke
WD-158-007



Palmer
WD-158-008



WD-159-011
Complete Set Of 3



Elastic Separating Plier
WD-159-012



Elastic Separating Plier
WD-159-013



Mershon
WD-159-014
14.5 cm / 5 3/4"
Band Pusher



WD-159-015
15 cm / 6"
Band Setter With
Tin Insert



WD-160-016
Molar Clamps - Winged



WD-160-017



WD-160-018



WD-160-019



WD-160-020



WD-160-021



WD-160-022



WD-160-023



WD-160-024
Set Of 8 Pieces Stainless Steel Clamps
Stainless Steel Bicuspid Clamps-Wingless



WD-160-025
Molar Clamps - Winged




WD-160-026



WD-160-027



WD-160-028



WD-160-029



WD-160-030



WD-160-031
Labial Clamps - Winged



WD-160-032



WD-160-033
Cervical Clamps





WD-162-052
Fig. 28



WD-162-053
Fig. 200



WD-162-054
Fig. 201



WD-162-055
Fig. 202



WD-162-056
Fig. 203



WD-162-057
Fig. 204



WD-162-058
Fig. 205



WD-162-059
Fig. 3



WD-162-060
Fig. W3



WD-162-061
Fig. 4



WD-162-062
Fig. 56A



WD-162-063
Fig. 7



WD-162-064
Fig. 7A



WD-162-065
Fig. 8



WD-162-066
Fig. 8A



WD-162-067
Fig. W8A



WD-162-068
Fig. 10



WD-162-069
Fig. 11




WD-162-070
Fig. 12A



WD-162-071
Fig. 13A

Rubber Dam Clamp / Bicuspid Clamps With Dam-Engaging Projection (Winged)







WD-164-094
Rubber Dam Instrument Kit

Contents

- 1 Instrument Tray Base-Stainless Steel
- 1 Instrument Tray Lid-Stainless Steel
- 1 Insert Frame-Stainless
- 1 Rubber Dam Frame-Stainless Steel
- 1 Clamp-Forceps-Brewer
- 1 Punch - Forceps - Ivory
- 1 Rubber Dam Clamp Tray - Stainless
- 1 Bicuspid Clamp No 60 - Wingless
- 1 Molar Clamp No 26 - Wingless
- 1 Molar Clamp No 51 - Wingless
- 1 Bicuspid Clamp No 29 - Wingless
- 1 Molar Clamp No 202 - Winged
- 1 Molar Clamp No 205 - Winged
- 1 Cervical Clamp No 212
- 1 Labial Clamp No 210 - Winged



WD-165-095
Stainless Steel Rubber Dam Clamp Tray With 8 Diferent Clamps



WD-165-096
Stainless Steel Rubber Dam Clamps Tray (Without Clamps)



WD-166-097
Autoclavable Plastic Amalgam Holders



WD-166-098
Rubber Dam Vacuum
6" x 6" Cut



WD-166-099
2 mm Ø
Metal Tip



Amalgam Well
WD-166-100
3.5 cm / 1½" Ø
3.5 cm / 1½" Height



Amalgam Carriers
WD-166-101
2 mm Ø



Amalgam Carriers
WD-166-102
2 mm Ø



Amalgam Carriers
WD-166-103
2 mm Ø

2.0 mm Ø



1.5 mm Ø

WD-167-001
Medium/Small

2.0 mm Ø



2.5 mm Ø

WD-167-002
Medium/Large

2.0 mm Ø



3.0 mm Ø

WD-167-003
Medium/Jumbo

2.5 mm Ø



3.0 mm Ø

WD-167-004
Large/Jumbo

1.5 mm Ø



WD-167-005
Small

2.0 mm Ø



WD-167-006
Medium

2.5 mm Ø




WD-167-007
Large


3.0 mm Ø



WD-167-008
Jumbo




Amalgam Instruments



Westcott
WD-169-017-SD
WD-169-017-HH
Fig. 1
Burnishers

Amalgam Instruments


Serrated



WD-170-018-SD WD-170-019-SD WD-170-020-SD WD-170-021-SD
WD-170-018-HH WD-170-019-HH WD-170-020-HH WD-170-021-HH
Fig. 1 Fig. 2 Fig. 3 Fig. 4

WD-170-026-SD WD-170-027-SD
WD-170-026-HH WD-170-027-HH
Fig. 1 Fig. 2


Smooth



WD-170-022-SD WD-170-023-SD WD-170-024-SD WD-170-025-SD
WD-170-022-HH WD-170-023-HH WD-170-024-HH WD-170-025-HH
Fig. 1L Fig. 2L Fig. 3L Fig. 4L

WD-170-028-SD
WD-170-028-HH
Fig. 3


Amalgam Instruments



Amalgam Instruments



WD-172-046-SD
WD-172-046-HH
Set Of 6 Pieces



WD-172-040-SD
WD-172-040-HH
Fig. 1

WD-172-041-SD
WD-172-041-HH
Fig. 2

WD-172-042-SD
WD-172-042-HH
Fig. 3


WD-172-043-SD
WD-172-043-HH
Fig. 4

WD-172-044-SD
WD-172-044-HH
Fig. 5

WD-172-045-SD
WD-172-045-HH
Fig. 6

Burnishers





WD-174-057-SD WD-174-058-SD WD-174-059-SD
WD-174-057-HH WD-174-058-HH WD-174-059-HH
Fig. 1 Fig. 2 Fig. 3

WD-174-064-SD WD-174-065-SD WD-174-066-SD
WD-174-064-HH WD-174-065-HH WD-174-066-HH
Fig. 1 Fig. 2 Fig. 3

WD-174-056-SD WD-174-060-SD WD-174-061-SD WD-174-062-SD WD-174-063-SD WD-174-067-SD WD-174-068-SD
WD-174-056-HH WD-174-060-HH WD-174-061-HH WD-174-062-HH WD-174-063-HH WD-174-067-HH WD-174-068-HH
Fig. 1 Fig. 2 Fig. 1 Fig. 2 Fig. 1 Fig. 2 Fig. 3



WD-175-078-SD
WD-175-078-HH
Set of 5 Instruments



WD-175-069-SD
WD-175-069-HH
Fig. 1

WD-175-070-SD
WD-175-070-HH
Fig. 2

WD-175-071-SD
WD-175-071-HH
Fig. 3



WD-175-075-SD
WD-175-075-HH
Fig. 2

WD-175-076-SD
WD-175-076-HH
Fig. 3




WD-175-072-SD
WD-175-072-HH
Fig. 1


WD-175-073-SD
WD-175-073-HH
Fig. 2

WD-175-074-SD
WD-175-074-HH
Fig. 3



WD-175-077-SD
WD-175-077-HH
Fig. 1





Glick
WD-177-091A-SD
WD-177-091A-HH
Fig. 1


Glick
WD-177-091B-SD
WD-177-091B-HH
Fig. 2

WD-177-091-SD
WD-177-091-HH
Fig. 4

WD-177-092-SD
WD-177-092-HH
Fig. 5

WD-177-093-SD
WD-177-093-HH
Fig. W-3

Gregg
WD-177-094-SD
WD-177-094-HH
Fig. 4/5



CF 5A
WD-178-095A-SD
WD-178-095A-HH

CF 4A
WD-178-095B-SD
WD-178-095B-HH

Frahm
WD-178-095-SD
WD-178-095-HH
Fig. 1

USC
WD-178-096-SD
WD-178-096-HH
Fig. 2

Interproximal
WD-178-097-SD
WD-178-097-HH
Fig. 1

Interproximal
WD-178-098-SD
WD-178-098-HH
Off-Angle
Fig. 2



WD-179-099-SD
WD-179-099-HH
Fig. SHO-A



WD-179-100-SD
WD-179-100-HH
Fig. 3-SHO



WD-179-101-SD
WD-179-101-HH
Fig. 104



Wiland
WD-179-102-SD
WD-179-102-HH
Fig. 8

Cleoids-Discoids



WD-180-103-SD
WD-180-103-HH
Fig. 1




WD-180-104-SD
WD-180-104-HH
Fig. 2




WD-180-105-SD
WD-180-105-HH
Fig. 3



WD-180-106-SD
WD-180-106-HH
Fig. 4



Behrman
WD-181-107-SD
WD-181-107-HH



Black
WD-181-108-SD
WD-181-108-HH



Back-Action
WD-181-109-SD
WD-181-109-HH
1/12 Plugger


WD-182-110-SD
WD-182-110-HH


Black
WD-182-111-SD
WD-182-111-HH
1-2 (1.5 mm - 2.5 mm)


Bodkin
WD-182-112-SD
WD-182-112-HH
1 (0.7 mm - 1 mm)

Bodkin
WD-182-113-SD
WD-182-113-HH
1 (0.7 mm - 1 mm)

Bodkin
WD-182-114-SD
WD-182-114-HH
2 (1.5 mm - 2 mm)

WD-182-115-SD
WD-182-115-HH
2 mm - 3 mm





Instruments For Subgingival Application Of Tinctures & Medicaments



Orban Loop
WD-185-001-SD
WD-185-001-HH
Wire Retractor

WD-185-002-SD
WD-185-002-HH

WD-185-003-SD
WD-185-003-HH


WD-185-004-SD
WD-185-004-HH

WD-185-005-SD
WD-185-005-HH

WD-185-006-SD
WD-185-006-HH

Luniatschek
WD-185-007-SD
WD-185-007-HH
2 mm
WD-185-008-SD
WD-185-008-HH
3 mm
Gauze Packers

Gingival Margin Trimmers



Black
WD-186-010-SD
WD-186-010-HH
Fig. 1, 53-54

Black
WD-186-011-SD
WD-186-011-HH
Fig. 2, 51-52

Wedelstaet
WD-186-012-SD
WD-186-012-HH
Fig. 1/2

Wedelstaet
WD-186-013-SD
WD-186-013-HH
Fig. 3/4

Wedelstaet
WD-186-014-SD
WD-186-014-HH
Fig. 5/6



Black
WD-186-009-SD
WD-186-009-HH
Fig. 1, 26




Black
WD-186-015-SD
WD-186-015-HH
Fig. 2, 27



Black
WD-186-016-SD
WD-186-016-HH
Fig. 3, 28




Black
WD-186-017-SD
WD-186-017-HH
Fig. 4, 29




Black
WD-186-018-SD
WD-186-018-HH
Fig. 1, 77-78



Black
WD-186-019-SD
WD-186-019-HH
Fig. 2, 79-80



WD-187-027-SD
WD-187-027-HH
Set Of 5 Pieces



WD-187-020-SD
WD-187-020-HH
Fig. 0

WD-187-021-SD
WD-187-021-HH
Fig. 1

WD-187-022-SD
WD-187-022-HH
Fig. 2

WD-187-023-SD
WD-187-023-HH
Fig. 3

WD-187-024-SD
WD-187-024-HH
Fig. 4

WD-187-025-SD
WD-187-025-HH
Fig. 1

WD-187-026-SD
WD-187-026-HH
Fig. 2



WD-188-028-SD
WD-188-028-HH
Fig. 6



WD-188-029-SD
WD-188-029-HH
Fig. 14



WD-188-030-SD
WD-188-030-HH
Fig. 17



Oval Spoon
WD-188-031-SD
WD-188-031-HH
Fig. 17W



WD-189-032-SD
WD-189-032-HH
Fig. 18



Oval Spoon
WD-189-033-SD
WD-189-033-HH
Fig. 18W




WD-189-034-SD
WD-189-034-HH
Fig. 19



Oval Spoon
WD-189-035-SD
WD-189-035-HH
Fig. 19W

Root Canal Excavators



WD-190-036-SD
WD-190-036-HH
Fig. 38/39

WD-190-037-SD
WD-190-037-HH
Fig. 63/64

WD-190-038-SD
WD-190-038-HH
Fig. 65/66

WD-190-039-SD
WD-190-039-HH
Fig. 31L

WD-190-040-SD
WD-190-040-HH
Fig. 32L

WD-190-041-SD
WD-190-041-HH
Fig. 33L

Bone Chisels & Gouges



Ochsenbein
WD-191-001
Fig. 1



Ochsenbein
WD-191-002
Fig. 2



Wakefield
WD-191-003
Fig. 1W



Wakefield
WD-191-004
Fig. 2W



Buckley
WD-191-005
Fig. 3



Gardner
WD-192-006
Fig. 1



Gardner
WD-192-007
Fig. 2



Gardner
WD-192-008
Fig. 5



Gardner
WD-192-009
Fig. 5R



Gardner
WD-192-010
Fig. 6




WD-193-011
3 mm



WD-193-012
4 mm



WD-193-013
5 mm



WD-193-014
6 mm

Partsch
13.5 cm / 5 1/4"
Bone Chisels



WD-193-015
3 mm



WD-193-016
4 mm




WD-193-017
5 mm



WD-193-018
6 mm

Partsch
13.5 cm / 5 1/4"
Gouges



WD-193-019
3 mm



WD-193-020
4 mm



WD-193-021
5 mm




WD-193-022
6 mm

Partsch
17 cm / 6 3/4"
Bone Chisels



WD-193-023
3 mm



WD-193-024
4 mm



WD-193-025
5 mm



WD-193-026
6 mm

Partsch
17 cm / 6 3/4"
Gouges




WD-193-027
18 cm / 7"
Bone Chisel



Gerzog
WD-194-028
19 cm / 7½"
Weight 10 oz. / 285 gr.
Head Weight 7 oz. / 200 gr.
Head Diameter 25 mm



Partsch
WD-194-029
19 cm / 7½"
Mallets, Lead Filled



Mead
WD-194-030
19.5 cm / 7¾"
With Plastic Inserts




Mead
WD-194-031
17 cm / 6¾"
With Interchangeable Ferrozell
Faces

WD-193-031R
Ferozell Faces Only, In Pairs



Tomati
WD-194-032
19 cm / 7½"
Gum Guillotine Forceps



Mini Friedman
WD-194-033
12 cm / 4¾"
Small Beaks For Use In Limited-Access Areas



Goldman-Fox
WD-194-034
10 cm / 4"





Periosteals & Elevators



König
WD-197-044
15.5 cm / 6"

Partsch
WD-197-045
14.5 cm / 5¾"

Seldin Retractor
WD-197-046
Fig. 22


Seldin Retractor
WD-197-047
Fig. 23

Williger
WD-197-048
17 cm / 6¾"

Williger
WD-197-049
16 cm / 6¼"

Williger
WD-197-050
16 cm / 6¼"

Williger
WD-197-051
13 cm / 5"



Sugarman
WD-198-052-SD
WD-198-052-HH
Fig. 1S/2S
Files

Sugarman
WD-198-053-SD
WD-198-053-HH
Fig. 3S/4S
Files

Schluger
WD-198-054-SD
WD-198-054-HH
Fig. 9/10
Files

Buck
WD-198-055-SD
WD-198-055-HH
Fig. 11/12
Files

Ochsenbein
WD-198-056-SD
WD-198-056-HH
Fig. 3

Ochsenbein
WD-198-057-SD
WD-198-057-HH
Fig. 4

Rhodes Back-Action
WD-198-058-SD
WD-198-058-HH
Fig. 36/37



TG Chisel
WD-199-059-SD
Fig. 1/2

TG Chisel
WD-199-060-HH
Fig. 1/2

TG Chisel
WD-199-061-SD
WD-199-061-HH
Fig. 1

TG Chisel
WD-199-062-SD
WD-199-062-HH
Fig. 2

Gingivectomy Knives



Kirkland
WD-200-063-SD
WD-200-063-HH
Fig. 1
(15K)

Kirkland
WD-200-064-SD
WD-200-064-HH
Fig. 2
(16K)

Orban
WD-200-065-SD
WD-200-065-HH
Fig. 1

Orban
WD-200-066-SD
WD-200-066-HH
Fig. 2

Kirkland
WD-200-067-SD
WD-200-067-HH
Fig. 1/2
(15K / 16K)

Orban
WD-200-068-SD
WD-200-068-HH
Fig. 1/2

Gingivectomy Knives



Goldman-Fox
WD-201-069-SD
WD-201-069-HH
Fig. 7

Goldman-Fox
WD-201-070-SD
WD-201-070-HH
Fig. 8

Goldman-Fox
WD-201-071-SD
WD-201-071-HH
Fig. 9

Goldman-Fox
WD-201-072-SD
WD-201-072-HH
Fig. 11

Buck
WD-201-073-SD
WD-201-073-HH
Fig. 3/4

Buck
WD-201-074-SD
WD-201-074-HH
Fig. 5/6

Gingivectomy Knives




Sanders
WD-202-075-SD
WD-202-075-HH
Fig. 1/2

Sanders
WD-202-076-SD
WD-202-076-HH
Fig. 4/5

(University Of South California)
WD-202-077-SD
WD-202-077-HH
Fig. 19/20
USC Towner

Crane-Kaplan
WD-202-078-SD
WD-202-078-HH
Fig. CK3

Curettes & Scalers



WD-203-086-SD
WD-203-086-HH
Set Of 7 Pieces
Columbia University Set



Curette
WD-203-079-SD
WD-203-079-HH
Fig. HU.1/2



Columbia
WD-203-080-SD
WD-203-080-HH
Fig. 2R/2L



Columbia
WD-203-081-SD
WD-203-081-HH
Fig. 4R/4L



Columbia
WD-203-082-SD
WD-203-082-HH
Fig. 13/14



WD-203-083-SD
WD-203-083-HH
Fig. 411/412



WD-203-084-SD
WD-203-084-HH
Fig. 413/414



Zerfing-Chisel
WD-203-085-SD
WD-203-085-HH

Curettes & Scalers



Langer
WD-204-087-SD
WD-204-087-HH
Fig..1/2

Langer
WD-204-088-SD
WD-204-088-HH
Fig..3/4


Barnhart
WD-204-089-SD
WD-204-089-HH
Fig..1/2

Barnhart
WD-204-090-SD
WD-204-090-HH
Fig..5/6

WD-204-091-SD
WD-204-091-HH
Fig. ME-1
[10/11]

WD-204-092-SD
WD-204-092-HH
Fig. 3/4
Rule [23/24]


Curettes & Scalers



WD-205-100-SD
WD-205-100-HH
Set Of 7 Pieces



Curettes & Scalers



Gracey
WD-206-101-SD
WD-206-101-HH
Fig. 11/14
Mesial-Distal

Gracey
WD-206-102-SD
WD-206-102-HH
Fig. 12/13
Mesial-Distal


De Marco
WD-206-103-SD
WD-206-103-HH
Fig. 1

De Marco
WD-206-104-SD
WD-206-104-HH
Fig. 2


Indiana University
WD-206-105-SD
WD-206-105-HH
Fig. 13/14
Anterior

Indiana University
WD-206-106-SD
WD-206-106-HH
Fig. 17/18
Posterior

Curettes & Scalers



WD-207-114-SD
WD-207-114-HH
Set Of 7 Pieces



Goldman-Fox
WD-207-107-SD
WD-207-107-HH
Fig. 1

Goldman-Fox
WD-207-108-SD
WD-207-108-HH
Fig. 2

Goldman-Fox
WD-207-109-SD
WD-207-109-HH
Fig. 3


Goldman-Fox
WD-207-110-SD
WD-207-110-HH
Fig. 4

Goldman-Fox
WD-207-111-SD
WD-207-111-HH
Fig. 5

Goldman-Fox
WD-207-112-SD
WD-207-112-HH
Fig. 6

Goldman-Fox
WD-207-113-SD
WD-207-113-HH
Fig. 21

Curettes & Scalers



Mc Call
WD-208-115-SD
WD-208-115-HH
Fig. 11/12

Mc Call
WD-208-116-SD
WD-208-116-HH
Fig. 11A/12A


Mc Call
WD-208-117-SD
WD-208-117-HH
Fig. 13/14

Mc Call
WD-208-118-SD
WD-208-118-HH
Fig. 13S/14S

Mc Call
WD-208-119-SD
WD-208-119-HH
Fig. 17S/18S

Mc Call
WD-208-120--SD
WD-208-120-HH
Fig. 19/20

Curettes & Scalers



Taylor
WD-209-121-SD
WD-209-121-HH
Fig. 2/3

WD-209-122-SD
WD-209-122-HH
Fig. 204

WD-209-123-SD
WD-209-123-HH
Fig. 204S

WD-209-124-SD
WD-209-124-HH
Fig. 2/3

Crane-Kaplan
WD-209-125-SD
WD-209-125-HH
Fig. 6

Younger-Good
WD-209-126-SD
WD-209-126-HH
Fig. 7/8

Curettes & Scalers



Hirschfeld
WD-210-127-SD
WD-210-127-HH
Fig. 3/7
File

Jaquette
WD-210-128-SD
WD-210-128-HH
Fig. 30/33

Jaquette
WD-210-129-SD
WD-210-129-HH
Fig. 31/32

Jaquette
WD-210-130-SD
WD-210-130-HH
Fig. 34/35

Jaquette
WD-210-131-SD
WD-210-131-HH
Fig. 37/38



WD-211-132
Periodontal Instrument Set

Contents


- Williams Periodontal Pocket Probe
- Goldman-Fox Universal Scaler Fig. 1
- Mc Call Scaler Fig. 13/14
- Crane-Kaplan Scaler Fig. 6
- Towner Scaler Fig. U 15
- Gracey Curette Fig. 5/6
- Gracey Curette Fig. 11/12
- Gracey Curette Fig. 13/14
- Kirkland Gingivectomy Knife Fig. 15K/16K
- Orban Gingivectomy Knife Fig. 1/2
- Buck Gingivectomy Fig. 5/6
- Molt Periosteal Elevator Fig M9
- Goldman-Fox Gum Scissors-Curved, Saw Edge
- Goldman-Fox Periodontal Pocket Maker-Right Side
- Goldman-Fox Periodontal Pocket Maker-Left Side



WD-212-133
Periodontal Set

Contents


- Instrument Tray
- Ochsenbein
- Rhodes Back-Action
- Knife Kirkland
- Knife Buck
- Knife Goldman-Fox
- Knife Sanders
- Curette Langer
- Scaler Bernhart
- Scaler Rule
- Scaler De Marco
- Scaler Goldman-Fox
- Curette Mccall
- Curette Taylor
- Curette Jaquette
- Scaler Towner Jaquette
- File Hirschfeld
- File Sugarmann
- File Buck
- Chisel
- De Marco



WD-212-134
Periodontal Set

Content

- Periosteals Molt
- Periosteals Prichard
- Periosteals Goldman-Fox
- Knives Kirkland
- Mallet-Mead
- Needle Holder Castroviejo
- Needle Holder Mathieu
- Scalpel Handle
- Instrument Tray



WD-212-135
Periodontal Set

Content

- Scalpel Handle
- Scalpel Blades (100 Pcs)
- Retractor Sternberg
- Scissor Goldman-Fox
- Periosteals Molt
- Periosteals Prichard
- Needle Holder Crile Wood
- Tissue Plier Semken Taylor
- Explorer
- Handle
- Mouth Mirror
- Instrument Tray



Towner
WD-213-001-SD
WD-213-001-HH
Fig. U15


Towner Jaquette
WD-213-002-SD
WD-213-002-HH
Fig. U15/33

Towner Jaquette
WD-213-003-SD
WD-213-003-HH
Fig. U15/30

WD-213-004-SD
WD-213-004-HH
Fig. H5/33

WD-213-005-SD
WD-213-005-HH
Fig. H6/H7

WD-213-006-SD
WD-213-006-HH
Fig. 1/2





Hygienist
WD-214-024
Fig. 4



Black
WD-214-028
Fig. 3



Black
WD-214-032
Fig. 4




WD-214-020



WD-214-021



WD-214-022



Zerfing
WD-214-023



Hygienist



Hygienist
WD-214-025
Fig. 5



Hygienist
WD-214-026
Fig. 8



Hygienist
WD-214-027
Fig. 9



Black



Black
WD-214-029
Fig. 7



Black
WD-214-030
Fig. 9



Black
WD-214-031
Fig. 11



Black
WD-214-033
Fig. 8



Black
WD-214-034
Fig. 10



Black
WD-214-035
Fig. 12



McCall
WD-214-036
Fig. 1



McCall
WD-214-037
Fig. 2



McCall
WD-214-038
Fig. 3



McCall
WD-214-039
Fig. 4



McCall
WD-214-040
Fig. 5



McCall
WD-214-041
Fig. 6



McCall
WD-214-042
Fig. 7



McCall
WD-214-043
Fig. 8



McCall
WD-214-044
Fig. 9



McCall
WD-214-045
Fig. 10



McCall
WD-214-046
Fig. 11



McCall
WD-214-047
Fig. 12



McCall
WD-214-048
Fig. 13



McCall
WD-214-049
Fig. 14



McCall
WD-214-050
Fig. 17



McCall
WD-214-051
Fig. 18



McCall
WD-214-052
Fig. 19



McCall
WD-214-053
Fig. 20



McCall
WD-214-054
Fig. 21




McCall
WD-214-055
Fig. 22




McCall
WD-214-056
Fig. 25





McCall
WD-214-057
Fig. 26






McCall











Hemingway
WD-219-027-SD
WD-219-027-HH
Fig. 1

Hemingway
WD-219-028-SD
WD-219-028-HH
Fig. 2

Hemingway
WD-219-029-SD
WD-219-029-HH
Fig. 3

Hemingway
WD-219-030-SD
WD-219-030-HH
Fig. 4




Lucas
WD-219-031-SD
WD-219-031-HH
Fig. 85

Lucas
WD-219-032-SD
WD-219-032-HH
Fig. 86

Lucas
WD-219-033-SD
WD-219-033-HH
Fig. 87

Lucas
WD-219-034-SD
WD-219-034-HH
Fig. 88




WD-220-001
Crown Remover, Handle Only
WD-220-005
Complete Set With 3 Attachments

Attachments



WD-220-006
Crown Spreader, Crown Remover.
Special Instrument For Crown Spreading
After Separating, Simple & Safe Handling.







WD-223-018
Fig. 1




WD-223-019
Fig. 2




Universal
WD-223-020
Fig. 3




WD-223-017
Fig. 0
Crown Spreaders




a	b	c	d
4.4	5.0	3.4	3.8 mm




a	b	c	d
3.3	4.0	2.3	3.0 mm



Böhm
WD-224-021
 14.5 cm / 5 3/4"
 Telescope Crown Pliers



Böhm
WD-224-022
 11.5 cm / 4 1/2"



WD-225-023
15 cm / 6"
Paper Forceps



Miller
WD-225-024
15 cm / 6"
Paper Forceps



Ivory 8
WD-226-001 5 mm
Fig. 5
WD-226-002 6 mm
Fig. 6
WD-226-003 7 mm
Fig. 7



Ivory 8 N
WD-226-004



Didner 8 P
WD-226-005
Adjustable



Tofflemire-Universal
WD-226-006
Fig. 1




Tofflemire-Junior
WD-226-007
Fig. 2




Tofflemire-Senior
WD-226-008
Fig. 3

For Molars



For Molars & Premolars




Standard Package = 1 Dozen



WD-227-016 5 mm x 100 cm
WD-227-017 6 mm x 100 cm
WD-227-018 7 mm x 100 cm
For Matrix Retainers



Tofflemire
WD-227-019
Assorted In Dozen Package



Tofflemire
WD-227-020
Fig. 1

Tofflemire
WD-227-021
Fig. 2

Tofflemire
WD-227-022
Fig. 3
Standard Package = 1 Dozen



3 mm Ø



WD-228-001
18 cm / 7 $\frac{1}{8}$ "
Suction Tubes
Small

4 mm Ø



WD-228-002
18 cm / 7 $\frac{1}{8}$ "
Suction Tubes
Medium

5 mm Ø



WD-228-003
18 cm / 7 $\frac{1}{8}$ "
Suction Tubes
Large

2 mm Ø



Frazier
WD-228-004
15 cm / 6"
Fr. 6

2.67 mm Ø



Frazier
WD-228-005
15 cm / 6"
Fr. 8

3.33 mm Ø



Frazier
WD-228-006
15 cm / 6"
Fr. 10

4 mm Ø



Frazier
WD-228-007
15 cm / 6"
Fr. 12



WD-229-008



WD-229-009



WD-229-010



WD-229-012
Autoclavable Adaptador

WD-229-011
Aspirator Tips Specially
Shaped For Surgical Use.
Autoclavable

Disposable Saliva Ejectors With Metal Wire, Smooth Edges, Will Not Aspirate Tissue, Flexible, Comfortable, Bonded & Non Bonded Tip.



Instruments Trays - Stainless Steel | Bandejas De Acero Inoxidable
Norm-Tray Aus Rosfreiem Edeltah - Normgröße



WD-230-001
Perforated Instrument Tray Complete with Insert Frame For 16 Instruments

WD-230-002
Lid Perforated: Dimensions. 288 x 187 x 29 mm

WD-230-003
Perforated Base: Dimensions. 284 x 183 x 17 mm

WD-230-004
Perforated Instrument Tray Complete with Insert Frame For 23 Instruments.



WD-230-005
Instrument Tray Complete with Perforated Base, Lid
& Insert Frame For 16 Instruments

WD-230-006
Lid: Dimensions. 288 x 187 x 29 mm

WD-230-003
Perforated Base. Dimensions: 284 x 183 x 17 mm

WD-230-007
Instrument Tray Complete with Perforated Base, Lid
& Insert Frame For 23 Instruments




WD-231-008
Instrument Tray Complete with Base, Lid & Insert
Frame For 16 Instruments

WD-231-008L
Lid: Dimensions. 288 x 187 x 29 mm


WD-231-008B
Base: Dimensions. 284 x 183 x 17 mm

WD-231-009
Instrument Tray Complete With Base, Lid & Insert Frame For 23 Instruments



WD-231-010
Insert Frame:
For 16 Instruments

WD-231-011
Insert Frame:
For 23 Instruments



WD-232-012
Instrument Tray Complete with Base, Lid & Insert Frame

Lid : WD-232-013
Base : WD-232-014
Insert Frame : WD-232-015

WD-232-016
Instrument Tray Complete with Base, Lid & Insert Frame

Lid : WD-232-017
Base : WD-232-018
Insert Frame : WD-232-015



WD-232-019
Instrument Tray Complete with Base, Lid & Insert Frame

Lid : WD-232-020
Base : WD-232-021
Insert Frame : WD-232-022



WD-233-001
Lid-Dimensions. 204 x 105 x 54 mm



WD-233-002
Lid-Dimensions. 204 x 105 x 54 mm



WD-233-003
Lid-Dimensions. 204 x 105 x 54 mm



WD-233-004
Lid-Dimensions. 204 x 105 x 54 mm



WD-233-005
Aluminium Endodontic MAXI Boxes, Complete With
Lid-Dimensions. 204 x 105 x 54 mm



WD-233-006
Lid-Dimensions. 125 x 105 x 54 mm



WD-233-007
Aluminium Endodontic MIDI Boxes, Complete With
Lid-Dimensions. 125 x 105 x 54 mm



WD-234-008
Dimensions. 104 x 44 x 54 mm



WD-234-009
Dimensions. 104 x 44 x 54 mm



WD-234-010
Dimensions. 104 x 44 x 54 mm



WD-234-011
Dimensions. 104 x 44 x 54 mm



WD-234-012
Dimensions. 104 x 44 x 54 mm



WD-234-013
Dimensions. 104 x 44 x 54 mm



WD-234-014
Dimensions. 104 x 44 x 54 mm

Aluminium Endodontic MINI Boxes, Complete With Lid-
Dimensions. 104 x 44 x 54 mm




WD-234-015
Round Aluminium Endodontic Boxes, Complete With Lid
Diameter = 48 mm
Height = 62 mm



WD-234-016
Round Aluminium Endodontic Boxes, Complete With Lid
Diameter = 94 mm
Height = 70 mm



Spirit Lamp
WD-235-017
120 ccm



Spirit Lamp
WD-235-018
120 ccm




WD-235-021
15 cm / 6"
Spare Wick



Spirit Lamp
WD-235-019
120 ccm




WD-235-020
Alcohol Burner



WD-236-022
30 cm / 11¾"



WD-236-023
20 cm / 8"




Stainless Steel Water-Cooled Impression Trays With Retention Rim.



WD-237-001
MAXI Set Of 10 Pieces

Dimension : a x b mm



WD-238-002
68 x 52 mm
U5, S



WD-238-003
69 x 54 mm
U4, NM



WD-238-004
73 x 53 mm
U3, M



WD-238-005
75 x 60 mm
U2, L



WD-238-006
83 x 62 mm
U1, XL

Dimension : a x b mm



WD-238-007
71 x 53 mm
L5, S




WD-238-008
74 x 60 mm
L4, NM




WD-238-009
78 x 57 mm
L3, M




WD-238-010
81 x 61 mm
L2, L



WD-238-011
82 x 63 mm
L1, XL




Stainless Steel Impression Trays With Retention Rim.



WD-239-012
MAXI Set Of 10 Pieces

Dimension : a x b mm



WD-240-013
65 x 51 mm
U5, S



WD-240-014
67 x 53 mm
U4, NM



WD-240-015
71 x 55 mm
U3, M




WD-240-016
74 x 60 mm
U2, L



WD-240-017
81 x 63 mm
U1, XL

Dimension : a x b mm



WD-240-018
69 x 50 mm
L5, S



WD-240-019
72 x 57 mm
L4, NM



WD-240-020
76 x 57 mm
L3, M




WD-240-021
78 x 59 mm
L2, L




WD-240-022
80 x 63 mm
L1,XL




Stainless Steel Impression Trays With Retention Rim.



WD-241-023
MAXI Set Of 12 Pieces




WD-241-024
MIDI Set Of 8 Pieces



WD-241-025
MINI Set Of 6 Pieces

Dimension : a x b mm



WD-242-026
61 x 47 mm
U1, XXS



WD-242-027
64 x 51 mm
U2, XS



WD-242-028
66 x 53 mm
U3, S



WD-242-029
69 x 57 mm
U4, M



WD-242-030
71 x 60 mm
U5, L



WD-242-031
73 x 62 mm
U6, XL

Dimension : a x b mm



WD-242-032
55 x 47 mm
L1, XXS



WD-242-033
60 x 49 mm
L2, XS



WD-242-034
72 x 54 mm
L3, S



WD-242-035
75 x 57 mm
L4, M



WD-242-036
76 x 59 mm
L5, L



WD-242-037
78 x 61 mm
L6, XL

Impression Trays For Crown & Bridge-Work.



WD-242-038
U1P



WD-242-039
L2D



WD-242-040
L4D



WD-242-041
L6D



WD-242-042
L1P




WD-242-043
L2S




WD-242-044
L4S




WD-242-045
L6S




Solid Stainless Steel Impression Trays With Retention Rim.



WD-243-046
MAXI Set Of 12 Pieces




WD-243-047
MIDI Set Of 8 Pieces




WD-243-048
MINI Set Of 6 Pieces

Dimension : a x b mm





Dimension : a x b mm








Solid Stainless Steel Depressed Anterior Impression Trays With Retention Rim.

Dimension : a x b mm




Dimension : a x b mm





Solid, Complete Denture Stainless Steel Impression Trays With Rim.

Dimension : a x b mm



WD-246-068
73 x 60 mm
U64, S

WD-246-069
70 x 64 mm
U63, NM

WD-246-070
76 x 67 mm
U62, M

WD-246-071
79 x 70 mm
U61, L

Dimension : a x b mm




WD-246-072
70 x 57 mm
L69, S


WD-246-073
73 x 60 mm
L68, M

WD-246-074
70 x 64 mm
L67, NL


WD-246-075
86 x 64 mm
L66, L




WD-246-076
MIDI Set of 8 pieces




Solid Regular Stainless Steel Impression Trays Without Retention Rim.



WD-247-077
MAXI Set Of 12 Pieces




WD-247-078
MIDI Set Of 8 Pieces



WD-247-079
MINI Set Of 6 Pieces

Dimension : a x b mm



WD-248-080
61 x 47 mm
U1, XXS



WD-248-081
64 x 51 mm
U2, XS



WD-248-082
66 x 53 mm
U3, S



WD-248-083
69 x 57 mm
U4, M




WD-248-084
71 x 60 mm
U5, L



WD-248-085
73 x 62 mm
U6, XL

Dimension : a x b mm



WD-248-086
55 x 47 mm
L1, XXS



WD-248-087
60 x 49 mm
L2, XS



WD-248-088
72 x 54 mm
L3, S



WD-248-089
74 x 57 mm
L4, M



WD-248-090
76 x 59 mm
L5, L



WD-248-091
78 x 61 mm
L6, XL

Impression Trays For Crown & Bridge-Work



WD-248-092
L2D



WD-248-093
L4D



WD-248-094
L6D




WD-248-095
L2S




WD-248-096
L4S




WD-248-097
L6S




Solid Stainless Steel Impression Trays Without Retention Rim.



WD-249-098
MAXI Set Of 12 Pieces




WD-249-099
MIDI Set Of 8 Pieces



WD-249-100
MINI Set Of 6 Pieces

Dimension : a x b mm



WD-250-101
61 x 47 mm
U1A, XXS



WD-250-102
64 x 51 mm
U2A, XS



WD-250-103
66 x 53 mm
U3A, S



WD-250-104
69 x 57 mm
U4A, M




WD-250-105
71 x 60 mm
U5A, L



WD-250-106
73 x 62 mm
U6A, XL

Dimension : a x b mm



WD-250-107
55 x 47 mm
L1, XXS



WD-250-108
60 x 49 mm
L2, XS



WD-250-109
72 x 54 mm
L3, S




WD-250-110
74 x 57 mm
L4, M




WD-250-111
76 x 59 mm
L5, L



WD-250-112
78 x 61 mm
L6, XL




Perforated Regular Stainless Steel Impression Trays With Retention Rim.




WD-251-113
MAXI Set Of 12 Pieces

U6
U5
U4
U3
U2
U1



WD-251-114
MIDI Set Of 8 Pieces

U6
U5
U4
U3



WD-251-115
MINI Set Of 6 Pieces

U5
U4
U3

Dimension : a x b mm



WD-252-116
63 x 47 mm
U1, XXS




WD-252-117
64 x 51 mm
U2, XS



WD-252-118
66 x 53 mm
U3, S



WD-252-119
69 x 57 mm
U4, M




WD-252-120
71 x 60 mm
U5, L



WD-252-121
73 x 62 mm
U6, XL

Dimension : a x b mm



WD-252-122
55 x 47 mm
L1, XXS



WD-252-123
60 x 49 mm
L2, XS



WD-252-124
72 x 54 mm
L3, S



WD-252-125
74 x 57 mm
L4, M



WD-252-126
76 x 59 mm
L5, L



WD-252-127
78 x 61 mm
L6, XL

Impression Trays For Crown & Bridge-Work



WD-252-128
U1P



WD-252-129
L2D



WD-252-130
L4D



WD-252-131
L6D



WD-252-132
L1P




WD-252-133
L2S




WD-252-134
L4S




WD-252-135
L6S




Perforated Stainless Steel Impression Trays With Retention Rim.



WD-253-136
MAXI Set Of 12 Pieces



WD-253-137
MIDI Set Of 8 Pieces






WD-253-138
MINI Set Of 6 Pieces

Dimension : a x b mm




Dimension : a x b mm







Perforated Stainless Steel Depressed Anterior Impression Trays With Retention Rim.

Dimension : a x b mm



Dimension : a x b mm





Perforated-Complete Denture Stainless Steel Impression Trays With Rim.

Dimension : a x b mm



WD-256-158
73 x 60 mm
U64, S



WD-256-159
70 x 64 mm
U63, NM




WD-256-160
76 x 67 mm
U62, M



WD-256-161
79 x 70 mm
U61, L

Dimension : a x b mm



WD-256-162
70 x 57 mm
L69, S




WD-256-163
73 x 60 mm
L68, M




WD-256-164
70 x 64 mm
L67, NL




WD-256-165
86 x 64 mm
L66, L




WD-256-166
MIDI Set Of 8 Pieces




Perforated Regular Stainless Steel Impression Trays Without Retention Rim.



WD-257-167
MAXI Set Of 12 Pieces



WD-257-168
MIDI Set Of 8 Pieces



WD-257-169
MINI Set Of 6 Pieces

Dimension : a x b mm



WD-258-170
61 x 47 mm
U1, XXS



WD-258-171
64 x 51 mm
U2, XS



WD-258-172
66 x 53 mm
U3, S



WD-258-173
69 x 57 mm
U4, M



WD-258-174
71 x 60 mm
U5, L



WD-258-175
73 x 62 mm
U6, XL

Dimension : a x b mm



WD-258-176
55 x 47 mm
L1, XXS



WD-258-177
60 x 49 mm
L2, XS



WD-258-178
72 x 54 mm
L3, S



WD-258-179
74 x 57 mm
L4, M



WD-258-180
76 x 59 mm
L5, L



WD-258-181
78 x 61 mm
L6, XL

Impression Trays For Crown & Bridge-Work



WD-258-182
L2D



WD-258-183
L4D



WD-258-184
L6D




WD-258-185
L2S




WD-258-186
L4S




WD-258-187
L6S




Perforated Stainless Steel Impression Trays Without Retention Rim.



WD-259-188
MAXI Set Of 12 Pieces



WD-259-189
MIDI Set Of 8 Pieces



WD-259-190
MINI Set Of 6 Pieces

Dimension : a x b mm



WD-260-191
61 x 47 mm
U1A, XXS



WD-260-192
64 x 51 mm
U2A, XS



WD-260-193
66 x 53 mm
U3A, S



WD-260-194
69 x 57 mm
U4A, M




WD-260-195
71 x 60 mm
U5A, L



WD-260-196
73 x 62 mm
U6A, XL

Dimension : a x b mm



WD-260-197
55 x 47 mm
L1, XXS



WD-260-198
60 x 49 mm
L2, XS



WD-260-199
72 x 54 mm
L3, S



WD-260-200
74 x 57 mm
L4, M



WD-260-201
76 x 59 mm
L5, L



WD-260-202
78 x 61 mm
L6, XL



WD-261-215
MAXI Set Of 24 Pieces



WD-261-203
L2S

WD-261-204
L2D

WD-261-205
61 x 49 mm
U1

WD-261-206
65 x 53 mm
U2

WD-261-207
70 x 59 mm
U3

WD-261-208
73 x 63 mm
U4

WD-261-209
L1S

WD-261-210
L1D

WD-261-211
61 x 50 mm
L1


WD-261-212
67 x 52 mm
L2

WD-261-213
72 x 56 mm
L3

WD-261-214
74 x 58 mm
L4



WD-261-228
Set Of 12 Pieces
Plastic Impression Trays
Autoclavable 130° C



WD-261-216
L2S

WD-261-217
L2D

WD-261-218
61 x 49 mm
U1

WD-261-219
65 x 53 mm
U2

WD-261-220
70 x 59 mm
U3

WD-261-221
73 x 63 mm
U4

WD-261-222
L1S

WD-261-223
L1D

WD-261-224
61 x 50 mm
L1

WD-261-225
67 x 52 mm
L2

WD-261-226
72 x 56 mm
L3

WD-261-227
74 x 58 mm
L4



WD-261-237
MIDI Set Of 8 Pieces
Orthodontic Impression Trays
Autoclavable 130° C



WD-261-229
63 x 59 mm
U4




WD-261-230
59 x 55 mm
U3



WD-261-231
56 x 51 mm
U2



WD-261-232
50 x 46 mm
U1



WD-261-233
66 x 56 mm
L4



WD-261-234
62 x 53 mm
L3




WD-261-235
58 x 49 mm
L2





WD-261-236
52 x 44 mm
L1





Dimension : a x b mm





Perforated Regular Stainless Steel Impression Trays Without Retention Rim.



Dimension : a x b mm



WD-262-238
Set Of 11 Pieces



1. WD-263-002 Mounting Plates
2. WD-263-005 Occlusal Table
3. WD-263-006 Incisal Table
4. WD-263-007 Incisal Guide Pin
5. WD-263-008 Incisal Table & Incisal Guide Pin Lock Screws
6. WD-263-009 Mounting Plates & Occlusal Table Lock Screws
7. WD-263-010 Adjustable lateral Excursion Lock Screws
8. WD-263-011 Incisor Indicator Pin
9. WD-263-012 Upper Jaw Member Screws



HI-Top Model
WD-264-013
Chrome



Crown & Bridge
WD-264-014
Chrome



Apex II Type
WD-264-015
Chrome



HI-Top Model
WD-264-016
Brass



Crown & Bridge
WD-264-017
Brass



Apex II Type
WD-264-018
Brass



WD-265-019
Stainless Steel
Finishing



WD-265-020
Black Metallic
Golgen




WD-265-021
Orange Metallic
Golgen



WD-265-022
Brown




WD-265-023
Royal Blue




WD-265-024
Silver Blue Metallic



WD-265-025
Silver Blue Metallic



WD-265-026
Silver Green Metallic



WD-265-027
Copper Black Metallic




WD-265-028
Articulators




WD-265-029
Articulators



WD-266-030



WD-266-031



WD-266-032



Hollenback
WD-267-001-SD
WD-267-001-HH
Fig. 1/2



Hollenback
WD-267-002-SD
WD-267-002-HH
Fig. 3



Hollenback
WD-267-003-SD
WD-267-003-HH
Fig. 3S



Ward
WD-267-004-SD
WD-267-004-HH
Fig. 1



Ward
WD-267-005-SD
WD-267-005-HH
Fig. 1S



Ward
WD-267-006-SD
WD-267-006-HH
Fig. 2



Frahm
WD-267-007-SD
WD-267-007-HH
Fig. 1



USC
WD-267-008-SD
WD-267-008-HH
Fig. 2

Waxing Instruments With Replaceable Points



WD-268-009



WD-268-010



WD-268-011



WD-268-009A




WD-268-010A



WD-268-011A



Maehlum WD-269-012	Hylin WD-269-013	Hylin-Le Cron WD-269-014	Hylin-Ceramica WD-269-015	Hylin WD-269-016	WD-269-017	Le Cron WD-269-018	Le Cron WD-269-019	Le Cron WD-269-020
-----------------------	---------------------	-----------------------------	------------------------------	---------------------	------------	-----------------------	-----------------------	-----------------------



Hylin Carver
WD-270-021
Garded Handle

Le Cron
WD-270-022

Zehle
WD-270-023

Zehle
WD-270-024

Zehle
WD-270-025


Gritman
WD-270-026

Beale
WD-270-027

Roach
WD-270-028

Vehe
WD-270-029

Vehe
WD-270-030



WD-271-036-SD
Set Of 5 Pieces



P.K Thomas
WD-271-031-SD
Fig. 1



P.K Thomas
WD-271-032-SD
Fig. 2




P.K Thomas
WD-271-033-SD
Fig. 3



P.K Thomas
WD-271-034-SD
Fig. 4



P.K Thomas
WD-271-035-SD
Fig. 5



WD-272-042-HH
Colour Coded
Set Of 5 Pieces



P.K. Thomas
WD-272-037-HH
Fig. 1



P.K. Thomas
WD-272-038-HH
Fig. 2



P.K. Thomas
WD-272-039-HH
Fig. 3



P.K. Thomas
WD-272-041-HH
Fig. 5



WD-273-011
Set Of 10 Pieces



WD-273-004 WD-273-005 WD-273-006 WD-273-007

WD-273-001 WD-273-002 WD-273-003

WD-273-008 WD-273-009 WD-273-010



WD-274-012



WD-274-013



WD-274-014



WD-274-015

WD-274-016

WD-274-017

Flexible




WD-274-018

WD-274-019

WD-274-020

Super Flexible



Fahnenstock
WD-275-022
12.5 cm / 5"
Fig. 1

Fahnenstock
WD-275-023
17 cm / 6¾"
Fig. 2

Lessmann
WD-275-024
12.5 cm / 5"
Fig. 1

Lessmann
WD-275-025
17 cm / 6¾"
Fig. 2




WD-276-026
17 cm / 6¾"
Fig. 1



Gritman
WD-276-027
17 cm / 6¾"
Fig. 2



WD-276-028
16 cm / 6¼"
Fig. 3



WD-276-029
16 cm / 6¼"
Fig. 4



WD-277-030
18.5 cm / 7¼"
Fig. 1
Rigid

WD-277-031
18 cm / 7⅞"
Fig. 2
Rigid

WD-277-032
18.5 cm / 7¼"
Fig. 3
Flexible

WD-277-033
19.5 cm / 7¾"
Fig. 4
Flexible

WD-277-034
21.5 cm / 8½"
Fig. 5
Flexible

WD-277-035
20.5 cm / 8"
Fig. 6
Flexible



WD-278-036
19.5 cm / 7¾"
Fig. 7
Flexible

WD-278-037
20 cm / 8"
Fig. 8
Rigid

WD-278-038
22.5 cm / 9"
Fig. 9
Rigid




WD-279-039-R WD-279-039-G WD-279-039-B WD-279-039-F WD-279-039-V WD-279-039-A

WD-279-039
19 cm / 7½"
Spatula Flexible
Complete Set




WD-279-040
16 cm / 6¼"
Plaster Cutting Plier



WD-279-041
19.5 cm / 7¾"
Plaster Cutting Plier



Bone Morselizer
WD-279-042
The Bone Morselizer Is A Hand
Instrument Which Helps In The
Crushing Of Bone Material.



Bone Grafting Instruments, Compactor / Tampers



WD-281-007
Graft Packer, Concave
Tips Diameter 4/6 mm



WD-281-008
Bone Spoon / Packer
Tips Diameter 3 mm



WD-281-009
Graft Packer, Serrated
Tips Diameter 6/8 mm



WD-281-010
Graft Packer, Serrated
Tips Diameter 3/4 mm

Bone Grafting Instruments, Compactor / Tampers


WD-282-011
Graft Packer, Smooth
Tips Diameter 5/8 mm



WD-282-012
Titanium Graft Packer / Placement
Site Depth Gauge 3.3 mm & 4.0 mm Diameter
Tips Marked at 8, 11, 13, 15 & 18 mm



WD-282-013
Graft Material Spoon With
3 mm Serrated Graft Paker
Concave Spoon Tip For Large Grafts



WD-282-014
Graft Material Spoon
Concave Tip For Large Grafts

Stainless Steel Tissue Punch



WD-283-015
Straight



WD-283-016
Straight



WD-283-017
Straight



WD-283-018
Offset

WD-283-019
Offset

WD-283-020
Offset

WD-283-021
Offset

WD-283-022
Offset

WD-283-023
Offset




WD-284-001
22 cm / 8½"
Riveting Mallet



WD-284-002
22 cm / 8½"
Riveting Mallet



WD-284-003
20 cm / 8"
Horn Mallet




WD-285-004
25 cm / 9 7/8"



WD-285-005
Fret Saw Blades For Plaster

WD-285-006
Fret Saw Blades For Metal



WD-285-007
18 cm / 7 1/8"



WD-285-008
0.25 x 3 mm

WD-285-009
0.25 x 5 mm

WD-285-010
0.12 x 5 mm



Vernier
WD-286-001
0-120 mm / 0-8"



WD-286-002
0-200 mm / 0-8"



Jameson
WD-286-003
0-80 mm / 0-3"



Townley
WD-286-004
0-80 mm / 0-3"




WD-287-005
0-100 mm / 0-4"



WD-287-006
0-80 mm / 0-3"



WD-287-007
0-80 mm / 0-3"




Vernier
WD-287-010
0-200 mm / 0-8"



Boley
WD-287-008
0-100 mm / 0-4"




Boley
WD-287-009
0-100 mm / 0-4"




Iwanson
WD-288-011
12 cm / 4 3/4"
For Crowns




Iwanson
WD-288-012
12 cm / 4 3/4"
With Drop Tips For Wax



Castroviejo
WD-288-013
8 cm / 3"
0-25 mm



Castroviejo
WD-288-014
8 cm / 3"
0-25 mm




Castroviejo
WD-288-015
8 cm / 3"
0-25 mm



Castroviejo-Epker
WD-289-016
9 cm / 3 1/2"
0-40 mm



Castroviejo-Epker
WD-289-017
18 cm / 7 1/8"
0-40 mm




WD-289-018 15 cm / 6"
WD-289-019 20 cm / 8"



WD-289-020
0-80 mm



WD-289-021
0-60 mm



WD-289-022
150 mm
Measuring Range & Hole Gauge




WD-290-001
31 cm / 12½"



WD-290-002
43 cm / 16"



WD-290-003
40 cm / 15¾"



WD-290-004
51 cm / 20"



WD-290-005
64 cm / 25"



WD-291-009
26 cm / 10 1/4"
Fig. 6



WD-291-006
16.5 cm / 6 1/2"
Fig. 1



WD-291-007
16 cm / 6 3/4"
Fig. 2



WD-291-008
12 cm / 4 3/4"
Fig. 2



WD-291-010
18 cm / 7 1/8"
Fig. 4



WD-291-011
15 cm / 6"
Fig. 3




WD-291-012
11.5 cm / 4 1/2"
Fig. 3




WD-292-001
Light Wire Bird Beak
Pliers With Groove



Tweed
WD-292-002
Loop Forming Pliers



WD-292-003
Arch Bending Pliers



Jarabak
WD-292-004



WD-292-005
Hollow Chop Contouring Pliers



WD-292-006
Lingual Arch Forming Pliers



WD-293-007
Rectangular Arch Bending Pliers



WD-293-008
Posterior Band Removing Pliers-Long



WD-293-009
Posterior Band Removing Pliers-Short



WD-293-010
Force Module Separating Pliers



WD-293-011
Tying Pliers Coon Style Ligature



WD-293-012
Bird Beak Pliers

Pliers For Orthodontics & Prosthetics | Alicates Para Ortodoncia Y Prótesis

Zangen Für Kieferorthopädie Und Prothetik



WD-294-013
Omega Loop Forming Pliers



WD-294-014
Three Jaw Pliers




WD-294-015
Cutter




Weingart
WD-294-016



Weingart
WD-294-017
Slim Pliers




How
WD-294-018
Straight



How
WD-295-019
Pliers Offset




WD-295-020
Band Seating Pliers



WD-295-021
Adhesive Removing Pliers



WD-295-022
Hard Wire Cutter
15 cm



WD-295-023
Distal End Cutter
Flush Cut



WD-295-024
Distal End Cutter
Universal Cut & Hold

Pliers For Orthodontics & Prosthetics | Alicates Para Ortodoncia Y Prótesis
Zangen Für Kieferorthopädie Und Prothetik



WD-296-025
Pin & Ligature Cutter



WD-296-026
Mini Pin & Ligature Cutter



WD-296-027
Hard Wire Cutter Straight



WD-296-028
Micro-Mini Pin & Ligature Cutter



WD-297-029
E-Z Ligature
Placer



WD-297-030
Angled Direct Bone
Bracket Remover
TC



WD-297-031
Distal End Cutter Safety Hold Is
Idea For Cutting Wires
Within The Patient's Mouth
TC




WD-297-032
Direct Bonding Bracket
Remover (Wide) Plier
TC



WD-297-033
Pin & Ligature Cutter
TC




WD-297-034
Direct Bracket Remover



Nance
WD-298-035
 Loop Closing Plier, For Bending
 Of Stop Loops & Adams' Clamps
 Wire Up To 0.7 mm / 22" x 0.28"



WD-298-036
 Plier For Double Bars




WD-298-037
 Hard Wire Cutter Up To
 56 mm / 0.22" x 71 mm / 0.28"



WD-298-038
 Heavy Duty Wire Cutter Unique Joint
 Wires Up To .080"
 This Cutter Is Not Autoclavable




WD-298-039
 Extra-Heavy Wire Cutter Heavy Duty
 To Cut All Archwires & Retainer
 Wires Up To .060"



WD-299-040
Hard Wire Cutters
0.82 mm / 0.32" To 1.50 mm / 0.57"



WD-299-041
Hard Wire Cutters
0.80 mm / 0.32" To 1.50 mm / 0.57"




Univrsal
WD-299-042
Face Bow Bending & Cutting Pliers
Max 1.80 mm / 0.71"




Magnm
WD-299-043
Face Bow Bending & Cutting Pliers
Max 1.80 mm / 0.71"




WD-300-044
Micro-Minature Light Wire Cutter (With Stop Screw). Stop Screw Prevents Tip Damage Due To Overstressing Cutter. Caution: Maximum Cutting Capacity: Soft Ligature Wire Up To .012\"



WD-300-045
Distal End Cutter - Safety Hold
Maximum Cutting Capacity Archwires
Up To .022\" x .028\". Rectangular
Round Wires Up To .020\".



WD-300-046
Micro-Distal End Cutter - Safety Hold




Nance
WD-300-047
Loop Closing Plier. For Bending Of Stop Loops & Special Plier For Adams' Clamps. Wires Up To 0.7 mm / 0.22\" x 0.28\"




Universal
WD-300-048
Universal Pin Roughing, Cutting & Bending Plier. For Wires Up To 0.035\" . Will Cut Wires Up To 0.028\"



Universal
WD-300-049
Universal Pin, Roughing, Cutting & Bending Plier. For Wires Up To 0.035\". Will Cut Wires Up To 0.028\".



Universal
WD-301-050
Cutting, Max 0.7 mm




Waldsachs
WD-301-051
Wire, Bending & Cutting Forming, Wires Up
To 1 mm / 0.39", Will Cut Wirer Up To 0.7 mm / 0.28"




WD-301-052
Round & Concave Plier
For Loop Forming & Bending
Wires Up To .020".





Reynolds
WD-301-053
Contouring Plier




WD-301-054
Double Rounded Jaw Wire Bending Plier
Wires Up To 0.30".




WD-301-055
Round Nose Wire Bending Plier.
Wire Up To 0.20".




WD-302-056
Angled
Bird Beak Plier For Bending Round, Square Or
Rectangular Wire Up To .030".



Adams
WD-302-057
Adams Plier For Bending Wire Up To .030".



Adams
WD-302-058




WD-302-059
Light Wire Bird Beak Plier For Bending Round,
Square Or Rectangular Wire Up To .016"



WD-302-060
Flat Nose Plier With Smooth Beaks For
Bending Wire Up To .030"




WD-302-061
Bending Wire Up To 0.76 mm / 0.30"



WD-303-062
Bending Wire Up To 0.76 mm / 0.30"



WD-303-063
Bending Wire Up To 0.76 mm / 0.30"



Johnson
WD-303-064


3 mm
4 mm
5 mm
6 mm



Nance
WD-303-065




Young
WD-303-066
Wires Up To 0.7 mm / 0.28"



WD-303-067




Makellops
WD-304-068
Wires Up To 1.5 mm / 0.59"



Angle
WD-304-069
Bird Beak Plier For Bending Round,
Square Or Rectangular Wires Up To .030". With Cutter




Tweed (Angle)
WD-304-070
Tweed (Angle) Arch Bending Plier - Short
- Plier With Blades (10 mm) Designed For
Convenient Handling Of Square Or
Rectangular Wire, The Blades Are Parallel At
0.20" Plier Opening.



Tweed
WD-304-071
Wire Up To 0.7 mm / 0.28"




Weingart
WD-304-072
Utility Plier "Sure-Grip" Serrated
Pads Are Near The Plier's Center Of Rotation
For Better Balance.



Aderer
WD-304-073
Three Jaw Wire Bending Plier. The
Precision-Aligned Tips Insure Consistent
Bends. Gradually Tapered For Delicate Bending.



Aderer
WD-305-074
Three Jaw Wire Bending Plier.
For Wire & Clamps Adjusting & Contouring.



How
WD-305-075
How Plier - Regular - Well Balanced For
Optimum Control. Completely Lip-Safe Beaks.



How
WD-305-076
How Plier - Offset - Utility Plier
With 45° Offset Angle Of Tips.



WD-305-077
Lingual Arch Plier Used To Form Double &
Triple Back Bends In .030" Or .036" Wire for
Insertion in Lingual Sheaths.



De La Rosa
WD-305-078
De La Rosa Arch Forming Plier - Grooved.
The Beak Has Grooves Of
.016", .018" & .022".



De La Rosa
WD-305-079
Arch Forming Plier-Non Grooved -
Modified With A Non Grooved Contouring
Surface For The Forming Of Precise Arches.




WD-306-080
Direct Bond Attachment Remover (Angled).



Tweed
WD-306-081
Tweed Loop Forming Plier. For Making Loops In Square Or Rectangular Wires Up To 0.022" x 0.028". The Replaceable Cylindrical Beak Has Three Diameter Sections: 0.045", 0.060", 0.075"



Tweed
WD-306-082
Omega Loop Forming Plier. For making Perfect Omega Loops In Square Or Rectangular Wires Up To .022" x .028".



WD-306-083
Single Diameter Loop Forming Plier. The Replaceable Cylindrical Beak Has Only One Diameter Section Of .060" & .270" In Length.



WD-306-084
Direct Bonding Bracket Remover (Wide)
The Tips Wedge Between Edges Of Bracket Base & Tooth Surface



WD-306-085
Direct Bonding Bracket Remover With Pad (Wide).




WD-307-086A
4.8 mm (3/16")

WD-307-086B
4.4 mm (1/4")

Replacement Pads For Pliers
WD-307-086A
WD-307-086B

WD-307-086
Direct Bonding Bracket Remover With Pad
The Tip Is Narrow Tip



WD-307-087
Posterior Band Remover
Short



WD-307-088
Posterior Band Remover
Long



WD-307-089



WD-307-090



WD-307-091
Wire Cutter. It Cuts Ligature Wires & Pins Of
Any Size. It Is Strong Enough To Bend
Archwire Ends After Distal Cutting.
Max. Archwire Up To .022" x .028".

Pliers For Orthodontics & Prosthetics | Alicates Para Ortodoncia Y Prótesis
Zangen Für Kieferorthopädie Und Prothetik



Schwarz
WD-308-092
Arrow Clamps Forming & Bending Pliers



Schwarz
WD-308-093
Arrow Clamps Forming & Bending Pliers
Max 0.70 mm / 0.28"



Nance
WD-308-094
Adams Clamps Forming Pliers
Max 0.70 mm / 0.28"




Aderer
WD-308-095
Facebow Loop Forming Pliers
Max 1.60 mm / 0.63"




Mckellops
WD-308-096
Loop Forming Pliers
Max 1.80 mm / 0.32"




Aderer Clamps
WD-309-097
Three Jaw Wire Bending Plier For Contouring
& Bending Wire Up To 76 mm / 30"




Bird Beak Plier
WD-309-098
For Bending Round Square & Rectangle
Wire Up To 76 mm / 30"




Bird Beak Plier
WD-309-099
For Light Wire Bending Round
Square & Rectangular
Wire Up To 0.41 mm / 0.16"



WD-309-100
Double Round Jaw Bending Pliers Hooks &
Wire Up To 0.51 mm / 0.20"



Adams Pliers
WD-309-101
For Hooks & Bending
Wire Up To 0.76 mm / 0.30"



Johnson Pliers
WD-309-102
For Crown Band Contouring



WD-310-103
Falcon Beak Contouring Plier



WD-310-104
For Placements & Removal
Of Arch Wires, Pins &
Other Auxiliaries




WD-310-105
For Placements & Removal
Of Arch Wires, Pins &
Other Auxiliaries In Angles Design



WD-310-106
Contouring Plier



WD-310-107
Bending Plier. In Order To Skirt
Molares Bands. With Grooved
Recommended Wire
0.4 mm / 0.16" x 0.46 mm / 0.18" x 0.56 mm / 0.22"




WD-310-108
Loop Forming In Four Length
3,4,5 & 6 mm Wire Capacity




WD-311-109
Arch Bending Plier. Design For
Convenient Handling Of Square Or
Rectangular Wire, 0.51 mm / 0.20"




WD-311-110
In Order To Skirt Bands & Crowns




WD-311-111
Wire Up To
0.7 mm / 0.20"



WD-311-112
Round & Concave Plier For Loop
Forming & Banding
Wire Up To 0.51 mm / 0.20"




WD-311-113
Banding For Round Square
or Rectangular
Wire Up To 0.76 mm / 0.30"



WD-311-114
For Double Back & Triple Back
In 76 mm / 0.30" x 91 mm / 0.36"
Lingual Arch Wires

Pliers For Orthodontics & Prosthetics | Alicates Para Ortodoncia Y Prótesis


Zangen Für Kieferorthopädie Und Prothetik



WD-312-115
Detailing Step Plier Double-Sided
Offset Beak Permits Bayonet Bends
Creates Right Hand Or Left Hand
Step-Up Or Stepdown At $\frac{1}{2}$, $\frac{3}{4}$, Or
1 mm Increments




WD-312-116
"V" Bend Stop Special Design Makes
1 mm "V" Bends To Shorten Archwires
or Provide A Positive Stop. Excellent
For Placing Stops In Nitium Wires.
For Round Or Rectangular Wires Up To 0.25"




How
WD-312-117
Angle How Plier Efficient Utility Plier
With Serrated Beaks For A Firm Grip
45° Angled Tips To Access
Hard To Reach Areas



WD-312-118
Weingart Plier, Serrated Tips,
Hold Wires Firmly At A Convenient Working Angle
Tapered Beaks Fit Easily
Between Brackets & Hard To Reach
Areas. Our Most Popular Utility Plier



WD-312-119
Multi-Purpose Pliers, Bird Beak,
De La Rosa & Cutting In One
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.46 mm / 180° x 0.56 mm / 0.22"



WD-312-120
Cinch Back / Tie Back
Distal End Bending Pliers
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"




Hammer Head
WD-313-121
Distal End Bending Pliers
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



Nance
WD-313-122
Loop Forming Pliers
0.36 mm / 0.14" To 0.71 mm / 0.28"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



WD-313-123
Detailing Step Pliers / Bayonet
Bending Pliers - 1 mm
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"




WD-313-124
Stop Pliers V-Bebed 1 mm
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.43 mm / 0.17" x 0.56 mm / 0.28"



WD-313-125
Lingua / Palata Arch Forming Plier
0.76 mm / 0.30" x 0.90 mm / 36"




WD-313-126
Torquing Pliers Set



WD-314-127
Tweed / Angel
Rectangular Arch Forming
& Torquing Pliers
0.30 mm / 0.12" To 0.71 mm / 0.28"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



WD-314-128
Surgical Ball Hook
Crimping Pliers, Angled



WD-314-129
Surgical Ball Hook
Crimping Pliers, Straight




Admas
WD-314-130
Wire Bending Pliers
0.30 mm / 0.14" To 0.81 mm / 0.32"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"




Admas
WD-314-131
Wire Bending Pliers
0.30 mm / 0.14" To 0.81 mm / 0.32"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



WD-314-132
Universal Small
Wire Bending & Cutting Pliers
Bend 0.9 mm / 0.36"
Cuts 0.7 mm / 0.28"



WD-315-133
Universal Dolphin Style
Wire Bending & Cutting Pliers
Bend 0.9 mm / 0.36"
Cuts 0.7 mm / 0.28"



WD-315-134
Universal Taper Style
Wire Bending & Cutting Pliers
Bend 0.9 mm / 0.36"
Cuts 0.7 mm / 0.28"



Waldsachs
WD-315-135
Wire Bending & Cutting Pliers
Bend 0.9 mm / 0.36"
Cuts 0.7 mm / 0.28"



WD-315-136
Wire Bending & Cutting
0.70 mm / 0.28" To 1.00 mm / 0.40"
TC




WD-315-137
Hard Wire Cutter
0.70 mm / 0.28" To 1.00 mm / 0.40"
TC



WD-315-138
Hard Wire Cutter
0.60 mm / 0.24" To 1.80 mm / 0.32"
TC




WD-316-139
Hard Wire Cutters
0.60 mm / 0.24" To 1.80 mm / 0.32"
TC




WD-316-140
Hard Wire Cutters
0.87 mm / 0.28" To 1.1 mm / 0.45"




WD-316-141
Hard Wire Cutters
0.87 mm / 0.28" To 1.1 mm / 0.45"



WD-316-142
Band Removing Pliers




WD-316-143
Bracket Removing Pliers
Angled



WD-316-144
Bracket Removing Pliers
Straight



WD-317-145
Bracket Removing Pliers
With Plastic Tip



WD-317-146
Adhesive Removing Pliers



Weingart-Mini
WD-317-147
Utility Pliers
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"




Weingart
WD-317-148
Utility Pliers
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"




How
WD-317-149
Utility Pliers, Straight
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"




WD-317-150
Utility Pliers, Angled
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" X 0.71 mm / 0.28"



WD-318-151
Wire Bending & Loop Forming Pliers
Angled (Bird Beak)
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"




WD-318-152
Light Wire Bending & Loop
Forming Pliers
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.41 mm / 0.16" x 0.56 mm / 0.22"




WD-318-153
Light Wire Bending Pliers
Fitted With Hard Metal Tips Cutter
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.46 mm / 0.18" x 0.56 mm / 0.22"




WD-318-154
Light Wire Bending &
Loop Forming Pliers With Cutter
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.46 mm / 0.18" x 0.56 mm / 0.22"



WD-318-155
Light Wire Bending &
Loop Forming Pliers, Grooved Tip
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.46 mm / 0.18" x 0.56 mm / 0.22"



Tweed
WD-318-156
Loop Forming Pliers With Cutter
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.53 mm / 0.21" x 0.64 mm / 0.25"




Jarabak
WD-319-157
Light Wire Bending Pliers
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.46 mm / 0.18" x 0.56 mm / 0.22"




Aderer
WD-319-158
3 Jaw Wire Bending Pliers
0.30 mm / 0.18" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



WD-319-159
Ligature Forming Pliers
25 mm / 0.10" To 38 mm / 0.15"



Hollow-Chop
WD-319-160
Archwire Forming Pliers
Concave & Convex Tips
0.30 mm / 0.12" To 0.71 mm / 0.28"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"




WD-319-161
Mini Head Distal End Cutter With
Safety Hold, Long Handles Fitted With
Un-Breakable, Corrosion Resistant &
Durable Stainless Steel Tips
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.53 mm / 0.21" x 0.64 mm / 0.25"
0.38 mm / 0.155" To 0.53 mm / 0.215"
Max 0.53 mm / 0.21" x 0.64 mm / 0.25"




WD-319-162
Hard Wire Cutter, 15° Angled
Fitted With Un-Breakable, Corrosion
Resistant & Durable Stainless Steel Tips
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"

Pliers For Orthodontics & Prosthetics | Alicates Para Ortodoncia Y Prótesis

Zangen Für Kieferorthopädie Und Prothetik



WD-320-163
Soft Wire & Ligature Cutter, Mini
Fitted With Un-Breakable, Corrosion
Resistant & Durable Stainless Steel Tips
25 mm / 0.10" To 38 mm / 0.15"



WD-320-164
Band Removing Pliers



WD-320-165
Band Removing Pliers With Plastic Tip



Weingart
WD-320-166
Utility Pliers
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



WD-320-167
Standard Distal End Cutter With
Safety Hold. Fitted With
Un-Breakable, Corrosion Resistant &
Durable Stainless Steel Tips
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



WD-320-168
Mini Head Distal End Cutter With
Safety Hold. Fitted With Un-Breakable,
Corrosion Resistant &
Durable Stainless Steel Tips
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.53 mm / 0.18" x 0.64 mm / 0.25"
0.38 mm / 0.155" To 0.53 mm / 0.215"
Max 0.53 mm / 0.21" x 0.64 mm / 0.25"



WD-321-169
Distal End Cutter With Safety Hold
Mini Head, Long Handles
Fitted With Un-Breakable, Corrosion Resistant &
Durable Stainless Steel Tips
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.46 mm / 0.18" x 0.64 mm / 0.25"
0.38 mm / 0.155" To 0.53 mm / 0.215"
Max 0.53 mm / 0.21" x 0.64 mm / 0.25"



WD-321-170
Flush-Cut Distal End Cutter With Safety Hold
Fitted With Un-Breakable, Corrosion Resistant &
Durable Stainless Steel Tips
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.46 mm / x 0.64 mm / 0.25"



WD-321-171
Soft Wire & Ligature Cutter, Mini
Fitted With Un-Breakable, Corrosion Resistant &
Durable Stainless Steel Tips
25 mm / 0.10" To 38 mm / 0.15"



WD-321-172
Soft Wire & Ligature Cutter, Standard
Fitted With Un-Breakable, Corrosion Resistant &
Durable Stainless Steel Tips
25 mm / 0.10" To 38 mm / 0.15"



WD-321-173
Soft Wire & Ligature Cutter, 15° Angled
Fitted With Un-Breakable, Corrosion Resistant &
Durable Stainless Steel Tips
25 mm / 0.10" To 38 mm / 0.15"




WD-321-174
Soft Wire & Ligature Cutter, 45° Angled
Fitted With Un-Breakable, Corrosion Resistant &
Durable Stainless Steel Tips
25 mm / 0.10" To 38 mm / 0.15"



WD-322-175
Hard Wire Cutter, Straight
Fitted With Un-Breakable, Corrosion
Resistant & Durable Stainless Steel Tips
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



WD-322-176
Hard Wire Cutter, 15° Angled
Fitted With Un-Breakable, Corrosion
Resistant & Durable Stainless Steel Tips
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" X 0.71 mm / 0.28"



WD-322-177
Hard Wire Cutter, 15° Angled
Non-Autoclavable
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



Weingart Narrow
WD-322-178
Slim & Sleek Design With Long Handles
For Better Access To Posterior Area,
Also Suitable For Lingual Technique



Micro
WD-322-179
Ligature Cutters
23 mm / 0.09" To 30 mm / 0.12"




Mini
WD-322-180
Ligature Cutters
23 mm / 0.10" To 38 mm / 0.15"



WD-323-181
Wire Bending & Loop Forming Pliers, Angled
0.36 mm / 0.14" To 0.71 mm / 0.28"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



WD-323-182
Wire Bending & Loop Forming Pliers, Angled
0.36 mm / 0.14" To 0.50 mm / 0.20"
Max 0.41 mm / 0.16" x 0.56 mm / 0.22"



How
WD-323-183
Utility Pliers Standard, Straight
4 mm Wide Tips
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



How
WD-323-184
Utility Pliers Standard, Angled
4 mm Wide Tips
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



WD-323-185
Mini
2.4 mm Wide Tips, Straight
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



How
WD-323-186
Utility Pliers Mini
2.4 mm Wide Tips, Angled
0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



Tweed
WD-324-187
 Loop Forming Pliers
 0.36 mm / 0.14" To 0.64 mm / 0.25"
 Max 0.53 mm / 0.21" x 0.64 mm / 0.25"



Tweed
WD-324-188
 Loop Forming Pliers
 0.36 mm / 0.21" To 0.50 mm / 0.20"
 Max 0.46 mm / 0.18" x 0.56 mm / 0.22"



Nance
WD-324-189
 Loop Forming Pliers
 0.36 mm / 0.14" To 0.71 mm / 0.28"
 Max 0.56 mm / 0.22" x 0.71 mm / 0.28"




De La Rose
WD-324-190
 Arch Forming & Contouring Pliers
 0.36 mm / 0.12" To 0.50 mm / 0.20"
 Max 0.53 mm / 0.21" x 0.64 mm / 0.25"




Jarabak
WD-324-191
 Light Wire Bending Pliers
 0.36 mm / 0.18" To 0.50 mm / 0.22"
 Max 0.46 mm / 0.12" x 0.56 mm / 0.20"




WD-324-192
 Light Wire Bending & Loop Forming Pliers
 0.36 mm / 0.12" To 0.50 mm / 0.20"
 Max 0.41 mm / 0.16" x 0.56 mm / 0.22"



Aderer
WD-325-193
3-Jaw Wire Bending Pliers
0.36 mm / 0.14" To 0.76 mm / 0.30"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



WD-325-194
Lingual / Palatal Arch Forming Pliers
0.76 mm / 0.091" To 0.30 mm / 0.036"



WD-325-195
Slim 0.64 mm (0.025") Wide Tips
Arch Forming & Torquing Pliers
0.36 mm / 0.012" To 0.50 mm / 0.020"
Max 0.41 mm / 0.016" x 0.56 mm / 0.022"



WD-325-196
Angel 1.27 mm / 0.05" Wide Tips
Arch Forming & Torque Pliers
0.36 mm / 0.012" To 0.56 mm / 0.020"
Max 0.53 mm / 0.021" x 0.64 mm / 0.025"




WD-325-197
Tweed / Angel 1.78 mm / 0.07" Wide Tips
Arch Forming & Torque Pliers
0.36 mm / 0.012" To 0.71 mm / 0.028"
Max 0.56 mm / 0.022" x 0.71 mm / 0.028"



WD-325-198
Bracket Removing Pliers

Pliers For Orthodontics & Prosthetics | Alicates Para Ortodoncia Y Prótesis
Zangen Für Kieferorthopädie Und Prothetik



WD-326-199
Bracket Removing Pliers



WD-326-200
Band Removing Pliers



WD-326-201
Band Removing Curved, Left & Right



WD-326-202
Left



WD-326-203
Right, Pliers For Bands



Johnson
WD-326-204
Contouring Pliers



Reynolds
WD-327-205
Contouring Pliers





WD-327-206
For Band Seating



WD-327-207
For Incisors With Convex Beaks



WD-327-208
For Bicuspsids & Molars
With Convex Beaks



WD-327-209
Straight
Pliers Ball Hook Crimping



WD-327-210
Angled
Pliers Ball Hook Crimping



WD-327-211
Micro Distal End Cutters
With Safety Hold
Max 0.30 mm / 0.12" To 0.50 mm / 0.20"
Max 0.46 mm / 0.18" To 0.64 mm / 0.25"
Max 0.38 mm / 0.155" To 0.53 mm / 0.215"
Max 0.53 mm / 0.21" To 0.64 mm / 0.25"
TC


WD-328-212

Mini

Max 0.30 mm / 0.12" To 0.50 mm / 0.20"

Max 0.46 mm / 0.18" To 0.64 mm / 0.25"

Max 0.38 mm / 0.155" To 0.53 mm / 0.215"

Max 0.53 mm / 0.21" To 0.64 mm / 0.25"

TC


WD-328-213

Mini With Long Handles

Max 0.30 mm / 0.12" To 0.50 mm / 0.20"

Max 0.46 mm / 0.18" To 0.64 mm / 0.25"

Max 0.38 mm / 0.155" To 0.53 mm / 0.215"

Max 0.53 mm / 0.21" To 0.64 mm / 0.25"

TC

Distal End Cutters With Safety Hold


WD-328-214

Mini With Extra Long Off-Set Handles

Max 0.30 mm / 0.12" To 0.50 mm / 0.20"

Max 0.46 mm / 0.18" To 0.64 mm / 0.25"

Max 0.38 mm / 0.155" To 0.53 mm / 0.215"

Max 0.53 mm / 0.21" To 0.64 mm / 0.25"

TC


WD-328-215

Distal End Bending Pliers

Max 0.30 mm / 0.12" To 0.50 mm / 0.20"

Max 0.56 mm / 0.22" To 0.71 mm / 0.28"


Weingart Mini
WD-328-216

Slim Design For Better Accessibility
In Posterior & Lingual Areas


Weingart Narrow
WD-328-217

Narrow Long Tips For Better
Access To Hard To Reach Areas



Weingrat
WD-329-218
Utility Pliers
Standard



Aderer
WD-329-219
Small
0.36 mm / 0.014" To 0.50 mm / 0.020"
Max 0.53 mm / 0.021" x 0.64 mm / 0.025"



Aderer
WD-329-220
Medium
0.36 mm / 0.014" To 0.71 mm / 0.028"
Max 0.56 mm / 0.022" x 0.71 mm / 0.028"



Schwarz
WD-329-221
Small
Max 0.50 mm / 0.020"



Schwarz
WD-329-222
Medium
Max 0.70 mm / 0.028"
Wire Bending & Loop Forming Pliers



WD-329-223
Angle
0.36 mm / 0.014" To 0.71 mm / 0.028"
Max 0.56 mm / 0.022" x 0.71 mm / 0.028"

Pliers For Orthodontics & Prosthetics | Alicates Para Ortodoncia Y Prótesis

Zangen Für Kieferorthopädie Und Prothetik




WD-330-224
Angle
0.36 mm / 0.14" To 0.71 mm / 0.28"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



Adams
WD-330-225
0.36 mm / 0.14" To 0.71 mm / 0.28"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



Adams
WD-330-226
Small
0.36 mm / 0.14" To 0.71 mm / 0.28"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



WD-330-227
0.36 mm / 0.14" To 0.71 mm / 0.28"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"




Adams
WD-330-228
Large
0.36 mm / 0.14" To 0.81 mm / 0.32"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"
Wire Bending & Loop Forming Pliers



Young
WD-330-229
0.36 mm / 0.14" To 0.71 mm / 0.28"
Max 0.56 mm / 0.22" x 0.71 mm / 0.28"



Adams
WD-331-230
Large
0.36 mm / 0.014" To 0.81 mm / 0.032"
Max 0.56 mm / 0.022" x 0.71 mm / 0.028"




Aderer
WD-331-231
Large
0.50 mm / 0.020" To 0.91 mm / 0.036"
Max 0.56 mm / 0.022" x 0.71 mm / 0.028"



Schwarz Graduated
WD-331-232
Max 0.70 mm / 0.028"




Schwarz Graduated
WD-331-233
Large
Max 0.70 mm / 0.028"



Peeso
WD-331-234
Max 0.70 mm / 0.028"
Wire Bending & Loop Forming Pliers



Peeso
WD-331-235
Max 0.70 mm / 0.028"



WD-331-236
0.36 mm / 0.014" To 0.71 mm / 0.028"
Max 0.56 mm / 0.022" x 0.71 mm / 0.028"


WD-332-001

WD-332-002


WD-332-003

WD-332-004

WD-332-005

WD-332-006

**Methot Jar
WD-332-007**
50 mm Ø x 75 mm
Without Lid

Kidney Tray
WD-332-008 150 mm
WD-332-009 200 mm
WD-332-010 250 mm
WD-332-011 300 mm

Mayo Tray Perforate
WD-332-012 150 x 100 x 25 mm
WD-332-013 200 x 100 x 25 mm
WD-332-014 220 x 120 x 25 mm
WD-332-015 240 x 140 x 25 mm
WD-332-016 200 x 150 x 25 mm
WD-332-017 250 x 200 x 25 mm
WD-332-018 300 x 250 x 25 mm
WD-332-019 300 x 200 x 25 mm
WD-332-020 350 x 250 x 25 mm
WD-332-021 400 x 300 x 25 mm
WD-332-022 450 x 300 x 25 mm
WD-332-023 450 x 350 x 25 mm

Instruments Tray
WD-332-024 150 x 100 x 50 mm
WD-332-025 200 x 100 x 50 mm
WD-332-026 220 x 120 x 50 mm
WD-332-027 240 x 140 x 50 mm
WD-332-028 200 x 150 x 50 mm
WD-332-029 250 x 200 x 50 mm
WD-332-030 300 x 250 x 50 mm
WD-332-031 300 x 200 x 50 mm
WD-332-032 350 x 250 x 50 mm
WD-332-033 400 x 300 x 50 mm
WD-332-034 450 x 300 x 50 mm
WD-332-035 450 x 350 x 50 mm

Scalars Tray
WD-332-036
200 x 100 x 25 mm

Instrument Boxes with Lid
WD-332-037 170 x 70 x 35 mm
WD-332-038 180 x 80 x 40 mm
WD-332-039 200 x 100 x 40 mm
WD-332-040 200 x 100 x 50 mm
WD-332-041 200 x 100 x 60 mm
WD-332-042 220 x 120 x 60 mm
WD-332-043 250 x 120 x 60 mm
WD-332-044 300 x 150 x 60 mm
WD-332-045 300 x 200 x 80 mm
WD-332-046 400 x 160 x 80 mm
WD-332-047 400 x 200 x 100 mm
WD-332-048 500 x 200 x 120 mm


Direct Bonding Bracket Holders - Self Holding



WD-333-049
13.5 cm / 5¼"



WD-333-050
13.5 cm / 5¼"



WD-333-051
14 cm / 5½"

Band Seating & Removing Instrument




WD-334-052-SD
WD-334-052-HH

WD-334-053-SD
WD-334-053-HH

WD-334-054-SD
WD-334-054-HH

WD-334-055-SD
WD-334-055-HH



WD-335-001



WD-335-002



WD-335-003
Folding Type, Cork Screw Point



WD-335-004



WD-335-005



WD-335-006



Euro
WD-336-007
Breech Loading Type
1.8 ml



WD-336-008
T-Handle Type
1.8 ml



WD-336-009
Euro Type
1.8 ml



WD-336-010
3-Ring Type
1.8 ml



WD-336-011
Euro Type
1.8 ml
Manual Aspitating Syringes



WD-336-012
Type A
1.8 ml



WD-337-013
Type A
1.8 ml



WD-337-014
Type C
1.8 ml



WD-337-015
Euro Type
1.8 ml



WD-337-016
3-Ring Type
1.8 ml



WD-337-017
2-Ring Type
1.8 ml



Euro
WD-337-018
Breech Loading Type
1.8 ml



Euro
WD-338-019
Breech Loading Type
1.8 ml



WD-338-020
Disc System
Self Aspiring Syringes
1.8 ml



WD-338-021
Disc System
1.8 ml



WD-338-022
Disc System
1.8 ml



WD-338-023
Projection System
1.8 ml



WD-338-024
Projection System
1.8 ml



WD-339-025
Type A



WD-339-026
Qyatro Type
1.8 ml



WD-339-027
Trio Type
1.8 ml



WD-339-028
Type A
1.8 ml



WD-339-029
Type C
1.8 ml



Euro
WD-339-030
Breech Loading Type
1.8 ml




Pistol-Ject
WD-340-031
1.8 ml



WD-340-032
Plastic Protection Tube

WD-340-033
Anesthesia Ayringe
In Black Colour



WD-340-034
Swedish Type
1.8 ml



Falcon-Ject
WD-340-035
Angled
1.8 ml



Falcon-Ject
WD-340-036
Straight
1.8 ml

Index English

A			
Abscess Knives	98		
Acorn Filling Instruments	184		
Adams Plier	302, 309, 314, 330, 331		
Aderer Clamps	309		
Aderer Facebow loop Forming Pliers	308		
Aderer Jaw Wire Bending Pliers	319		
Aderer Pliers	329, 331		
Aderer Three Jaw Wire Bending Pliers	305, 318		
Aderer Wire Bending Pliers	304, 305, 325		
Adhesive Removing Pliers	295, 317		
Admas Wire Bending Pliers	314		
Adson Needle Holders	102		
Adson Tissue Forceps	135		
Adson-Brown Tissue Forceps	135		
Ainsworth Rubber Dam Punch Forceps	158		
Alcohol Burners	235		
Alm Retractors	130		
Aluminium Endodontic Boxes	233, 234		
Amalgam Carriers	166, 167		
Amalgam Carvers	177 - 180		
Amalgam Pluggers	181 - 183		
Amalgam Well	166		
Angle Bird Beak Plier	304		
Apexo Root Elevators	64, 67, 74		
Apical Root Elevators	64		
Arch Bending Pliers	292		
Arruga Needle Holders	105, 112		
Articulators	263 - 266		
Ash Abscess Knives	98		
Autoclavable Adaptador Saliva Ejectors	229		
Autoclavable Plastic Amalgam Holders	166		
Awls	99		
B			
Back-Action Amalgam Pluggers	181		
Backhaus Towel Clamps	116		
Ball Filling Instruments	173		
Band Removing Pliers	316, 320, 326		
Band Seating & Removing Instrument	334		
Band Seating Pliers	295		
Bane Bone Rongeurs	56		
Barnhart Curettes	204		
Barraquer Micro Scissors	88		
Barraquer Needle Holders	106, 113		
Barry Root Elevators	75		
Barsky Delicate Retractors	122		
Baumgartner Needle Holders	100, 107, 110		
Baynton Needle Holders	105, 112		
Beale Periosteal Instrument	195		
Beale Wax & Modelling Carvers	270		
Bee Bee Stitch & Wire Cutting Scissors	85		
Behrman Amalgam Pluggers	181		
Bein Root Elevators	63, 71, 72		
Beyer Bone Rongeurs	56, 57		
Bird Beak Pliers	292, 293, 302, 309		
Black Amalgam Pluggers	181 - 183		
Black Cavity Preparation Instruments	186		
Black Scalars	214B		
Blumenthal Bone Rongeurs	54		
Bodkin Amalgam Pluggers	182		
Böhler Bone Cutting Forceps	61		
Böhler Bone Rongeurs	57		
Böhm Telescope Crown Pliers	224		
Boley Measuring Instruments	287		
Bonding Bracket Remover Plier	297, 306		
Bone Bracket Remover Plier	297		
Bone Chisels	193		
Bone Curettes	217 - 219		
Bone Cutting Forceps	59 - 61		
Bone Files	215 - 216		
Bone Grafting Instruments	282		
Bone Morselizer	279		
Bone Rongeurs	53 - 58		
Bone Spoons / Packers	281		
Bone Tampers	280		
Bowman Antrum Probes	131		
Boxes For Instruments	235		
Boynton Needle Holders	105, 112		
Bracket Holders	333		
Bracket Removing Pliers	297, 316, 325, 326		
Bradley Root Elevators	66		
Brewer Rubber Dam Instruments	158		
Brünings Tongue Depressors	118		
Brushes	273		
Buchwald Tongue Depressors	118		
Buck Files	198		
Buck Gingivectomy Knives	201		
Buckley Bone Chisels & Gouges	191		
Burners	235		
Burnishers	169, 173, 184		
C			
Calipers	286 - 289		
Castro Micro Scissors	88		
Castroviejo (Barraquer) Needle Holders	106, 113		
Castroviejo Measuring Instruments	288		
Castroviejo Needle Holders	106, 113		
Castroviejo-Epker Measuring Instruments	289		
Cattoni Double - Ended Explorers	146, 147		
Cavity Preparation Instruments	186		
Cawood Minnesota Tongue Depressors	118		
Cement Spatulas	150		
Ceramic Sculpturing Instruments	274		
CF5A Filling Instruments	178		
Chisels	193		
Chompret Syndesmotomes	79		
Cleveland Bone Cutting Forceps	60		
Cleveland Bone Rongeurs	55		
Cogswell Root Elevators	67		
Coleman Root Elevators	62		
College Cotton & Dressing Pliers	141		
Collier Needle Holders	100, 107		
Collin Tongue Forceps	119		
Columbia Tongue Depressors	118		
Columbia University Curettes & Scalars	203		
Combination Explorer - Probes	149		
Compactors	280, 282		
Cone Socket - Mouth Mirror Handles	142		
Conservative Treatment Instruments Set	143		
Contouring Plier	310		
Cottle-Kazanjian Bone Cutting Forceps	61		
Cotton & Dressing Pliers	138 - 141		
Coupland Root Elevators	69, 70		
Crane Root Elevators	70		
Crane-Kaplan Curettes	209		
Crane-Kaplan Gingivectomy Knives	202		
Crile Haemostatic Forceps	115		

Crile-Wood Needle Holders	102, 109
Crimping Pliers	314
Crown Instruments	220 - 225
Crown Removers	220
Crown Spreaders	220, 223
Crucible Tongs	290
Cryer Extracting Forceps	32, 33, 51
Cryer Root Elevators	68, 71, 74
Curettes	203 - 210, 217 - 219
Curtis Root Elevators	66
D	
De La Rosa Pliers	305
De La Rose Arch Forming & Contouring Pliers	324
De Marco Curettes	206
Dean Gum Scissors	90, 94
Delicate Retractors	122 - 125
Denhart Mouth Gags	119
Dental Syringes	335 - 340
Derf (Wright) Needle Holders	100, 107
Didner Matrix Retainers	226
Dieffenbach Operating Knives	97
Distal End Cutter	295, 297, 300
Double Rounded Jaw Wire Bending Plier	301
Double-Ended Explorers	146
Doyen-Collin Mouth Gags	120
Doyen-Jansen Mouth Gags	120
Dressing Forceps	117, 137 - 141
E	
Eiselberg-Mathieu Needle Holders	103
Elastic Separating Pliers	159
Elevators	62 - 76
Endodontic Boxes	233, 234
Endodontic Instruments	154 - 157
Euro Dental Syringes	336, 337, 338, 339
Excavators	187 - 190
Explorers	144 - 146
Extracting Forceps	1 - 52
Extracting Forceps (American Pattern)	27 - 34
Extracting Forceps For Children	35 - 38
F	
Face Bow Bending & Cutting Pliers	299
Fahnenstock Waxing Instruments	275
Falcon Beak Contouring Plier	310
Falcon-Ject Dental Syringes	340
Farabeuf Retractors	125
Farrell Cotton Applicator	131
Ferguson-Ackland Mouth Gags	120
Files	215 - 216
Filling Instruments	168 - 184
Fine Operating Scissors	89
Flagg Cotton & Dressing Pliers	140
Flexible Spatulas	279
Flohr Root Elevators	64, 71, 72
Force Module Separating Pliers	293
Frahm Amalgam Carvers	178
Frahm Wax & Modelling Carvers	267
Frazier Suction Tubes	228
Freer Periodontia Instrument	196
Fret Saws	285
Friedman Root Elevators	70
Friedmann Bone Rongeurs	53
Furrer Copper Ring Removing Pliers	222
G	

Gardner Periodontia Instruments	192
Gerzog Mallet	194
Gillies Delicate Retractors	122
Gillies Needle Holders	103, 110
Glick Endo Pluggers & Excavators	177
Goldman-Fox (Crane-Kaplan) Dressing Pliers	138, 194
Goldman-Fox / Williams Probes	147 - 149
Goldman-Fox Curettes	207
Goldman-Fox Gingivectomy Knives	201
Goldman-Fox Gum Scissors	89, 90, 93
Goldman-Fox Periodontia Instrument	194, 195
Goldman-Fox Probes	148
Gracey Curettes	205, 206
Graft Material Spoon	282
Graft Packers	281, 282
Gregg Filling Instruments	177
Gritman Plaster Knives	276
Gritman Wax & Modelling Carvers	270
Gross-Maier Dressing Forceps	117
Gum Scissors	89 - 92
Gum Scissors With Tungsten Carbide Inserts	93, 94
Guthrie Delicate Retractors	123
H	
Haemostatic Forceps	114, 115
Halsey Needle Holders	100, 107
Halsted - Mosquito Haemostatic Forceps	114
Hammer Head Distal End Bending Pliers	313
Handles For Crown Instruments	221
Handles For Explorers	144
Handles For Mouth Mirrors	142
Handles For Root Elevators	62A
Harris Extracting Forceps	29
Hartmann Bone Rongeurs	54
Hartmann Haemostatic Forceps	114
Hartmann-Halle Nasal Speculas	131
Heath-Stitch & Wire Cutting Scissors	85
Hegar-Baumgartner Needle Holders	100, 107
Heidbrink Root Elevators	67, 74
Heidbrink Root-Tip Picks	77
Heidemann Separating Spatulas	151 - 153
Heister Tongue Forceps	119
Hemingway Bone Curettes	219
Henahan Extracting Forceps	50
Hirschfeld Curettes	210
Hollenback Filling Instruments	171, 176
Hollenback Wax & Modelling Carvers	267
Hollow Chop Contouring Pliers	292
Hollow Handles For Explorers	144
Hollow-Chop Archwire Forming Pliers	319
Hopkins Periodontia Instruments	196
Horn Mallet	284
How Pliers	294, 295, 305, 312
How Utility Pliers	317, 323
Howard Bone Files	216
Hull Extracting Forceps	32
Hygienist Scalers	214A
Hyllin Root Elevators	62
Hyllin Wax & Modelling Carvers	269
Hyllin Wax & Modelling Carvers	269, 270
I	
Implantology Instruments	280 - 283
Impression Trays	237 - 262
Impression Trays Plastic - Disposable	261

Indiana University Curettes	206	Loop Forming Plier	306
Instruments Boxes	233 - 234, 332	Lorna (Edna) Towel Clamps	116
Instruments Set For Inplant Dentistry	143	Lucas Bone Curettes	219
Instruments Set For Oral Surgery	80	Luer Bone Rongeurs	55
Instruments Trays	230 - 232, 332	Luks Root-Canal Pluggers	155
Interproximal Amalgam Carvers	178	Luniatschek Operative Instruments	185
Iris Gum Scissors	89, 93	M	
Ivory Matrix Retainers	226	Mac Millan Root Elevators	70
Iwanson Measuring Instruments	288	Machatt Delicate Retractors	123
J		Maehlum Modeling Carver	269
Jameson Measuring Instruments	286	Magnm Orthodontics Pliers	299
Jansen Bone Rongeurs	54	Makellops Plier	304
Jaquette Curettes	210	Mallets	194, 284
Jarabak Light Wire Bending Pliers	319, 324	Marquardt Bone Rongeurs	58
Jarabak Orthodontic Pliers	292	Marquette Filling Instruments	171
Jennings Mouth Gags	121	Martin Ergo Extracting Forceps	51
Johnson Pliers	303, 309, 326	Martini Bone Curettes	217
Jones Towel Clamps	116	Mathieu Needle Holders	104, 111
K		Matrix Bands	227
Kells Extracting Forceps	31	Matrix Retainers	226
Kelly Gum Scissors	91, 94	Mayfield Bone Rongeurs	57
Kelly Haemostatic Forceps	114	Mayo Surgical Scissors	86
Kerrison Bone Rongeurs	58	Mayo Trays	332
Kidney Trays	332	Mayo-Hegar Needle Holders	102, 109
Kilner Delicate Retractors	125	Mc Call Curettes	208
Kirkland Delicate Retractors	125	Mc Call Scalers	214B
Kirkland Gingivectomy Knives	200	Mc Intosh Laryngoscopes	132, 133
Klein Extracting Forceps	47, 48	Mckellops Loop Forming Pliers	308
Kleinert-Kutz Bone Rongeurs	53	Mead Bone Rongeurs	53
Kliner Reverdin Needles & Awls	99	Mead Mallets	194
Knapp Delicate Retractors	123	Mead Periodontia Instruments	195
Knives	97, 98, 275, 276	Measuring Instruments	286 - 289
Kocher Haemostatic Forceps	115	Meriam Cotton & Dressing Pliers	140
Kocher Retractors	127	Merrit Periodontal Pocket Probes	148
König Periodontia Instrument	197	Mershon Rubber Dam Instruments	159
Kopp Root Elevators	62B	Methot Jar	332
L		Metzenbaum Surgical Scissors	87
La Grange Gum Scissors	90	Micro Ligature Cutters	322
Laboratory Tweezers	291	Micro Scissors	88
Ladmore Filling Instruments	173	Micro-Adson Tissue Pliers	135
Langenbeck Nasal Specula	131	Micro-Vascular Needle Holders	101, 108
Langenbeck Retractors	128	Middeldorp Retractors	128
Langer Curettes	204	Miller Bone Files	215
Laryngoscopes	132, 133	Miller Crown Instruments	225
Lawrence-Read Extracting Forceps	40, 47	Miller Paper Forceps	225
Le Cron Wax & Modelling Carvers	269, 270	Miller-Apexo Root Elevators	65, 73
Lecluse Root Elevators	75	Mini Ligature Cutters	322
Lempert Bone Rongeurs	56	Mini Pin & Ligature Cutter	296
Lessmann Waxing Instruments	275	Mini-Friedman Rongeur	53, 194
Lichtenberg Needle Holders	105, 112	Mini-Ryder Needle Holders	101, 108
Ligature Forming Pliers	319	Modelling Carvers	267 - 272
Ligature Placer	297	Molt Bone Curettes	217
Lingual Arch Forming Pilers	292	Molt Mouth Gags	120
Lip Retractors	118	Molt Periodontia Instrument	195
Lister Bandage Scissors	87	Montfort Crown Slitter	222
Liston Bone Cutting Forceps	59	Mosquito Haemostatic Forceps	114
Littauer Bone Cutting Forceps	59	Mouth Gags	119 - 121
Littauer Stitch & Wire Cutting Scissors	85	Mouth Mirrors	142
Locking Cotton & Dressing Pliers	138	N	
Locklin Gum Scissors	91, 94	Nabers Probes	147
London College Dressing Pliers	139	Nance Loop Closing Plier	300
London Hospital Root Elevators	62	Nance Loop Forming Pliers	313, 324
Loop Closing Plier	298	Nance Plier	303, 308

Nasal Speculas	131
Needle Holders	100 - 113
Needle Holders With Tungsten Carbide Inserts	107 - 113
Neumann Gum Scissors	90
Nevius Extracting Forceps	31
Niro Bone Cutting Forceps	59
Noyes (Castroviejo) Gum Scissors	92
O	
Ochsenbein Bone Chisels & Gouges	191
Ochsenbein Periodontia Instruments	198
Olsen-Hegar Needle Holders	103, 109
Omega Loop Forming Pliers	294
Operating Knives	97
Operating Scissors	81 - 84
Operative Instruments	185
Orban Gingivectomy Knives	200
Orban Loop Operative Instruments	185
Ortho Boxes	332
Orthodontic Instruments	333, 334
Oval Spoon Excavators	188, 189
P	
P.K Thomas Wax & Modelling Carvers	271, 272
Palmer Rubber Dam Instruments	158
Paper Forceps	225
Parmly Extracting Forceps	29, 50
Partsch Bone Chisels & Gouges	193
Partsch Mallet	194
Partsch Periodontia Instrument	197
Pean Haemostatic Forceps	115
Peeso Pliers	331
Periodontal Instruments Set	211, 212
Periodontal Pocket Probes	147 - 148
Periodontia Instruments	191 - 211
Perry Cotton & Dressing Pliers	139
Physick Extracting Forceps	29
Pin & Ligature Cutting Pliers	296, 297
Pistol-Ject Syringe	340
Plaster Cutting Pliers	279
Plaster Knives	276
Pliers	292 - 331
Pliers For Removing Broken Broaches	157
Polokoff Bone Files	215
Porcelain Sculpturing Brushes	273
Posterior Band Removing Plier	293, 307
Potts Root Elevators	62
Potts-Smith Dressing Pliers	138
Prichard Periodontia Instrument	195
Probes	147 - 149
Q	
Quinby Gum Scissors	90
R	
Ragnell Delicate Retractors	122
Ramus Retractors	127
Rectangular Arch Bending Pliers	293
Retractors	122 - 130
Reverdin Needles & Awls	99
Reynolds Contouring Pliers	301, 327
Rhodes Back-Action Periodontia Instruments	198
Riveting Mallet	284
Roach Wax & Modelling Carvers	270
Rongeurs	53 - 58
Root Elevators	62 - 76
Root Elevators-Handles	62A

Root-Canal Excavators & Explorers	156
Root-Tip Picks	77
Roser-Koenig Mouth Gags	119
Round & Concave Plier	301
Round Nose Wire Bending Plier	301
Rubber Dam Clamp Trays	165
Rubber Dam Clamps	160 - 163
Rubber Dam Instruments	158, 159, 166
Rubber Dam Instruments Kit	164
Ruskin Bone Rongeurs	58
Ruskin-Liston Bone Cutting Forceps	61
S	
Saliva Ejectors	228, 229
Sanders Gingivectomy Knives	202
Scaler Trays	230 - 232, 332
Scalers	213 - 214B
Scalpel Handles - Scalpel Blades	95, 96
Schaedle Towel Clamps	116
Schluger Files	198
Schmeckebier Root Elevators	64, 72
Schroeder Abscess Knives	98
Schwarz Graduated Pliers	331
Schwarz Pliers	308, 329
Scissors	81 - 94
Seldin Retractor	197
Seldin Root Elevators	65 - 66, 67, 73
Semken-Taylor Forceps	134
Senn-Miller Delicate Retractors	122
Spatulas Flexible Complete Set	279
Spatulas For Plaster & Alginate	277 - 278
Spencer Stitch & Wire Cutting Scissors	85
Spencer-Wells Haemostatic Forceps	115
Spirit Lamps	235
Spreader Endodontic Instruments	154
Stellbrink Bone Rongeurs	56
Sterilizing Forceps	236
Sternberg Tongue Depressors	118
Stevens Dissecting Scissors	86
Stieglitz Forceps & Pliers	157
Stitch & Wire Cutting Scissors	85
Stocke Rubber Dam Instruments	158
Stout Root Elevators	66
Straight Delicate Retractors	122
Suction Tubes	228
Sugarman Files	198
Sullivan Gum Scissors	91
Surgical Scissors	86, 87
Suture Pliers	134, 135
Swedish Needle Holders	101, 108
Syndesmotomes	78, 79
Syringes	335 - 340
T	
Tampers	282
Tanner Filling Instruments	171
Taylor Curettes	209
TG Chisels	199
Three Jaw Pliers	294
Tissue Forceps	136 - 137
Tissue Pliers - Suture Pliers	134 - 137
Tissue Punches	283
Toennis Needle Holders	103, 110
Tofflemire Matrix Bands	227
Tofflemire Matrix Retainers	226

Tomati Gum Guillotine Forceps	194	Woodson Filling Instruments	176
Tomes Extracting Forceps	31	Woodward Root Elevators	69, 70
Tongue Depressors	118	Y	
Tongue Forceps	119	Young Arch	159
Torquing Pliers Set	313	Young Plier	303, 330
Towel Clamps	116	Young Tongue Forceps	119
Towner Jaquette Scalers	213	Younger-Good Curettes	209
Towner Scalers	213	Z	
Townley Measuring Instruments	286	Zaufal-Jansen Bone Rongeurs	57
Trays For Instruments	230 - 234	Zehle Wax & Modelling Carvers	270
Tweed / Angel Rectangular Pliers	314	Zerfing Chisels	203
Tweed Loop Forming Pliers	324	Zerfing Scalers	214A
Tweed Loop Forming With Cutter Plier	318		
Tweed Pliers	292, 304, 306		
Tying Pliers	293		
U			
U.S.A. Model Bone Curettes	218		
Universal Crown Instruments	223		
Universal Pin Raughing, Cutting & Bending Plier	300, 301		
Universal Stitch & Wire Cutting Scissors	85		
Universal Surgical Scissors	87		
University Of Michigan Probes	148, 149		
University Of Michigan Root Elevators	70		
University Of South California Gingivectomy Knives	202		
USC Amalgam Carvers	178		
USC Towner Gingivectomy Knives	202		
USC Wax & Modelling Carvers	267		
V			
Vehe Wax & Modelling Carvers	270		
Vernier Calipers	286 - 289		
Vertical Condenser Endodontic Instruments	155		
W			
Wahl Bone Files	216		
Wakefield Bone Chisels & Gouges	191		
Waldsachs Pliers	301, 315		
Walter F. Barry Root Elevators	75		
Ward Filling Instruments	171, 176		
Ward Wax & Modelling Carvers	267		
Warwick-James Root-Tip Picks	77		
Wax & Modelling Carvers	267 - 272		
Waxing Instruments	275		
Wedelstaet Cavity Preparation Instruments	186		
Weingar Plier Serated Tips Hold Wires	312		
Weingart Pliers	294, 304, 317, 320, 322, 328		
Weingart-Mini Utily Pliers	317, 329		
Weitlaner - Wullstein Retractors	129		
Weitlaner Retractors	129, 130		
Westcott Burnishers	169, 173		
Westcott Gum Scissors	92		
White Burnishers	184		
Whitehead Mouth Gags	121		
Wiland Amalgam Carvers	179		
Williams Probes	148, 149		
Williger Bone Curettes	217		
Williger Periosteals & Elevators	197		
Williger Retractors	126		
Winter Root Elevators	76		
Wire Bending Pliers	299 - 315, 318, 319, 323 - 325, 329 - 331		
Wire Bird Beak Pliers	292		
Wire Cutters	294, 298 - 300		
Wire Cutting Pliers	295, 296		
Wire Cutting Scissors	85		
Witzel Extracting Forceps	51, 52		

Indice Español

A			
Abaja Lenguas	118		
Abre Bocas	119 - 121		
Abre coronas	220, 223		
Acorn Bruñidores	184		
Adams Alicates Para Ortodoncia Y Prótesis	302, 309, 314, 330, 331		
Aderer Campos Y Alixates Para Ortodoncia Y Prótesis	304, 305, 308, 309, 319, 325, 329, 331		
Adson Pinzas De Disección	135		
Adson Porta-Agujas	102		
Adson-Brown Pinzas De Disección	135		
Ainsworth Instrumentos Para Dique	158		
Alcance De Medición	286 - 289		
Alicates Para Alambres	294, 295, 296, 298 - 300		
Alicates Para Bandes	295		
Alicates Para Brocas	157		
Alicates Para Cortar Yeso	279		
Alicates Para Ligaduras	296		
Alicates Para Ortodoncia Y Prótesis	292 - 331		
Alm Separadores	130		
Angle Bird Alicates Para Ortodoncia Y Prótesis	304		
Apexo Elevadores Para Raíces	64, 67, 74		
Apical Elevadores Para Raíces	64		
Arranca-Coronas	220		
Arruga Porta-Agujas	105, 112		
Articuladores	263 - 266		
Ash Bisturios Para Abscesos	98		
B			
Backhaus Pinzas De Campo	116		
Ball Instrumentos Para Obturaciones	173		
Bandejas	230 - 232, 332		
Bandejas De Acero Inoxidable	230 - 232		
Bandejas Escalador	230 - 232, 332		
Bandejas Para Instrumentos	230 - 234		
Bandejas Porta-Clamp	165		
Bane Pinzas Gubias	56		
Barnhart Curetas	204		
Barraquer Micro Tijeras	88		
Barraquer Porta-Agujas	106, 113		
Barry Elevadores Para Raíces	75		
Barsky Separadores Finos	122		
Baumgartner Porta-Agujas	100, 107, 110		
Baynton Porta-Agujas	105, 112		
Beale Instrumentos Para Periodoncia	195		
Beale Talladores Para Cera Y Para Modelar	270		
Bee Bee Tijeras Para Hilos Y Alambre	85		
Behrman Picos De Amalgama	181		
Bein Elevadores Para Raíces	63, 71, 72		
Beyer Pinzas Gubias	56, 57		
Bird Alicates Para Ortodoncia Y Prótesis	292, 293, 302, 309		
Bisturios Para Abscesos	98		
Bisturios Para Operaciones	97, 98		
Black Instrumentos Para Esmalte	186		
Black Limpiadores De Sarro	214A		
Black Picos De Amalgama	181 - 183		
Blumenthal Pinzas Gubias	54		
Bodkin Picos De Amalgama	182		
Bohler Pinzas Cortantes Para Huesos	61		
Bohler Pinzas Gubias	57		
Böhm Instrumentos Para Coronas	224		
Boley Instrumentos De Medición	287		
Bowman Sondas Para El seno Maxilar	131		
Boynton Porta-Agujas	105, 112		
Brewer Instrumentos Para Dique	158		
Bruñidores	169, 173, 184		
Brünings Abaja Lenguas	118		
Buchwald Abaja Lenguas	118		
Buck Bisturios Para La Gingivectomia	201		
Buck Limas	198		
Buckley Cinceles Y Gubias	191		
C			
Cajas De Endodoncia	233, 234		
Cajas De Ortho	332		
Cajas Para Instrumentos	233, 234, 332		
Calibres	286 - 289		
Castro Micro Tijeras	88		
Castroviejo (Barraquer) Porta-Agujas	106, 113		
Castroviejo Instrumentos De Medición	288		
Castroviejo Porta-Agujas	106, 113		
Castroviejo-Epker Instrumentos De Medición	289		
Cattoni Exploradores Doubles	146, 147		
Cawood Minnesota Abaja Lenguas	118		
CF5A Instrumentos Para Obturaciones	178		
Chompret Sindesmótomos	79		
Cinceles	193		
Cinceles Para Huesos	193		
Cintas De Matrices	227		
Clamps Para Dique	160 - 163		
Cleveland Pinzas Cortantes Para Huesos	60		
Cleveland Pinzas Gubias	55		
Cogswell Elevadores Para Raíces	67		
Coleman Elevadores Para Raíces	62		
College Pinzas Para Algodón Y Curaciones	141		
Collier Porta-Agujas	100, 107		
Collin Abre Bocas Pinzas - Tira Lengua	119		
Columbia Abaja Lenguas	118		
Columbia Instrumentos Para Periodoncia	203		
Cottle-Kazanjian Pinzas Cortantes Para Huesossetzt	61		
Coupland Elevadores Para Raíces	69, 70		
Crane Elevadores Para Raíces	70		
Crane-Kaplan Bisturios Para La Gingivectomia	202		
Crane-Kaplan Curetas	209		
Crile Pinzas Hemostáticas	115		
Crile-Wood Porta-Agujas	102, 109		
Cryer Elevadores Para Raíces	68, 71, 74		
Cryer Fórceps De Extracción (Forma Americana)	32, 33, 51		
Cubetas De Plástico Un Solo Uso	261		
Cubetas Para Impresiones	237 - 262		
Cuchillos	97, 98, 275, 276		
Cuchillos Para Yeso	276		
Curetas	203 - 210, 217 - 219		
Curetas Cortantes	217 - 219		
Curtis Elevadores Para Raíces	66		
D			
De La Rosa Alicates Para Ortodoncia Y Prótesis	305		
De Marco Curetas	206		
Dean Tijeras Para Encías	90, 94		
Denhart Abre Bocas Pinzas - Tira Lengua	119		
Derf (Wright) Porta-Agujas	100, 107		
Didner Porta Matrices	226		
Dieffenbach Bisturios Para Operaciones	97		
Doyen-Collin Abre Bocas	120		
Doyen-Jansen Abre Bocas	120		
E			
Eiselberg-Mathieu Porta-Agujas	103		
Elevadores	62 - 76		

Elevadores Para Fragmentos	77
Elevadores Para Raíces	62 - 76
Escaladores higienist	214A
Espátula Flexible	279
Espátulas Para Cemento	150
Espátulas Para Yeso Y Alginatos	277 - 278
Espéculos Nasaes	131
Espejos De Boca	142
Excavadores	187 - 190
Excavadores Y Exploradores	156
Exploradores	144 - 146
Exploradores Doubles	146
F	
Falcon Alicates Para Ortodoncia Y Prótesis	310
Falcon-Ject Jeringas Dentales	340
Farabeuf Separadores	125
Farrell Espéculos Nasaes	131
Farrell Titular De Algodón	131
Ferguson-Ackland Abre Bocas	120
Flagg Pinzas Para Algodón Y Curaciones	140
Flohr Elevadores Para Raíces	64, 71, 72
Force Alicates Para Ortodoncia Y Prótesis	293
Fórceps De Extracción	1 - 52
Fórceps De Extracción (Forma Americana)	27-34
Fórceps De Extracción, Para Niños	35 - 38
Frahm Instrumentos Para Modelar Las Obturaciones	178
Frahm Talladores Para Cera Y Para Modelar	267
Frazier Tubos Para Saliva	228
Freer Instrumentos Para Periodoncia	196
Fret Sierra Para Modelos Y Hojas	285
Friedman Elevadores Para Raíces	70
Friedmann Pinzas Gubias	53
Furrer Instrumentos Para Coronas	222
G	
Gardner Instrumentos Para Periodoncia	192
Gerzog Martillo	194
Gillies Porta-Agujas	103, 110
Gillies Separadores	122
Glick Excavadores	177
Goldman-Fox Bisturios Para La Gingivectomia	201
Goldman-Fox Curetas	207
Goldman-Fox Instrumentos Para Periodoncia	194, 195
Goldman-Fox Pinzas Para Algodón Y Curaciones	138, 194
Goldman-Fox Sondas	147 - 149
Goldman-Fox Tijeras Para Encías	89, 90, 93
Gracey Curetas	205, 206
Gregg Instrumentos Para Modelar Las Obturaciones	177
Gritman Cuchillos Para Yeso	276
Gritman Talladores Para Cera Y Para Modelar	270
Gross-Maier Pinzas Portataponos	117
Gubias	53 - 58
Guthrie Separadores Finos	123
H	
Halsey Porta-Agujas	100, 107
Halsted - Mosquito Pinzas Hemostáticas	114
Harris Fórceps De Extracción (Forma Americana)	29
Hartmann Espéculos Nasaes	131
Hartmann Pinzas Gubias	54
Hartmann Pinzas Hemostáticas	114
Heath Tijeras Para Hilos Y Alambre	85
Hegar-Baumgartner Porta-Agujas	100, 107
Heidbrink Elevadores Para Fragmentos	77
Heidbrink Elevadores Para Raíces	67, 74

Heidemann Espátulas Para Separar	152, 153
Heister Abre Bocas Pinzas	119
Hemingway Curetas Cortantes	219
Henahan Fórceps De Extracción	50
Hirschfeld Curetas	210
Hollenback Instrumentos Para Obturaciones	171, 176
Hollenback Talladores Para Cera Y Para Modelar	267
Hopkins Instrumentos Para Periodoncia	196
Horn Martillas	284
How Alicates Para Ortodoncia Y Prótesis	295, 305, 312, 317, 323
Howard Limas Para Alveolotomia	216
Hull Fórceps De Extracción (Forma Americana)	32
Hyllin Talladores Para Cera Y Para Modelar	269, 270
I	
Implantologie	280 - 283
Indiana University Curetas	206
Instrumentos De Medición	286 - 289
Instrumentos Para Bandes	334
Instrumentos Para Cera	275
Instrumentos Para Coronas	220 - 225
Instrumentos Para Dique	158, 159, 166
Instrumentos Para Esmalte	186
Instrumentos Para Modelar La Ceramica	274
Instrumentos Para Modelar Las Obturaciones	177 - 180
Instrumentos Para Nervios	154 - 157
Instrumentos Para Obturaciones	168 - 184
Instrumentos Para Operatoria	185
Instrumentos Para Ortodoncia	333, 334
Instrumentos Para Periodoncia	191 - 212
Iris Tijeras Para Encías	89, 93
Ivory Porta Matrices	226
Iwanson Instrumentos De Medición	288
J	
Jameson Instrumentos De Medición	286
Jansen Pinzas Gubias	54
Jaquette Curetas	210
Jarabak Alicates Para Ortodoncia Y Prótesis	292, 319, 324
Jennings Abre Bocas	121
Jeringas Dentales	335 - 340
Johnson Alicates Para Ortodoncia Y Prótesis	303, 309, 326
Jones Pinzas De Campo	116
Juego De Instrumentos Para Dique	164
Juego De Instrumentos Para Implantología	143
Juego De Instruments Para Cirugía Oral	80
K	
Kells Fórceps De Extracción (Forma Americana)	31
Kelly Pinzas Hemostáticas	114
Kelly Tijeras Para Encías	91, 94
Kerrison Pinzas Gubias	58
Kilner Separadores Finos	125
Kirkland Bisturios Para La Gingivectomia	200
Kirkland Separadores	125
Klein Fórceps De Extracción	47, 48
Kleinert-Kutz Pinzas Gubias	53
Kliner Aguja De Reverdin Y Leznas	99
Knapp Separadores	123
Kocher Pinzas Hemostáticas	115
Kocher Separadores	127
König Instrumentos Para Periodoncia	197
Kopp Elevadores Para Raíces	62B
L	
La Grange Tijeras Para Encías	90
Ladmore Instrumentos Para Obturaciones	173

Lámparas De Alcohol	235	Micro-Adson Pinzas De Disección	135
Langenbeck Espéculos Nsales	131	Micro-Vascular Porta-Agujas	101, 108
Langenbeck Separadores	128	Middeldorpf Separadores	128
Langer Curetas	204	Miller Instrumentos Para Coronas	225
Laringoscopios	132, 133	Miller Limas Para Alveolotomia	215
Lawrence-Read Fórceps De Extracción	40, 47	Miller-Apexo Elevadores Para Raíces	65, 73
Le Cron Talladores Para Cera Y Para Modelar	269, 270	Mini-Friedman Pinzas Gubias	53, 194
Lecluse Elevadores Para Raíces	75	Mini-Ryder Porta-Agujas	101, 108
Lempert Pinzas Gubias	56	Molt Abre Bocas	120
Lessmann Instrumentos Para Cera	275	Molt Curetas Cortantes	217
Lezna Para Mandibula Superior	99	Molt Instrumentos Para Periodoncia	195
Lichtenberg Porta-Agujas	105, 112	Montfort Instrumentos Para Coronas	222
Limas Para Alveolotomia	215 - 216	Mosquito Pinzas Hemostáticas	114
Limpiadores De Sarro	213 - 214B	N	
Liston Pinzas Cortantes Para Huesos	59	Nabers Sondas	147
Littauer Pinzas Cortantes Para Huesos	59	Nance Alicates Para Ortodoncia Y Prótesis	300, 303, 308, 313, 324
Littauer Tijeras Para Hilos Y Alambre	85	Neumann Tijeras Para Encías	90
Locklin Pinzas Para Algodón Y Curaciones	138	Nevius Fórceps De Extracción (Forma Americana)	31
Locklin Tijeras Para Encías	91, 94	Niro Pinzas Cortantes Para Huesos	59
London Elevadores Para Raíces	62	Noyes (Castroviejo) Tijeras Para Encías	92
London Pinzas Para Algodón Y Curaciones	139	O	
Lorna (Edna) Pinzas De Campo	116	Obturadores De Amalgama	181
Lucas Curetas Cortantes	219	Ochsenbein Cinceles Y Gubias	191
Luer Pinzas Gubias	55	Ochsenbein Instrumentos Para Periodoncia	198
Luks Instrumentos Para Nervios	155	Olsen-Hegar Porta-Agujas	103, 109
Luniatschek Instrumentos Para Operatoria	185	Omega Alicates Para Ortodoncia Y Prótesis	294
M		Orban Bisturios Para La Gingivectomia	200
Mac Millan Elevadores Para Raíces	70	Orban Instrumentos Para Operatoria	185
Machatt Separadores Finos	123	P	
Maehlum Talladores Para Modelar	269	P.K Thomas Talladores Para Cera Y Para Modelar	271, 272
Magnm Alicates Para Ortodoncia Y Prótesis	299	Palmer Instrumentos Para Dique	158
Makellops Alicates Para Ortodoncia Y Prótesis	304	Parmly Fórceps De Extracción (Forma Americana)	29, 50
Mangos De Bisurí - Hojas De Bisturí	95, 96	Partsch Cinceles Y Gubias	193
Mangos De Elevadores Para Raíces	62A, 63	Partsch Instrumentos Para Periodoncia	197
Mangos Para Elevadores	63	Partsch Martillo	194
Mangos Para Elevadores Raices	62A	Pean Pinzas Hemostáticas	115
Mangos Para Espejos De Boca - Espejos De Boca	142	Peeso Alicates Para Ortodoncia Y Prótesis	331
Mangos Para Exploradores	144	Perry Pinzas Para Algodón Y Curaciones	139
Mangos Para Instrumentos Para Coronas	221	Physick Fórceps De Extracción (Forma Americana)	29
Marquardt Pinzas Gubias	58	Picos De Amalgama	181 - 183
Marquette Instrumentos Para Obturaciones	171	Pinceles Para Modelar La Ceramical	273
Martillos	284	Pinzas Abre Bocas	119
Martin Ergo Fórceps De Extracción	51	Pinzas Cortantes Para Huesos	59 - 61
Martini Curetas Cortantes	217	Pinzas De Campo	116
Mathieu Porta-Agujas	104, 111	Pinzas De Disección	136 - 137
Mayfield Pinzas Gubias	57	Pinzas De Los Brackets	333
Mayo Bandeja	332	Pinzas Gubias	53 - 58
Mayo Tijeras Para Cirugía	86	Pinzas Hemostáticas	114, 115
Mayo-Hegar Porta-Agujas	102, 109	Pinzas Para Algodón Y Curaciones	138 - 141
Mc Intosh Laringoscopios	132, 133	Pinzas Para Crisoles	290
McCall Curetas	208	Pinzas Para Esterilización	236
McCall Limpiadores De Sarro	214B	Pinzas Para Laboratorio	291
McKellops Alicates Para Ortodoncia Y Prótesis	308	Pinzas Para Papel De Articular	225
Mead Instrumentos Para Periodoncia	195	Pinzas Para Suture	134 - 137
Mead Instrumentos Para Periodoncia	193, 195	Pinzas Portatapones	117, 137 - 141
Mead Martillos	194	Polokoff Limas Para Alveolotomia	215
Mead Pinzas Gubias	53	Porta Aguja Con Carburo De Tungsteno	107 - 113
Meriam Pinzas Para Algodón Y Curaciones	140	Porta Matrices	226
Merrit Sondas Calibradas	148	Porta-Agujas	100 - 113
Mershon Instrumentos Para Dique	159	Porta-Amalgamas	166, 167
Metzenbaum Tijeras Para Cirugía	87	Potts Elevadores Para Raíces	62
Micro Curetas De Ligadura	322	Potts-Smith Pinzas Para Algodón Y Curaciones	138
Micro Tijeras	88	Prichard Instrumentos Para Periodoncia	195

Q		
Quinby Tijeras Para Encías	90	
R		
Ragnell Separadores	122	
Ramus Separadores	127	
Reverdin Agujas Y Leznas	99	
Reynolds Alicates Para Ortodoncia Y Prótesis	301, 327	
Rhodes Instrumentos Para Periodoncia	198	
Riñones	332	
Roach Talladores Para Cera Y Para Modelar	270	
Roser-Koenig Abre Bocas	119	
Ruskin Pinzas Gubias	58	
Ruskin-Liston Pinzas Cortantes Para Huesos	61	
S		
Sanders Bisturios Para La Gingivectomia	202	
Schaedle Pinzas De Campo	116	
Schluger Instrumentos Para Periodoncia	198	
Schmeckebier Elevadores Para Raíces	64	
Schroeder Bisturios Para Abscesos	98	
Schwarz Alicates Para Ortodoncia Y Prótesis	331	
Seldin Elevadores Para Raíces	65 - 66, 67, 73	
Seldin Instrumentos Para Periodoncia	197	
Semken-Taylor Pinzas De Disección	134	
Senn-Miller Separadores Finos	122	
Separadores	122 - 130	
Separadores De Los Labios	118	
Separadores Finos	122 - 125	
Sindesmótomos	78, 79	
Sondas	147 - 149	
Spencer Tijeras Para Hilos Y Alambre	85	
Spencer-Wells Pinzas Hemostáticas	115	
Spreader Instrumentos Para Nervios	154	
Stellbrink Pinzas Gubias	56	
Sternberg Abaja Lenguas	118	
Stevens Tijeras Para Cirugía	86	
Stieglitz Instrumentos Para Nervios	157	
Stocke Instrumentos Para Dique	158	
Stout Elevadores Para Raíces	66	
Straight Separador	122	
Sugarman Instrumentos Para Periodoncia	198	
Sullivan Tijeras Para Encías	91	
Swedish Porta-Agujas	101, 108	
T		
Talladores Para Cera Y Para Modelar	267 - 272	
Tanner Instrumentos Para Obturaciones	171	
Tarro Methot	332	
Taylor Curetas	209	
TG Cinceles	199	
Tijeras	81 - 94	
Tijeras Operaciones Fino	89	
Tijeras Para Cirugía	86, 87	
Tijeras Para Encías	89 - 92	
Tijeras Para Encías Con Carburo De Tungsteno	93, 94	
Tijeras Para Hilos Y Alambre	85	
Tijeras Para Vendajes	87	
Tijeras Quirúrgicas	81- 94, 158	
Toennis Porta-Agujas	103, 110	
Tofflemire Cintæ De Matrices	227	
Tofflemire Porta Matrices	226	
Tomati Pinzas Cortantes Para Encias	194	
Tomes Fórceps De Extracción (Forma Americana)	31	
Torquing Alicates Para Ortodoncia Y Prótesis	313	
Towner Jaquette Limpiadores De Sarro	213	
Towner Limpiadores De Sarro	213	
Townley Instrumentos De Medición	286	
Tubos Para Saliva	228, 229	
Tweed / Angel Alicates Para Ortodoncia Y Prótesis	314	
Tweed Alicates Para Ortodoncia Y Prótesis	292, 304, 306, 318, 324	
U		
Universal Alicates Para Ortodoncia Y Prótesis	300, 301	
Universal Instrumentos Para Coronas	223	
Universal Tijeras Para Cirugía	87	
Universal Tijeras Para Hilos Y Alambre	85	
University Of Michigan Elevadores Para Raíces	70	
University Of Michigan Sonidas Calibradas	148, 149	
University Of South California Bisturis Para La Gingivectomia	202	
USC Instrumentos Para Modelar Las Obturaciones	178	
USC Talladores Para Cera Y Para Modelar	267	
USC Towner Bisturis Para La Gingivectomia	202	
V		
Vaso Para Amalgama	166	
Vehe Talladores Para Cera Y Para Modelar	270	
Vernier Calibrador	286 - 289	
Vertical Condenser Instrumentos Para Nervios	155	
W		
Wahl Limas Para Alveolotomia	216	
Wakefield Cinceles Y Gubias	191	
Waldsachs Alicates Para Ortodoncia Y Prótesis	301, 315	
Walter F. Barry Elevadores Para Raíces	75	
Ward Instrumentos Para Obturaciones	171, 176	
Ward Talladores Para Cera Y Para Modelar	267	
Warwick-James Elevadores Para Fragmentos	77	
Wedelstaet Instrumentos Para Esmalte	186	
Weingar Alicates Para Ortodoncia Y Prótesis	294, 304, 312, 317, 320, 322, 328	
Weitlaner - Wullstein Separadores	129	
Weitlaner Separadores	129, 130	
Westcott Instrumentos Para Obturaciones	169, 173	
Westcott Tijeras Para Encías	92	
White Bruñidores	184	
Whitehead Abre Bocas	121	
Wiland Instrumentos Para Modelar Las Obturaciones	179	
Williams Sonidas	148, 149	
Williger Curetas Cortantes	217	
Williger Elevadores Para Raíces	76	
Williger Instrumentos Para Periodoncia	197	
Williger Separadores	126	
Witzel Fórceps De Extracción	51, 52	
Woodson Instrumentos Para Obturaciones	176	
Woodward Elevadores Para Raíces	69, 70	
Y		
Young Abre Bocas Pinzas	119	
Young Alicates Para Ortodoncia Y Prótesis	303, 330	
Young Arch	159	
Younger-Good Curetas	209	
Z		
Zaufal-Jansen Pinzas Gubias	57	
Zehle Talladores Para Cera Y Para Modelar	270	
Zerfing Cureta	203	
Zerfing Escaladores	214A	

Deutsch

Inhaltsverzeichnis

A		
Abdrucklöffel Aus Plastik, Einmaliger Gebrauch	261	
Abszeßmesser	98	
Acorn Kondensier Instrumente	184	
Adams Zangen	302, 309, 314, 330, 331	
Abdrucklöffel Aus Edelstahl	237 - 262	
Aderer Klemmen	304, 305, 308, 309	
	319, 325, 329, 331	
Adson Chirurgische Pinzetten	135	
Adson Nadelhalter	102	
Adson-Brown Chirurgische Pinzetten	135	
Ainsworth Koffferdam - Instrumente	158	
Alm Wundhaken	130	
Amalgam Stopfer	181 - 183	
Amalgambecher	166	
Amalgamträger	166, 167	
Angle Bird Zangen	304	
Apexo Wurzelheber	64, 67, 74	
Apical Wurzelheber	64	
Arruga Nadelhalter	105, 112	
Arterienklemmen	114, 115	
Artikulationspapierpinzetten	225	
Artikulatoren	263 - 266	
Ash Abszeßmesser	98	
B		
Backhaus Tuckklemmen	116	
Ball Füllungsinstrumente	173	
Bane Hohlmeißelzangen	56	
Barnhart Küretten	204	
Barraquer Mikro Scheren	88	
Barraquer Nadelhalter	106, 113	
Barry Wurzelheber	75	
Barsky Wundhaken Feine	122	
Baumgartner Nadelhalter	100, 107, 110	
Baynton Nadelhalter	105, 112	
Beale Instrumente Für Die Parodontologie	195	
Beale Wachs-Und Modellierinstrumente	270	
Bee Bee Fadenschere, Drahtschere	85	
Behrman Amalgam Stopfer	181	
Bein Wurzelheber	63, 71, 72	
Beyer Hohlmeißelzangen	56, 57	
Bird Zangen	292, 293, 302, 309	
Black Amalgam Stopfer	181 - 183	
Black Schmelzinstrumente	186	
Black Zahnerinigungs Instrumente	214A	
Blumenthal Hohlmeißelzangen	54	
Bodkin Amalgam Stopfer	182	
Bohler Hohlmeißelzangen	57	
Bohler Knochensplitterzangen, Übersetzt	61	
Böhm Kroneninstrumente	224	
Boley Meßbereich	287	
Boynton Nadelhalter	105, 112	
Brewer Koffferdam - Instrumente	158	
Brünings Zungenspatel	118	
Buchwald Zungenspatel	118	
Buck Feilen	198	
Buck Parodontalmesser	201	
Buckley Hohlmeißel	191	
C		
Castro Mikro Scheren	88	
Castroviejo (Barraquer) Nadelhalter	106, 113	
Castroviejo Meßbereich	288	
Castroviejo Nadelhalter	106, 113	
Castroviejo-Epker Meßbereich	289	
Cattoni Zahnsonden, Doppelendig	146, 147	
Cawood Minnesota Zungenspatel	118	
CF5A Füllungsinstrumente	178	
Chirurgische Pinzetten	136, 137	
Chirurgische Scheren	81- 94, 158	
Chompret Syndesmome	79	
Cleveland Hohlmeißelzangen	55	
Cleveland Knochensplitterzangen, Übersetzt	60	
Cogswell Wurzelheber	67	
Coleman Wurzelheber	62	
College Zahnpinzetten	141	
Collier Nadelhalter	100, 107	
Collin Mundsperrerr - Zungenzangen	119	
Columbia Instrumente Für Die Parodontologie	203	
Columbia Zungenspatel	118	
Cottle-Kazanjian Knochensplitterzangen, Übersetzt	61	
Coupland Wurzelheber	69, 70	
Crane Wurzelheber	70	
Crane-Kaplan Küretten	209	
Crane-Kaplan Parodontalmesser	202	
Crile Arterienklemmen	115	
Crile-Wood Nadelhalter	102, 109	
Cryer Wurzelheber	68, 71, 74	
Cryer Zahnzangen (Amerikanische Form)	32, 33, 51	
Curtis Wurzelheber	66	
D		
De La Rosa Zangen Für Kieferorthopädie Und Prothetik	305	
De Marco Küretten	206	
Dean Zahnfleischschere	90, 94	
Denhart Mundsperrerr - Zungenzangen	119	
Dental Spritzen	335 - 340	
Derf (Wright) Nadelhalter	100, 107	
Didner Matrizenspanner	226	
Dieffenbach Zahnanfängmesser	97	
Doyen-Collin Mundsperrerr	120	
Doyen-Jansen Mundsperrerr	120	
Drahtschneidezange	294 - 296, 298 - 300	
E		
Eiselberg-Mathieu Nadelhalter	103	
Endo-Container	233, 234	
Exkavatoren	187 - 190	
F		
Fadenschere, Drahtschere	85	
Falcon Zangen Für Kieferorthopädie Und Prothetik	310	
Falcon-Ject Dental Spritzen	340	
Farabeuf Wundhaken Feine	125	
Farrell Cotton Halter	131	
Farrell Nasenspekula	131	
Ferguson-Ackland Mundsperrerr	120	
Fine Swedish Nadelhalter	101, 108	
Flachmeißel	193	
Flachmeißel Für Knochen	193	
Flagg Zahnpinzetten	140	
Flexible Spachte	279	
Flohr Wurzelheber	64, 71, 72	
Force Zangen Für Kieferorthopädie Und Prothetik	293	
Frahm Instrumente Zum Modellieren Von Füllungen	178	
Frahm Wachs- Und Modellierinstrumente	267	
Frazier Speichelsauger	228	
Freer Instrumente Für Die Parodontologie	196	
Fret Säge Für Gipsmodelle Und Sägeblätter	285	
Friedman Wurzelheber	70	

Friedmann Hohlmeißelzangen	53	Instrumentenkästen	233 - 235, 332
Füllungsinstrumente	168 - 184	Instrumentensatz Für Implantation	143
Furrer Kroneninstrumente	222	Instrumentesatz Für Die Zahnärztliche Chirurgie	80
G		Instrumentos Para Periodoncia	215 - 216
Gardner Instrumente Für Die Parodontologie	192	Iris Zahnfleischscheren	89, 93
Gerzog Hammer	194	Ivory Matrizenspanner	226
Gillies Nadelhalter	103, 110	Iwanson Meßbereich	288
Gillies Wundhaken	122	J	
Gipsmesser	276	Jameson Meßbereich	286
Gipsschneidezange	279	Jansen Hohlmeißelzangen	54
Glas methot	332	Jaquette Küretten	210
Glick Exkavatoren	177	Jarabak Zangen Für Kieferorthopädie Und Prothetik	292, 319, 324
Goldman-Fox Instrumente Für Die Parodontologie	194, 195	Jennings Mundsperrerr	121
Goldman-Fox Küretten	207	Johnson Zangen Für Kieferorthopädie Und Prothetik	303, 309, 326
Goldman-Fox Parodontometer	147 - 149	Jones Tuckklemmen	116
Goldman-Fox Parodontalmesser	201	K	
Goldman-Fox Zahnfleischscheren	89, 90, 93	Kells Zahnzangen (Amerikanische Form)	31
Goldman-Fox Zahnpinzetten	138, 194	Kelly Arterienklemmen	114
Gracey Küretten	205, 206	Kelly Zahnfleischscheren	91, 94
Gregg Instrumente Zum Modellieren Von Füllungen	177	Keramikipinsel	273
Griffe für Aufzüge	62A	Kerrison Hohlmeißelzangen	58
Griffe Wurzelheber	62A	Kieferhöhlensonden	131
Gritman Gipsmesser	276	Kilderzahnzangen Nach Klein	35 - 38
Gritman Wachs- Und Modellierinstrumente	270	Kilner Wundhaken Feine	125
Gross-Maier Tupferzangen	117	Kirkland Parodontalmesser	200
Guthrie Wundhaken Feine	123	Kirkland Wundhaken Feine	125
H		Klammerbord	165
Halsey Nadelhalter	100, 107	Klein Zahnzangen	47, 48
Halsted - Mosquito Arterienklemmen	114	Kleinert-Kutz Hohlmeißelzangen	53
Hämmer	284	Klemmpinzetten, Zum Setzen Der Klebebrackets	333
Harris Zahnzangen (Amerikanische Form)	29	Kliner Reverdin Nadeln Und Ahlen	99
Hartmann Arterienklemmen	114	Knapp Wundhaken Feine	123
Hartmann Hohlmeißelzangen	54	Knochenfeilen	215 - 216
Hartmann Nasenspekula	131	Knochensplitterzangen, Übersetzt	59 - 61
Heath Fadenscheren, Drahtscheren	85	Kocher Arterienklemmen	115
Hegar-Baumgartner Nadelhalter	100, 107	Kocher Wundhaken	127
Heidbrink Wurzelheber	67, 74	Kofferdam-Klammern	160 - 163
Heidbrink Wurzelspitzenheber	77	Kofferdam - Instrumente	158, 159, 166
Heidemann Exkavatoren	151	Kondensier Instrumente	184
Heidemann Separierspatel	152, 153	Kondensier Instrumente	169, 173, 184
Heister Mundsperrerr - Zungenzangen	119	König Instrumente Für Die Parodontologie	197
Hemingway Kroneninstrumente	219	Kopp Wurzelheber	62B
Henahan Zahnzangen	50	Kronenaufbieger	220, 223
Hirschfeld Küretten	210	Kronenentferner	220
Hohlmeißelzangen	53 - 58	Kroneninstrumente	220 - 225
Hollenback Füllungsinstrumente	171, 176	Küretten	203 - 210, 217 - 219
Hopkins Instrumente Für Die Parodontologie	196	L	
Horn Hämmer	284	La Grange Zahnfleischscheren	90
How Zangen Für Kieferorthopädie Und Prothetik	295,	Laborpinzetten	291
	305, 312, 317, 323	Ladmore Füllungsinstrumente	173
Howard Knochenfeilen	216	Langenbeck Nasenspekula	131
Hull Zahnzangen (Amerikanische Form)	32	Langenbeck Wundhaken	128
Hygienist Scalers	214A	Langer Küretten	204
Hylin Wachs- Und Modellierinstrumente	269, 270	Laryngoskope	132, 133
I		Lawrence-Read Zahnzangen	40, 47
Implantoloia	280 - 283	Le Cron Wachs- Und Modellierinstrumente	269, 270
Indiana University Küretten	206	Lecluse Wurzelheber	75
Instrumente Für Bandaufsetzer	334	Lempert Hohlmeißelzangen	56
Instrumente Für Die Konservierende Behandlung	185	Lessmann Wachsinstrumente	275
Instrumente Für Die Orthodontie	333, 334	Lichtenberg Nadelhalter	105, 112
Instrumente Für Die Parodontologie	191 - 211	Lippenhalter	118
Instrumente Zum Modellieren Von Füllungen	177 - 180	Liston Knochensplitterzangen, Übersetzt	59
Instrumenteinsatz Für Die Parodontologie	211, 212	Littauer Fadenscheren, Drahtscheren	85

Littauer Knochensplitterzangen, Übersetzt	59
Locklin Zahnfleischscheren	91, 94
Locklin Zahnpinzetten	138
London Wurzelheber	62
London Zahnpinzetten	139
Lorna (Edna) Tuchklemmen	116
Lucas Kroneninstrumente	219
Luer Hohlmeißelzangen	55
Luks Nervinstrumente	155
Luniatschek Instrumente Für Die Konservierende Behandlung	185
M	
Mac Millan Wurzelheber	70
Machatt Wundhaken Feine	123
Maehlum Modellierinstrumente	269
Magnm Zangen Für Kieferorthopädie Und Prothetik	299
Makelops Zangen Für Kieferorthopädie Und Prothetik	304
Marquardt Hohlmeißelzangen	58
Marquette Füllungsinstrumente	171
Martin Ergo Zahnzangen	51
Martini Scharfe Löffel	217
Mathieu Nadelhalter	104, 111
Matrizenbänder	227
Matrizenspanner	226
Mayfield Hohlmeißelzangen	57
Mayo Chirurgische Scheren	86
Mayo Tablett	332
Mayo-Hegar Nadelhalter	102, 109
Mc Call Küretten	208
Mc Call Zahnerinigungs Instrumente	214B
Mc Intosh Laryngoskope	132, 133
Mckellops Zangen Für Kieferorthopädie Und Prothetik	308
Mead Hammers	194
Mead Hohlmeißelzangen	53
Mead Instrumente Für Die Parodontologie	195
Meriam Zahnpinzetten	140
Merrit Parodontometer	148
Mershon Koffferdam - Instrumente	159
Meßbereich	286 - 289
Messers	97, 98, 275, 276
Metzenbaum Chirurgische Scheren	87
Micro-Adson Chirurgische Pinzetten	135
Micro-Vascular Nadelhalter	101, 108
Middeldorp Wundhaken	128
Mikro Scheren	88
Miller Instrumentos Para Periodoncia	215
Miller Kroneninstrumente	225
Miller-Apexo Wurzelheber	65, 73
Mini-Friedman Hohlmeißelzangen	53, 194
Mini-Ryder Nadelhalter	101, 108
Modellierinstrumente Für Keramik	274
Molt Instrumente Für Die Parodontologie	195
Molt Mundsperrer	120
Molt Scharfe Löffel	217
Montfort Kroneninstrumente	222
Mosquito Arterienklemmen	114
Mundsperrer	119 - 121
Mundspiegel	142
Mundspiegelgriffe - Mundspiegel	142
Mundspiegelgriffe Kroneninstrumente	221
Mundspiegelgriffe Wurzelheber	63
Mundspiegelgriffe Zahnsonden	144
N	
Nabers Parodontometer	147

Nadelhalter	100 - 113
Nadelhalter Mit Hartmetalleinlagen	107 - 113
Nahtpinzette	134 - 137
Nance Zangen Für Kieferorthopädie Und Prothetik	300, 303, 308, 313, 324
Nasenspekula	131
Nervinstrumente	154 - 157
Neumann Zahnfleischscheren	90
Nevius Zahnzangen (Amerikanische Form)	31
Nieren-Tablett	332
Niro Knochensplitterzangen, Übersetzt	59
Norm-Tray Aus Rosfreiem Edeltah - Normgröße	230 - 232
Noyes (Castroviejo) Zahnfleischscheren	92
O	
Oberkiefer-Ahle	99
Ochsenbein Flachmeißel Und Hohlmeißel	191
Ochsenbein Instrumente Für Die Parodontologie	198
Olsen-Hegar Nadelhalter	103, 109
Omega Zangen Für Kieferorthopädie Und Prothetik	294
Operatin Schere Feine	89
Orban Instrumente Für Die Konservierende Behandlung	185
Orban Parodontalmesser	200
Ortho - Container	332
P	
P.K Thomas Wachs- Und Modellierinstrumente	271, 272
Palmer Koffferdam - Instrumente	158
Parmly Zahnzangen (Amerikanische Form)	29, 50
Parodontometer	147 - 149
Partsch Flachmeißel Und Hohlmeißel	193
Partsch Hammer	194
Partsch Instrumente Für Die Parodontologie	197
Pean Arterienklemmen	115
Peeso Zangen	331
Perry Zahnpinzetten	139
Physick Zahnzangen (Amerikanische Form)	29
Polokoff Instrumentos Para Periodoncia	215
Potts Wurzelheber	62
Potts-Smith Zahnpinzetten	138
Prichard Instrumente Für Die Parodontologie	195
Q	
Quinby Zahnfleischscheren	90
R	
Ragnell Wundhaken Feine	122
Ramus Wundhaken	127
Reverdin Nadeln Und Ahlen	99
Reynolds Zangen Für Kieferorthopädie Und Prothetik	301, 327
Rhodes Instrumente Für Die Parodontologie	198
Roach Wachs- Und Modellierinstrumente	270
Roser-Koenig Mundsperrer	119
Ruskin Hohlmeißelzangen	58
Ruskin-Liston Knochensplitterzangen, Übersetzt	61
S	
Sanders Parodontalmesser	202
Satz Zu Koffferdam-instrumente	164
Schaedle Tuchklemmen	116
Scharfe Löffel	217 - 219
Scheren	81 - 94
Schieblehren	286 - 289
Schluger Instrumente Für Die Parodontologie	198
Schmeckebier Wurzelheber	64
Schmelzinstrumente	186
Schroeder Abszeßmesser	98
Schwarz Zangen Für Kieferorthopädie Und Prothetik	331

Seldin Instrumente Für Die Parodontologie	197
Seldin Wurzelheber	65 - 67, 73
Semken-Taylor Chirurgische Pinzetten	134
Senn-Miller Wundhaken Feine	122
Skalpellgriffe - Skalpellklingen	95, 96
Sonden - Parodontometer - Kombination	149
Spatel Für Gips Und Alginate	277 - 278
Speichelsauger	228, 229
Spencer Fadenschere, Drahtschere	85
Spencer-Wells Arterienklemmen	115
Spirituslampe	235
Spreader Nervinstrumente	154
Stellbrink Hohlmeißelzangen	56
Sterilisierzangen	236
Sternberg Zungenspatel	118
Stevens Chirurgische Scheren	86
Stieglitz Nervinstrumente	157
Stocke Kofferdam - Instrumente	158
Stout Wurzelheber	66
Straight Wundhaken	122
Sugarman Instrumente Für Die Parodontologie	198
Sullivan Zahnfleischschere	91
Swedish Nadelhalter	101, 108
Syndesmotome	78, 79
T	
Tablett Escalador	230 - 232, 332
Tablett Für Instrumente	230 - 234
Tabletts Instrumente	230 - 232, 332
Tanner Füllungsinstrumente	171
Taylor Küretten	209
TG Meißel	199
Tiegelzangen	290
Toennis Nadelhalter	103, 110
Tofflemire Matrizenbänder	227
Tofflemire Matrizenspanner	226
Tomati Zahnfleischklappenstanze	194
Tomes Zahnzangen (Amerikanische Form)	31
Torquing Zangen Für Kieferorthopädie Und Prothetik	313
Towner Zahnreinigungsinstrumente	213
Townley Meßbereich	286
Tuchklemmen	116
Tupferzangen	117, 137 - 141
Tweed / Angel Zangen Für Kieferorthopädie Und Prothetik	314
Tweed Zangen Für Kieferorthopädie Und Prothetik	292, 304, 306, 318, 324
U	
Universal Chirurgische Scheren	87
Universal Fadenschere, Drahtschere	85
Universal Kroneninstrumente	223
Universal Zangen Für Kieferorthopädie Und Prothetik	300, 301
University Of Michigan Parodontometer	148, 149
University Of Michigan Wurzelheber	70
University Of South California Parodontalmesser	202
USC Instrumente Zum Modellieren Von Füllungen	178
USC Towner Parodontalmesser	202
USC Wachs- Und Modellierinstrumente	267
V	
Vehe Wachs- Und Modellierinstrumente	270
Verbandschere	87
Vernier Zehntelmaß	286 - 289
Vertical Condenser Nervinstrumente	155
W	
Wachs- Und Modellierinstrumente	267 - 272

Wachsinstrumente	275
Wahl Knochenfeilen	216
Wakefield Küretten	191
Waldsachs Zangen Für Kieferorthopädie Und Prothetik	301, 315
Walter F. Barry Wurzelheber	75
Ward Füllungsinstrumente	171, 176
Ward Wachs- Und Modellierinstrumente	267
Warwick-James Wurzelspitzenheber	77
Wedelstaet Schmelzinstrumente	186
Weingar Zangen Für Kieferorthopädie Und Prothetik	294, 304, 312, 317, 320, 322, 328
Weitlaner - Wullstein Wundhaken	129
Weitlaner Wundhaken	129, 130
Westcott Füllungsinstrumente	169, 173
Westcott Zahnfleischschere	92
White Kondensier Instrumente	184
Whitehead Mundsperrer	121
Wiland Instrumente Zum Modellieren Von Füllungen	179
Williams Parodontometer	148, 149
Williger Instrumente Für Die Parodontologie	197
Williger Scharfe Löffel	217
Williger Wundhaken	126
Williger Wurzelheber	76
Witzel Zahnzangen	51, 52
Woodson Instrumente Für Füllungen	176
Woodward Wurzelheber	69, 70
Wundhaken	122 - 130
Wurzelheber	62 - 76
Wurzelspitzenheber	77
Y	
Young Arch	159
Young Mundsperrer	119
Young Zangen Für Kieferorthopädie Und Prothetik	303, 330
Younger-Good Küretten	209
Z	
Zahnfleischmesser	97, 98
Zahnreinigungsinstrumente	213 - 214
Zahnfleischschere	89 - 92
Zahnfleischschere Mit Hartmetalleinlagen	93, 94
Zahnpinzetten	138 - 141
Zahnsonden	144 - 146
Zahnzangen	1 - 52
Zahnzangen (Amerikanische Form)	27 - 34
Zangen	157
Zangen Für Bandaufsetzer	295
Zangen Für Kieferorthopädie Und Prothetik	292 - 331
Zangen Für Ligatur	296
Zaufal-Jansen Hohlmeißelzangen	57
Zehle Wachs- Und Modellierinstrumente	270
Zementspateln	150
Zerfing Escaladores	214A
Zerfing Küretten	203
Zungenspatel	118

Instrument Care & Promotional Items

Most Common Instrument Problems.

Spotting (Light colored Spots)

Cause: Water droplets condensing on the instrument and evaporating very slowly. Usually traced to high mineral content in the water supply.

Solution: Follow the autoclave manufacturer's instructions exactly. Do not open the door until all the steam has been exhausted. Check for leaky gaskets or valves which may cause condensation during the drying cycle.

Corrosion

Cause: Insufficiently rinsed or linens containing residue of caustic detergent chemicals.

Solution: Check with the laundry to be sure linens are rinsed thoroughly.

Cause: Baked-on blood in the box locks, serrations and ratchets.

Solution: Open all box locks completely prior to cleaning

Rust Colored Film

Cause: Water softeners can cause rust colored, dust like film to occur under certain conditions. Also, if the water supply has a high iron content, mineral deposits on the instruments may result.

Solution: Consult with the hospital engineering staff for an evaluation of your water treatment program.

Pitting

Cause: Exposure to saline solutions, blood potassium chloride and a number of other compounds.

Solution: Rinse with distilled water as soon as possible after exposure.

Cause: Some detergent contain chloride or an acid base that, if not thoroughly rinsed from the instruments, unites with steam in the autoclave to form hydrochloric acid.

Solution: Avoid detergents with chloride bases and select a detergent with a pH range no higher than 8.5.

Staining (Purplish-Black Stain)

Cause: Amine chemicals, which are used in steam lines to clean lime deposits, can be deposited on stainless steel by electrolytic action.

Solution: Use distilled or de-mineralized water instead of the local water source to cycle through the autoclave.

Eye & E.N.T

Microsurgical instruments are most vulnerable to damage through misuse, abuse & rough handling.

Since their exacting performance is extremely crucial to microsurgery, extra attention should be given to the care of microsurgical instruments.

How to care for MicroSurgical Instruments

1. All microsurgical instruments should be carefully inspected when first purchased and after each use. Ideally a magnifying glass should be used to help in detecting flaws not visible to the naked eye. Look for burs on tips, nicks on cutting edges and proper alignment of jaws.

2. Extra care should be taken during surgery wipe off all blood and debris after each use. Non-fibrous sponge should be used for this purpose to prevent snagging and breaking delicate tips.

3. All microsurgical instruments should be clean by hand and thoroughly dried before packing or storing. Avoid using a washer-sterilize ultrasonic cleaner.
4. Metal to metal contact should always be avoided. Special racks and cases are available for many delicate instruments and should be used to protect and separate tips.
5. Preventative maintenance should be done on all micro instruments. Regular sharpening, realignment and adjustment can save costly repairs of irreparable damage.
6. Microsurgical instruments should only be used for their intended purpose. The delicate construction can be easily damaged if used in cases other than microsurgery.

Sterilization

All surgical instruments must be sterilized prior to surgery to prevent infection.

But even sterilization can leave contaminants behind if instruments are not properly cleaned.

To prevent this, sterilize instruments with the ratchets open. This allows for better steam penetration.

Plus it prevents the box locks (hinge area) from cracking.

If using a pan or tray, we recommend perforated. This will also enable better steam penetration, and aids in more effective drying as well. For efficiency, place heavy instruments at the bottom and lighter, more delicate instruments on top.

If sterilizing in paper or plastic pouches, do not stack pouches on top of one another during sterilization. When possible, use a spiral metal letter holder as a standing aid. This will permit proper steam flow.

Ultrasonic Cleaning

Tip : A heated tank does not significantly improve cleaning. Generally speaking, there is no measurable improvement in cleaning-- or timing -- when heat is applied to the solution. So, a heater is not a necessary feature in choosing an ultrasonic cleaner.

A method of cleaning that is growing in popularity is ultrasonic cleaning. This method is, by far, the most efficient and effective available today. Its ease of use and superior efficiency is quickly making ultrasonic cleaning the preferred choice for today's surgeons.

In fact, ultrasonic cleaning is 16 times more efficient than hand-cleaning. Place instruments in the ultrasonic unit for 10-15 minutes and use a neutral pH solution. Here are a few more tips for ultrasonic cleaning:

Before placing into the ultrasonic unit, clean instruments of all visible debris by washing them in an approved instrument cleaning solution.

Don't mix dissimilar metals (such as aluminum and stainless) in the same cycle. Make sure instruments have plenty of room. Don't overload your ultrasonic cleaner.

As with all types of cleaning, open all instruments so ratchets & box locks are fully exposed to the cleaning process.

Upon completion of the cycle, remove instruments immediately, and rinse them. Dry thoroughly with a towel, ensuring that no moisture is left on the instruments.

The use of water in an ultrasonic cleaner is not recommended. A neutral pH ultrasonic cleaner solution, when properly mixed, effectively reduces the surface tension of the solution and increases the ultrasonic cavitation process.

The solution should be changed at least daily, or sooner if the solution appears dirty or murky.

Steam Autoclave

Tip: Before autoclaving, lubricate all instruments that have moving parts. Use surgical lubricants rather than industrial oils, as they are superior for surgical instruments.

Most hospitals and practices use the Steam Autoclave for sterilizing instruments. If this is your method, there are a few important points to remember.

1. Only use distilled water for filling your sterilizer reservoir. Tap water contains minerals that will cause staining, if left to stand or dry on the instrument.
2. If your autoclave has a steam line filter, clean it regularly using the manufacturer's recommendations as a minimum guide.

Clean the inside of the chamber once a week. This will prevent build-up of scale and will allow the sterilizer to operate more efficiently. To clean the inside chamber, follow these steps:

Turn the unit off and allow chamber to cool. Then remove trays and racks.

Using towels, remove any water laying in the bottom of the autoclave.

Spray the inside walls of the unit with an approved solution, such as Spectrum's **Spectra-Clave™**. Allow the solution to stand five minutes (or as instructed on the bottle).

Using a stiff, nylon-bristle brush, brush the inside chamber very aggressively. After brushing, wipe the inside clean with clean towels.

Finally, using a moist towel, wipe the gasket clean on the inside of the door.

Lubrication

Tip: It is no longer recommended to use a lubricant bath because the container of the lubricant solution may contain certain bacteria from previous instruments dipped into this bath. A lubricant spray is advised. Sprays are safer, cost less, and take up less counter space as well.

One of the easiest, yet most effective ways to keep instruments in excellent condition is to lubricate them after every cleaning. Proper lubrication keeps instruments from rubbing and scraping, thus preventing dulling and strain to joints and hinges. Moving parts on instruments, such as joints, box locks,

ratchets, and screw joints, should be lubricated regularly. Before autoclaving, lubricate all instruments that have moving parts. Only use surgical lubricants because they are steam-penetrable.

Wrapping Materials, Surgical Towels and Drapes

Wrapping of surgical instruments is essential in the sterilization process. If your practice uses re-usable towels and drapes, please be certain to use as little laundry detergent as possible. Towels and drapes can retain soap particles. During the autoclave cycle, steam passes through the fabric, picking up these soap particles and depositing them on the surface of the instruments.

One suggestion for preventing this is to run an extra rinse cycle to remove excess soap particles. Also, if your instrument packs are coming out wet, we suggest a towel be placed inside the pack to absorb moisture.

Rust vs. Stain

Stains can be removed, whereas rust will leave permanent damage. To determine if a brown or orange discoloration is a stain or rust, use the eraser test. Rub a pencil eraser over the discoloration. If the discoloration is removed with the eraser and the metal underneath is smooth and clean, this is a stain. If a pit mark appears under the discoloration, this is corrosion or rust.

Trouble-Shooting Stain Guide for Surgical Instruments

Stain Color: Brown / Orange

Cause: A result of high-pH detergents, Chlorhexidine usage, or improper soaking of instruments. This color stain can also be caused by soaking in tap water.

Stain Color: Dark Brown

Cause: Low-pH instrument solutions. The brownish-colored film may also be caused by a malfunctioning sterilizer. Similar localized stain spots can also be a result of baked-on blood.

Stain Color: Bluish Black

Cause: Reverse plating, when instruments of different metal (e.g., chrome and stainless steel) are ultrasonically processed together. Additionally, exposure to saline, blood, or potassium chloride will cause this bluish black color.

Stain Color: Multicolor

Cause: Excessive heat by a localized hot spot in the sterilizer. The rainbow-colored stain can be removed.

Stain Color: Light and Dark Colored Spots

Cause: Water droplets drying on the instruments. With slow evaporation, the minerals sodium, calcium, and magnesium left behind can cause this spotting.

Stain Color: Bluish Gray

Cause: Liquid (cold) sterilization solutions being used beyond manufacturer's recommendations.

Stain Color: Black Stains

Cause: Contact with ammonia or a solution containing ammonia.

Stain Color: Gray

Cause: A liquid rust remover being used in excess of manufacturer's recommendations.

Stain Color: Rust

Cause: Dried blood that has become baked on the serrated or hinged areas of surgical instruments. This organic material, once baked on, may appear dark in color. Also can be caused by soaking in tap water.

Scissors Testing

Until recently, it's been difficult to gauge if scissors were performing optimally. To test them, surgeons often tried cutting a piece of material, but this yielded inconsistent results due to the variety of materials used. Today, there is a rubber-like material known as **Theraband** to test scissors' sharpness.

Medium-grade red **Theraband** is used for scissors measuring 4.5 inches to 12 inches in length.

Yellow **Theraband** is used for scissors measuring 3 inches to 4 inches in length.

Make several complete cuts through the material, cutting all the way to the distal tip (where scissors most often dull). This test will indicate satisfactory sharpness. The scissor test material is available from Spectrum.

Cold Sterilization

This is not a recommended method of disinfecting for high-quality surgical instruments. Immersing instruments for long periods of time in any solution can prove damaging. Never leave instruments in any solution for longer than 20 minutes. Do not immerse instruments with tungsten carbide inserts (gold handles) in solutions containing Benzyl Ammonium Chloride (BAC). This chemical is known to loosen the tungsten carbide.

Warnings & Precautions

Universal Precautions should be observed :

By all hospital personnel that work with contaminated or potentially contaminated medical devices. Caution should be exercised when handling devices with sharp points or cutting edges.

Personal Protective Equipment (PPE) should be worn :

When handling or working with contaminated or potentially contaminated materials, devices and equipment. PPE includes gown, mask, goggles or face shield, gloves and shoe covers.

Metal brushes or scouring pads must not be used :

During manual cleaning procedures. These materials will damage the surface and finish of instruments. Soft-bristled, nylon brushes and pipe cleaners should be used.

Cleaning agents with low foaming surfactant should be used during manual cleaning procedures to ensure that instruments are visible in the cleaning solution. Manual scrubbing with brushes should always be performed with the instrument below the surface of the cleaning solution to prevent formation of aerosols and splashing which may spread contaminants.

Cleaning agents must be easily and completely rinsed from device surfaces to prevent accumulation of detergent residue.

Do not place heavy instruments on top of delicate devices.

Do not allow contaminated devices to dry prior to reprocessing :

All subsequent cleaning and sterilization steps are facilitated by not allowing blood, body fluid, bone and tissue debris, saline, or disinfectants to dry on used instruments.

Saline and cleaning/disinfection agents containing aldehyde, mercury, active chlorine, chloride, bromine, bromide, iodine or iodide are corrosive and **should not** be used. Instruments **must not** be placed or soaked in **Ringers Solution**.

Mineral oil or silicone lubricants **should not** be used because they:

- 1) Coat microorganisms
- 2) Prevent direct contact of the surface with steam; and
- 3) are difficult to remove.

Only devices manufactured and/or distributed by Us should be included in Our instrument trays and cases. These validated reprocessing instructions are **not applicable** to Our trays and cases that include devices that are not manufactured and/or distributed by Us.

Descaling agents that include morpholine should not be used in steam sterilizers. These agents leave residue which can damage polymer instruments over time.



WD-364-001
Scalpel Handle



WD-364-002
Operating Scissors



WD-364-003
Dressing Forceps



Mathieu
WD-364-004
Needle Holder



Pean Forceps
WD-364-005



Kelly Forceps
WD-364-006



Mayo-Hegar
WD-364-007
Needle Holder




WD-364-008
Nasal Speculum



WD-365-009
Bone Curettes



Mead
WD-365-010
Bone Rongeur



WD-365-011
Cuticle Nipper



WD-365-012
Rongeur



WD-365-013



WD-365-014



WD-365-015



WD-365-016



WD-366-017



WD-366-018



WD-366-019



WD-366-020



WD-366-021



WD-366-022
Caliper

Weldon Instruments
Ganjay Wali Street
Rai Pur Stop, Wazirabad Road
Sialkot-51310, Pakistan.
Phone: +92 52 3571383 & 3571384
Fax : +92 52 3571382

E-mail: junaid@weldoninstruments.com
waqas@weldoninstruments.com



<https://www.facebook.com/weldoninstruments>

www.weldoninstruments.com